- EN**
3 **Ultrathane® Endoureterotomy Stent Set**
Instructions for Use
- CS**
5 **Souprava endoureterotomického stentu z materiálu Ultrathane®**
Návod k použití
- DA**
7 **Ultrathane® stentsæt til endoureterotomi**
Brugsanvisning
- DE**
10 **Endoureterotomie-Stentset aus Ultrathane®**
Gebrauchsanweisung
- EL**
13 **Σετ ενδοπρόσθεσης ενδοουρητηροτομής από Ultrathane®**
Οδηγίες χρήσης
- ES**
15 **Equipo de stent de endoureterotomía de Ultrathane®**
Instrucciones de uso
- FR**
18 **Set d'endoprothèse pour endo-urétérotomie en Ultrathane®**
Mode d'emploi
- HU**
21 **Ultrathane® endoureterotomiás sztenkészet**
Használati utasítás
- IT**
23 **Set con stent per endoureterotomia in Ultrathane®**
Istruzioni per l'uso
- NL**
26 **Ultrathane® endo-ureterostomie stentset**
Gebruiksaanwijzing
- NO**
28 **Stentsett av Ultrathane® for endoureterotomi**
Bruksanvisning
- PL**
31 **Zestaw stentu do endoureterotomii z Ultrathane®**
Instrukcja użycia
- PT**
33 **Conjunto de stent para endoureterotomia em Ultrathane®**
Instruções de utilização
- SV**
36 **Ultrathane®-stentset för endoureterotomi**
Bruksanvisning

ULTRATHANE® ENDOURETEROTOMY STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

This set contains the following components:

- Radiopaque double pigtail stent
- Coaxial positioner - consisting of an inner stent positioner and outer release sleeve

The stent has a segment with an expanded outer diameter. This segment is designed to maintain dilation in the associated section of the ureter.

NOTE: A .038" PTFE wire guide is recommended for use with these stents.

INTENDED USE

This device is intended for temporary internal drainage from the ureteropelvic junction to the bladder following dilation and/or incision of a stricture. Double pigtail ureteral stents have been used to relieve obstruction in a variety of benign, malignant, and post-traumatic conditions. These stents may be placed using endoscopic, percutaneous, or open surgical techniques.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- Complications of ureteral stent placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures.
- Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to Cook.
- Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.
- These stents must not remain indwelling more than two months. If the patient's status permits, the stent may be replaced with a new stent.
- These stents are not intended as permanent indwelling devices.
- Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered.
- A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the wire guide or stent should be avoided.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.
- If problems occur while using this device, please call your Cook sales representative or contact our Customer Service department.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Ultrathane Endoureterotomy Stent Set is **MR Conditional** according to ASTM F2503. A patient with this device can be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum magnetic field spatial gradient of 1900 gauss/cm (19.0 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of < 4.0 W/kg (Normal Operating or First Operating Mode)

Nonclinical testing demonstrated that the Ultrathane Endoureterotomy Stent is not expected to induce RF heating greater than that of biological tissue.

The image artifact extends approximately 1.5 mm from the Ultrathane Endoureterotomy Stent as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

POTENTIAL ADVERSE EVENTS

Complications of ureteral stent placement are documented. These complications include, but are not limited to:

- Extravasations
- Occlusion
- Migration
- Hemorrhage
- Sepsis
- Perforation of the urinary tract
- Peritonitis
- Encrustation
- Urinary tract infection
- Loss of renal function

INSTRUCTIONS FOR USE

Endoscopic / Percutaneous Placement

1. Pass the stent over a previously placed wire guide. **NOTE:** Ensure that the tapered enlarged segment will be placed in the desired location for treatment.
2. Pass the coaxial positioner over the wire guide and advance it until it touches the tip of the stent.
3. Loosen the compression fitting.
4. Stabilize the stent positioner and advance the release sleeve over the non-enlarged portion of the stent until it abuts the taper of the enlarged section.
5. Lock the compression fitting.
6. Advance the assembly into the desired position while monitoring via fluoroscopy.
7. Remove the wire guide.
8. Verify that the stent curl is properly positioned within the target location.
9. Release the compression fitting and retract the release sleeve over the stent positioner.
10. Remove the coaxial positioner.
11. Verify that the stent curl is properly positioned within the target location. **NOTE:** If necessary, final adjustment can be made with endoscopic forceps. The stent may be removed by gentle traction using endoscopic forceps. Fluoroscopy facilitates stent placement; however, standard radiography may be used.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

R. Clayman, M.D., Washington University School of Medicine, St. Louis, Missouri, U.S.A.

J. Hulbert, M.D., University of Minnesota, Minneapolis, Minnesota, U.S.A.

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "Self Retained Internal Ureteral Stents: A New Approach," *The Journal of Urology*, 119 (1978), 731.

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "The Pigtail Ureteral Stent in the Cancer Patient," *The Journal of Urology*, 121 (1979), 17.

M. F. Camacho, R. Pereira, H. Carrion, M. Bondhus, V. A. Politano: "Double-Ended Pigtail Ureteral Stent: Useful Modification To Single Ended Ureteral Stent," *Urology*, 13 (1979), 516-520.

M. A. Koyle, A. D. Smith: "Ureteral Stents," Chapter 31 in A. D. Smith, W. R. Castañeda-Zuñiga, J. G. Bronson: *Endourology: Principles and Practice*. Thieme Inc., New York, New York, 1986.

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J. Hulbert, M.D., University of Minnesota, Minneapolis, Minnesota, USA.

T.W. Hepperlen, H.K. Mardis, H. Kammandel: "Self Retained Internal Ureteral Stents: A New Approach," *The Journal of Urology*, 119 (1978), 731.

T.W. Hepperlen, H.K. Mardis, H. Kammandel: "The Pigtail Ureteral Stent in the Cancer Patient," *The Journal of Urology*, 121 (1979), 17.

M.F. Camacho, R. Pereiras, H. Carrion, M. Bondhus, V.A. Politano: "Double-Ended Pigtail Ureteral Stent: Useful Modification To Single Ended Ureteral Stent," *Urology*, 13 (1979), 516-520.

M.A. Koyle, A.D. Smith: "Ureteral Stents," Chapter 31 in A.D. Smith, W.R. Castaneda-Zuniga, J.G. Bronson: *Endourology: Principles and Practice*. Thieme Inc., New York, New York, 1986.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd