Endoureterotomy Stent Set

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T _ E N D O S _ R E V O

ENDOURETEROTOMY STENT SET

Used for temporary internal drainage from the ureteropelvic junction to the bladder following incision of a stricture. The segment with no sideports prevents postoperative narrowing of the ureteral lumen while preventing ingrowth of the ureteral wall to the stent. Double pigtail ureteral stents have been employed to relieve obstruction in a variety of benign, malignant, and post-traumatic conditions. These stents may be placed using endoscopic, percutaneous, or open surgical techniques.

CAUTION: Complications of ureteral stent placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures.

NOTE: Prior to use, immerse AQ® stent in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

NOTE: If no wire guide is provided with this set, the following is recommended:

5.0 French stents accept .035" PTFE wire guide,

6.0 French stents accept .038" PTFE wire guide,

7.0 French stents accept .038" PTFE wire guide.

SET CONTAINS *(set components may vary)*

- Radiopaque double pigtail stent
- Radiopaque release sleeve with lockable fitting
- Radiolucent stent positioner with radiopaque band

CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.

CONTRAINDICATIONS

There are no known contraindications.

SUGGESTED INSTRUCTIONS FOR USING ENDOURETEROTOMY STENT SET

ENDOSCOPIC PLACEMENT

1. Pass a flexible wire guide tip beyond the obstruction to the renal pelvis. Tortuosity in the obstructed ureter often can be resolved using a wire guide and an open-end ureteral catheter in combination.
2. Using a baseline pyelogram, estimate the proper stent length; add 1 cm to that estimated ureteral measurement. Accurate measurement enhances drainage efficiency and patient comfort.
3. Pass the stent over the wire guide through the cystoscope. Under direct vision, advance the stent into the ureter with the positioner/release sleeve assembly. **NOTE: If needed, loosen the release sleeve from the positioner and gently advance the release sleeve to push the stent into position.** Have an assistant hold the wire guide in position to prevent advancement of the wire guide into the renal parenchyma.
4. Watch for the distal end of the stent at the ureterovesical junction. At that point, halt advancement of the stent. As an assistant removes the wire guide, hold the stent in position with the positioner/release sleeve assembly. The stent pigtail will form spontaneously. Carefully remove the positioner/release sleeve assembly from the cystoscope. **NOTE: If necessary, final adjustment can be made with endoscopic forceps. The stent may be removed easily by gentle withdrawal traction using endoscopic forceps. Fluoroscopy facilitates stent placement; however, standard radiography may be used.**

CAUTIONS

Intended for one-time use only.

C-Flex® and Ultrathane® stents must not remain indwelling more than two (2) months. If the patient's status permits, the stent may be replaced with a new stent.

These stents are not intended as permanent indwelling devices.

Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered.

A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.

Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the wire guide or stent should be avoided. Use of a 0-degree scope lens is recommended. Scopes larger than 21.0 French are suggested.

Individual variations of interaction between stents and the urinary system are unpredictable.

Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.

If problems occur using this device, please call your Cook® sales representative or contact our Customer Quality Assurance department at the address/phone listed within this booklet.

REFERENCES

R. Clayman, M.D., Washington University School of Medicine, St. Louis, Missouri, U.S.A.

J. Hulbert, M.D., University of Minnesota, Minneapolis, Minnesota, U.S.A.

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "Self Retained Internal Ureteral Stents: A New Approach," *The Journal of Urology*, 119 (1978), 731.

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "The Pigtail Ureteral Stent in the Cancer Patient," *The Journal of Urology*, 121 (1979), 17.

M. F. Camacho, R. Pereira, H. Carrion, M. Bondhus, V. A. Politano: "Double-Ended Pigtail Ureteral Stent: Useful Modification To Single End Ureteral Stent," *Urology*, 13 (1979), 516-520.

M. A. Koyle, A. D. Smith: "Ureteral Stents," Chapter 31 in A. D. Smith, W. R. Castañeda-Zuñiga, J. G. Bronson: *Endourology: Principles and Practice*. Thieme Inc., New York, New York, 1986.

C-FLEX® is a registered trademark of Concept Polymer Technologies, Inc.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.