Feldman Aortic Stenosis Catheter

Instructions for Use

EN
4

Feldman aortastenosekateter

Brugsanvisning

DA
5

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T _ F A S C _ R E V . 2


Fig. 1

FELDMAN AORTIC STENOSIS CATHETERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Feldman Aortic Stenosis Catheter is a nylon catheter with stainless steel braiding for torque control. The catheter tip has two opposing sideports and a tapered endhole for use during ventriculography. The catheter shape design, coupled with the use of a straight movable core wire guide, facilitates crossing of stenotic aortic valves.

INTENDED USE

The Feldman Aortic Stenosis Catheter is intended for ventriculography.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

None known

WARNINGS

- **Do not exceed recommended flow rate.** Flow rates in excess of 10 ml/sec have been reported to cause ventricle perforation and intramuscular injection with improperly positioned catheters. Recommended flow rates of 6 ml/sec for a total volume of 25-35 ml is sufficient for most patients.

PRECAUTIONS

- **Carefully position catheter tip in the mid cavity of the left ventricle.** Obtain images in the RAO and LAO views to ensure proper positioning toward the apex of the ventricle. A test injection of contrast at a rate of 6 ml/sec for a total volume of 5-10 ml is recommended to prevent staining of the left ventricle prior to power injection.
- Manipulation of catheter requires fluoroscopic control.
- Due to thinwall construction, extreme care must be exercised during manipulation and withdrawal. Catheter insertion through synthetic grafts or severely scarred areas should be avoided. If necessary, use of an introducer sheath is recommended.
- Do not attempt to heat or reshape the catheter curve; the catheter tip is made from a heat-sensitive material.
- The possible whiplash effect of the long, soft catheter tip must be considered during selective angiography.
- Activate hydrophilic coating, if present, by wetting the distal end of the catheter with sterile water or saline. For best results, keep the catheter surface wet during placement.
- If an over-the-wire catheter design is used, where the wire is left in place during injections, the proximal end of the wire should be clamped to protect against accidental injection of the wire guide.

PRODUCT RECOMMENDATIONS

Recommended Flow Rates

Use of 6 ml/sec for a total volume of 25-35 ml is adequate for most patients. When the ventricle is extremely dilated, 10 ml/sec may be needed.

INSTRUCTIONS FOR USE

1. Over a straight wire guide, insert the catheter into the aortic root and rotate it clockwise to direct the catheter tip toward the center of the aortic root. No more than 90-120 degrees of rotation is necessary.
NOTE: A non-tapered, straight, movable-core wire guide may be used to change the angle of the catheter, allowing the operator to scan or interrogate the surface of the aortic valve. (**Fig.1**) By withdrawing the core of the movable core wire guide 3-4 inches, the tip of the wire guide will be made more flexible, thereby allowing the catheter to assume its formed curve completely. When the aortic root is small in diameter, the catheter can be straightened with the wire and formed into a less-curved shape which approximates a right Judkins curve. The angled design will reach the center of the aortic root in most patients unless the left ventricle and aortic root meet in an extremely acute angle.
2. Once the tip of the catheter and wire guide have passed through the stenotic valve, advance the catheter into the mid-portion of the left ventricle with the catheter tip positioned toward the apex.
3. Remove the wire guide and perform a test injection of 6 ml/sec for a total volume of 5-10 ml to confirm proper catheter positioning within the left ventricle. Both RAO and LAO projections should be checked prior to power injection. If the catheter tip is positioned against the ventricle wall, reposition the catheter so that it is pointing toward the apex and repeat test injection.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from

package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.