Fusion™ Biliary Dilation Balloon

Instructions for Use

**EN
4**

Fusion™ galdevejsdilatationsballon

Brugsanvisning

**DA
5**

Fusion™ Gallengang-Dilatationsballon

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**DE
7**

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FUSION™ BILIARY DILATION BALLOON

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Fusion™ Biliary Dilation Balloon is a multi-lumen catheter with a balloon near its distal tip. The manifold consists of two independent hubs: "Wire Guide" and "Inflation." The wire guide lumen, with the Zip Port, facilitates placement of wire guides. The inflation lumen is used to expand the balloon. (Fig. 1)

INTENDED USE

The Fusion™ Biliary Dilation Balloon is intended to dilate strictures of the biliary tree. The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for biliary access and dilation should be employed.

CONTRAINDICATIONS

Contraindications to dilation are those specific to ERCP and any procedures to be performed in conjunction with balloon dilation. They include, but are not limited to:

- Uncooperative patient
- Asymptomatic strictures
- Inability to advance the dilation balloon through the strictured area
- Coagulopathy
- Known or suspected perforation
- Severe inflammation or scarring near the dilation site

POTENTIAL ADVERSE EVENTS

Potential adverse events are those associated with ERCP. They include, but are not limited to:

- Pancreatitis
- Cholangitis
- Aspiration
- Perforation
- Hemorrhage
- Infection
- Sepsis
- Allergic reaction to contrast medium or medication
- Hypotension
- Respiratory depression or arrest
- Cardiac arrhythmia or arrest

WARNINGS

- Do not exceed the recommended balloon inflation pressure as listed on the product label or catheter tag of balloon.
- Do not advance the balloon if resistance is encountered. Stop and assess cause of resistance to determine if dilation should be re-attempted.

PRECAUTIONS

- Refer to package label for the minimum channel size required for this device.
- Do not inflate the balloon prior to introduction into the scope, as this may cause damage to the scope.
- The duodenoscope must remain as straight as possible with elevator in open position when advancing or withdrawing balloon.
- The entire balloon must be extended beyond the tip of the scope and be fluoroscopically visualized and positioned before inflation.
- The balloon is used in conjunction with an inflation device or manometer and a 1:1 mixture of contrast medium and normal saline. **DO NOT** use air or gaseous substances to inflate the balloon, as this will result in reduced balloon effectiveness.
- While withdrawing the balloon from the scope, negative pressure must be applied to maintain balloon deflation.
- A compromised balloon may prohibit removal from the accessory channel. Removal of the scope, along with the compromised balloon, may be required.

INSTRUCTIONS FOR USE

Illustrations

System Preparation

1. Attach the inflation device to the inflation hub.
2. Create and maintain a vacuum with the inflation device. The balloon is now ready for placement through the accessory channel.
3. Remove the protective sheath from the balloon. Apply a water-soluble lubricant to the balloon to allow easier passage through the accessory channel.

If Using Zip Port and Pre-Positioned Short Wire Guide:

1. Advance the tip of the balloon catheter over the pre-positioned wire guide, ensuring that the wire guide exits the catheter at the Zip Port.
2. Unlock the short wire guide from the Wire Guide Locking Device.
3. Introduce the deflated balloon into the accessory channel and relock the wire guide.
4. Continue to advance the device in short increments until it can be endoscopically visualized exiting the scope.
5. Monitoring fluoroscopically, continue advancing the device into the duct until the desired position is achieved.
6. Inflate balloon to specified inflation pressure, as listed on the tag attached to the device. **NOTE:** Position may be verified using fluoroscopic markers located at either end of the balloon.
7. To deflate the balloon, decrease pressure below 20 psi, then create and lock in negative pressure. **NOTE:** The indicator on the pressure gauge should be in the vacuum portion during deflation.
8. Maintain negative pressure and observe the balloon for complete deflation. If the balloon does not completely deflate, create negative pressure again, as indicated in step 7. Repeat this step as necessary until the balloon is completely deflated. **CAUTION:** The balloon must be completely deflated, with all fluid removed, before withdrawal.
9. **If wire guide is to remain in place during device withdrawal:**
 - a. To remove deflated balloon from accessory channel, straighten the scope tip and open the elevator. Maintain negative pressure while withdrawing the balloon, using a continuous twisting motion.
 - b. Once resistance is achieved, unlock the wire guide from the Wire Guide Locking Device and perform wire guide exchange.
 - c. Lock the wire guide into the Wire Guide Locking Device.

If Using Proximal Wire Port (PWP) and Pre-Positioned Long Wire Guide:

NOTE: For best results, the wire guide should be kept wet.

1. Remove stylet.
2. Advance the tip of the balloon catheter over the pre-positioned long wire guide, ensuring the wire guide exits the catheter at the PWP.
3. Advance the deflated balloon in short increments through the accessory channel, using standard long wire guide exchange technique, until endoscopically visualized exiting the scope.
4. Follow steps 5-8 in previous section.
5. To remove the deflated balloon from the accessory channel, straighten the scope tip and open the elevator. Maintain negative pressure while withdrawing the balloon, using a continuous twisting motion, via standard long wire guide exchange technique.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.