EN
3

Flow-Directed Balloon Catheter

Instructions for Use

DA
3

Flow-styret ballonkateter

Brugsanvisning

DE
4

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EL
5

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T _ F D B _ R E V 2

FLOW-DIRECTED BALLOON CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Flow-Directed Balloon Catheter is a double-lumen 5.3 French catheter with a balloon tip. The balloon has an inflation capacity of 1 ml and a maximum inflation diameter of 10 mm.

INTENDED USE

The Flow-Directed Balloon Catheter is intended for venous sampling and pressure monitoring. The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

None known

WARNINGS

- **This product contains NATURAL RUBBER LATEX, which may cause allergic reactions.**
- Use only bacteria-filtered carbon dioxide for balloon inflation.
- Do not use air for balloon inflation; doing so could produce air emboli in the event of balloon rupture.

PRECAUTIONS

- Manipulation of the product requires fluoroscopic control.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Using standard percutaneous technique, prepare and introduce an introducer sheath into the venous access site.
2. Introduce the Flow-Directed Balloon Catheter through the sheath into the vessel, and advance until the balloon tip has passed completely through the sheath.
3. Inflate the balloon for flow-assisted passage into the target vessel location.
4. Once positioned at the desired vessel location, the catheter may be used for monitoring and recording pressures in the target vessel, acquiring blood samples or delivering fluids.

NOTE: To allow free movement of the catheter within the vessel, ensure that the inflated balloon's diameter is less than the vessel diameter.

5. If required, a .035 inch or smaller wire guide can be used for catheter exchange.

NOTE: Deflate the balloon prior to retraction or removal.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2013
Bloomington, IN 47404 U.S.A.