Fascial Dilator Set

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Fasciedilatatorsæt

Fasziendilatator-Set

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T _ F D _ R E V 1

FASCIAL DILATOR SET

Used for fascial dilation. Intended for one-time use.

SET CONTAINS

- Radiopaque dilators 37 cm long
- .038 inch (0.97 mm) diameter stainless steel wire guide 80 cm long with 3 cm flexible tip

NOTE: Prior to use, immerse AQ® dilators in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

CONTRAINDICATIONS

There are no known contraindications.

SUGGESTED INSTRUCTIONS FOR USING FASCIAL DILATOR SET

1. Advance the wire guide directly from the enclosed holder.
2. Introduce the flexible tip of the wire guide into an 18 gage needle which has previously been positioned into the renal collecting system under fluoroscopic control. (The safety index mark, 39 cm from the wire guide's flexible tip, will be externally visible under normal circumstances.)
3. Once the wire guide is in position, remove the puncture needle from the wire guide prior to introduction of the dilators over the wire guide.
4. Progressing from the smallest to the largest appropriate size, pass the radiopaque dilators over the wire guide while maintaining the wire guide's position.

CAUTION: Do not advance dilators beyond the 39 cm safety index mark. When the end of the dilator is at the mark, a 2 cm length of wire protrudes from the tip. Further advancement of the dilator will override the wire guide.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd