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
Fig. 19


Fig. 20

GIANTURCO-ROEHM BIRD'S NEST® VENA CAVA FILTER FOR FEMORAL PLACEMENT


CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Gianturco-Roehm Bird's Nest Vena Cava Filter set consists of:

- 12 French introducer sheath assembly, appropriate length
- 11 French filter catheter, 45 cm long, with preloaded stainless steel filter and wire guide pusher

The stainless steel filter is introduced and placed in the inferior vena cava utilizing standard percutaneous entry (Seldinger) technique.


1. Safety Lock (Red)
2. Mark
3. Tuohy-Borst Valve
4. Wire Guide Pusher
5. Stopcock

Cross Section of Filter Catheter with Preloaded Filter

1. Filter Wires
2. Hook Wire Struts
3. Junction Point of Hook Wires
4. Hooks
5. Hook Stops
6. Wire Guide Pusher
7. Wire Guide Pusher Attachment Site

Bird's Nest® Filter

1. Filter Wires
2. Hook Wire Struts
3. Junction Point of Hook Wires
4. Hooks
5. Hook Stops
6. Wire Guide Pusher (not shown)
7. Wire Guide Pusher Attachment Site

INTENDED USE

The Gianturco-Roehm Bird's Nest Vena Cava Filter is intended for the prevention of recurrent pulmonary embolism via placement in the vena cava in the following situations:

1. Pulmonary thromboembolism when anticoagulants are contraindicated;
2. Failure of anticoagulant therapy in thromboembolic diseases;
3. Emergency treatment following massive pulmonary embolism where anticipated benefits of conventional therapy are reduced;
4. Prophylactically in patients with chronic, recurrent pulmonary embolism where anticoagulant therapy has failed or is contraindicated.

CONTRAINDICATIONS

- Vena cava filters should not be implanted in patients with risk of septic embolism.
- The Bird's Nest Filter should not be used in patients with a vena cava larger than 40 mm diameter.

WARNINGS

- Be certain the junction point of distal pair of hooks is in catheter tip when positioning these hooks into caval wall. If there is uncertainty regarding hook fixation, the catheter may be advanced up over the hook wires to the edge of the hooks. The catheter may then be repositioned, and the filter placement re-initiated in a slightly different caval orientation or level.
- No technique will completely eliminate the possibility of recurrent PTE. (With the Bird's Nest Vena Cava Filter, the observed incidence of recurrent PTE is clinically acceptable.)
- The Bird's Nest Vena Cava Filter is non-thrombogenic, but may occlude if it traps a large volume of embolic material in a short time period. Thrombolytic agents administered via femoral venipunctures have been successful in re-establishing filter and caval patency.
- If filter migration occurs, transcatheter retrieval of the filter is not recommended.
- Retroperitoneal hematomas or hemorrhage may occur when filter is used in patients receiving thrombolytic therapy.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.
- **Insertion requires the use of fluoroscopy.** Only those physicians skilled in the use of fluoroscopic images should attempt filter insertion. Read and follow the enclosed instructions carefully.
- **Do not remove the preloaded filter from the catheter. Attempts to reload may damage the catheter or the filter.**
- An effort should be made prior to filter placement to determine if emboli are coming from the lower extremities. The filter is not effective in filtering emboli from a superior or cardiac source if placed in the inferior vena cava.

MAGNETIC RESONANCE IMAGING

MRI-Safe: No additional risk to the patient, but may affect the quality of the diagnostic information if scanning in the area of the filter. MRI diagnostic procedures should be postponed for 6 weeks following filter implantation to assure device incorporation into the vessel wall.

POTENTIAL ADVERSE EVENTS

- Damage to the vena cava
- Pulmonary embolism
- Filter embolization
- Vena cava perforation
- Vena cava occlusion or thrombosis
- Hemorrhage
- Hematoma at vascular access site
- Infection at vascular access site
- Death

INSTRUCTIONS FOR USE

Illustrations

1. Prep and drape the access site.
2. After infiltration of local anesthesia, make a 3-4 mm skin nick at the access site, and perform a standard venipuncture using an 18 gage access needle. (**Fig. 1**)
3. Introduce a .038 inch (0.97 mm) diameter wire guide and advance it into the inferior vena cava. (**Fig. 2**)
4. Leaving the wire guide in place, withdraw the needle. (**Fig. 3**)
5. Introduce a standard diagnostic catheter and perform an inferior vena cavogram to localize the renal veins, typically located at the L1-L2 intervertebral space. Once localized, a lead marker taped to the skin at that level will help to serve as a radiographic landmark to assist precise filter placement.
6. Over the wire guide, exchange the diagnostic catheter for the filter introducer sheath assembly. Advance the introducer sheath assembly over the wire guide to its hub. (**Fig. 4**)

CAUTION: Do not advance distal end of introducer sheath assembly beyond distal tip of wire guide.

7. Prepare the filter catheter for introduction by priming its lumen with contrast medium via its sidearm. Use of a contrast-filled syringe and flexible connecting tube attached to the stopcock on the filter catheter sidearm is recommended.

8. Remove the introducer and wire guide from the sheath. (**Fig. 5**)
9. Insert the filter catheter into the introducer sheath and Luer lock it into position. The tip of the filter catheter will then extend approximately 1.5 cm beyond the tip of the sheath. (**Fig. 6**) While performing hand injections of contrast medium through the filter catheter sidearm, the filter catheter/introducer sheath assembly can be properly positioned to initiate filter placement. **When utilizing a femoral approach, position the filter catheter tip just below (caudad to) the renal veins. (Fig. 7)**
10. Loosen the Tuohy-Borst valve (**Fig. 8**) and while holding the wire guide pusher stationary, withdraw the filter catheter/introducer sheath assembly to the mark on the wire guide pusher. (**Fig. 9**) This movement will expose the distal anchoring hooks. Fluoroscopically, this corresponds to the position of the junction point of the hook wires remaining just within the tip of the filter catheter. (**Fig. 10**) If this can be clearly visualized with fluoroscopy, the mark can be ignored.
NOTE: At the end of the filter, the V-shaped hook wire struts are readily visible, but filter wires are too fine to be identified fluoroscopically.
11. **Gently** advance the filter catheter/introducer sheath assembly **one time**, 1-3 mm, to secure the hooks to the vena cava wall. (**Fig. 11**)
Be certain junction point of distal pair of hooks is in catheter tip when positioning these hooks into caval wall. If there is uncertainty regarding hook fixation, the catheter may be advanced up over the hook wires to the edge of the hooks. The catheter may then be repositioned, and the filter placement re-initiated in a slightly different caval orientation or level.
CAUTION: Overly forceful or multiple jabs with the filter catheter/introducer sheath assembly could result in perforation of the wall of the inferior vena cava by the exposed hooks and struts of the filter.
12. Again, hold the wire guide pusher stationary and withdraw the filter catheter/introducer sheath assembly 1-3 cm. (**Fig. 12**) This will facilitate subsequent passage of the filter wires through the tip of the filter catheter and provide room for filter formation within the vena cava.
CAUTION: Do not withdraw the filter catheter/introducer sheath assembly into the iliac vein. Repositioning into the inferior vena cava may prove difficult.
13. Holding the filter catheter/introducer sheath assembly stationary and with fluoroscopic control, advance the wire guide pusher with a steady, smooth motion until the junction point of the proximal hook wires is seen fluoroscopically at the filter catheter tip. (**Fig. 13**)
CAUTION: Do not advance too quickly or forcefully as wire kinking within filter catheter may occur, making advancement difficult.
14. Advance the filter catheter/introducer sheath assembly so that the proximal junction point is in close proximity to or overlaps the junction point of previously positioned hook wires. (**Fig. 14**)
NOTE: The approximation or slight overlap of the junction points ensures adequate compression and formation of the filter wires and places the second pair of hooks in the inferior vena cava. Placement of a hook in a renal vein or passage of filter wires distal to the first hook position (wire prolapse) may occasionally occur and is of no consequence.
15. While maintaining slight forward pressure on the wire guide pusher, slowly withdraw the filter catheter/introducer sheath assembly. (**Fig. 15**) This will permit the proximal hook wires to exit the catheter, springing open to engage with the caval wall. This is readily seen fluoroscopically.
To be sure the hooks are completely free from the filter catheter tip, the filter catheter/introducer sheath assembly should be withdrawn to within approximately 1 cm of the Bird's Nest Filter handle attached to the wire guide pusher.
CAUTION: Do not advance the wire guide pusher after the proximal hook wires have exited the filter catheter tip. Forceful advancement of wire guide pusher may alter the orientation of the filter attachment site, making disengagement difficult.
16. Gently tug or use a slight to-and-fro motion on the wire guide pusher to secure the second pair of hooks to the vena cava wall before disengaging the filter. (**Fig. 16**)
17. To release filter:
 - a. Maintaining wire guide pusher position, remove red safety lock from handle by grasping both sides and pulling up and away from the handle. (**Fig. 17**)
 - b. Holding the wire guide pusher steady, use the handle to fully depress the thumb grip of the wire guide pusher forward into the cannula. Maintaining this position will release the filter. (**Fig. 18**) Separation can be seen fluoroscopically.
 - c. Release pressure on the thumb grip of the handle. (**Fig. 19**)
 - d. Pull the wire guide pusher's distal end back into the filter catheter. (**Fig. 20**)
18. Remove the wire guide pusher and filter catheter, leaving the introducing sheath in place. Reposition the sheath for a final vena cavogram.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.