- EN**
3 **Firlit-Kluge Urethral Stent**
Instructions for Use
- CS**
4 **Uretrální stent Firlit-Kluge**
Návod k použití
- DA**
5 **Firlit-Kluge urethrastent**
Brugsanvisning
- DE**
6 **Firlit-Kluge Urethrastent**
Gebrauchsanweisung
- EL**
7 **Ουρηθρική ενδοπρόσθεση Firlit-Kluge**
Οδηγίες χρήσης
- ES**
8 **Stent uretral Firlit-Kluge**
Instrucciones de uso
- FR**
9 **Endoprothèse urétrale de Firlit-Kluge**
Mode d'emploi
- HU**
10 **Firlit-Kluge urethrális sztent**
Használati utasítás
- IT**
11 **Stent uretrale Firlit-Kluge**
Istruzioni per l'uso
- NL**
12 **Firlit-Kluge urethrastent**
Gebruiksaanwijzing
- NO**
13 **Firlit-Kluge uretrastent**
Bruksanvisning
- PL**
14 **Stent Firlit-Kluge do cewki moczowej**
Instrukcja użycia
- PT**
15 **Stent uretral Firlit-Kluge**
Instruções de utilização
- SV**
16 **Firlit-Kluge uretral stent**
Bruksanvisning

Firlit-Kluge Urethral Stent

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

This device is constructed out of silicone tubing and a silicone ball. The stent tubing comes in both 8 French and 10 French. The ball is located between the distal and proximal ends of the stent tubing. The 8 French stent comes in lengths of 31, 40, and 50 cm. The 10 French stent comes in lengths of 31 and 50 cm. The one variation of the 8 French stent is 31 cm in length and has sideports.

INTENDED USE

This device is used for stenting the urethra during hypospadias or epispadias repair and to allow postoperative drainage of the bladder in infants (29 days to less than 2 years old) and children (2 years old to less than 12 years old).

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions that create unacceptable risk during catheterization.

WARNINGS

None known

PRECAUTIONS

This device is intended for use by physicians trained and experienced in urological techniques. Standard techniques should be employed.

POTENTIAL ADVERSE EVENTS

- Urinary tract infection
- Perforation or laceration of the urethra
- Perforation or laceration of the bladder
- Abrasions

INSTRUCTIONS FOR USE

1. Measure the urethra and cut the long segment of the stent to the correct length.
2. At the appropriate time during the procedure, place the stent through the urethra to the bladder.
3. Position the molded silicone ball at the glans meatus.
4. Using absorbable suture, anchor the silicone ball to the glans meatus.
5. After the procedure, apply a dressing and, if necessary, trim the distal excess portion of the stent beyond the silicone ball.

NOTE: The stent must not remain indwelling for more than 21 days. During the healing process, the sutures will be absorbed, and the stent should be expelled spontaneously during micturition.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd