Formula 418 Stent

**Balloon-Expandable Stent Delivery System for Use
in Biliary Indication**

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**FR
11**

T _ F O R M B _ R E V 3

FORMULA 418 STENT BALLOON-EXPANDABLE STENT DELIVERY SYSTEM FOR USE IN BILIARY INDICATION

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Formula 418 Stent is a balloon-expandable stent made of 316L stainless steel with a slotted tube configuration. It is pre-mounted on a balloon catheter, which serves as the delivery system. The stent is positioned between two radiopaque marker bands, which are located inside the balloon at the proximal and distal tapers of the balloon. The cannula design of the Formula 418 Stent provides a low outside diameter profile and is pre-mounted on either an 80 cm or 135 cm balloon catheter delivery system. The stent is available in nominal expanded diameters ranging from 3-8 mm with lengths ranging from 12-30 mm.

INTENDED USE

The Formula 418 Stent is intended for use in palliation of malignant neoplasms in the biliary tree. The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of biliary duct access sheaths, guiding catheters/introducers and wire guides should be employed.

CONTRAINDICATIONS

The Formula 418 Stent is contraindicated for use in:

- Total biliary obstructions and postoperative strictures that cannot be dilated to permit passage of the sheath introducer.
- Stenting of a perforated duct where leakage from the duct could be exacerbated by placement of a stent.
- Patients with bleeding disorders.
- Severe ascites.

WARNINGS

- Persons allergic to 316L stainless steel may suffer an allergic reaction to this implant.
- The safety and effectiveness of this device for use in the vascular system have not been established.
- Stents bridging the papilla should extend no more than 5 mm beyond the papilla and into the duodenum after deployment.

PRECAUTIONS

Stent Placement

- The stent is preloaded on its delivery balloon and should not be removed. Do not attempt to place stent on another balloon catheter for deployment.
- DO NOT wipe or clean stent or catheter with organic solvent (i.e., isopropyl alcohol).
- DO NOT expose catheter to temperatures above 130° F.
- Special care must be taken not to handle or in any way disrupt the stent on the balloon. This is particularly important during removal of the catheter from packaging, placement over the wire guide and advancement through the large-bore Tuohy-Borst "Y" adapter and guiding catheter hub.
- Do not pre-inflate balloon prior to stent deployment. (See instructions for Preparation of Balloon Catheter)
- Use only the appropriate balloon inflation media. The standard inflation medium is a 1:1 mixture of contrast medium and normal saline. Do not use air or any gaseous substance as a balloon inflation medium.
- Expansion of the stent should not be undertaken if the stent is not properly positioned in the biliary duct. If the position of the stent is not optimal, it should not be expanded. (See instructions for Removal of Unexpanded Stent)

- Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure as indicated on the product label. Use of pressures higher than specified on the product label may result in a ruptured balloon, with possible biliary duct damage and dissection.
- Use prior to the expiration date on the package label.
- Intended for one-time use. (Do not re-sterilize)

MAGNETIC RESONANCE IMAGING

Non-clinical testing has demonstrated that the Formula 418 stent is MR Conditional. A patient with this stent may undergo MRI immediately after placement under the following conditions:

- Static magnetic field of 3 Tesla or less
- Maximum spatial gradient magnetic field of 720 Gauss/cm
- Maximum MR system reported whole body averaged specific absorption rate (SAR) of 3.0 W/kg for 15 minutes of scanning

In non-clinical testing, the Formula 418 produced a temperature rise of less than 1 degree C at an MR system reported maximum whole body averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of MR scanning in a 3-Tesla, Excite, Software G3.0-052B, General Electric Medical Systems MR scanner.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Formula 418. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

The effect of heating in the MRI environment for overlapping stents or stents with fractured struts is unknown.

POTENTIAL ADVERSE EVENTS

Adverse events that may be associated with the use of a biliary stent include:

- Death
- Stent migration
- Hemorrhage/hematoma
- Bile duct dissection, rupture and perforation
- Sepsis/infection/abscess
- Total occlusion/obstruction
- Peritonitis
- Cholangitis
- Stent embolization
- Tumor overgrowth at the stent ends
- Parenchymal hemorrhage
- Pancreatitis
- Hemobilia
- Stent obstruction secondary to tumor growth through the stent
- Stent restenosis of the stented bile duct
- Rupture into a neighboring organ or the retroperitoneum

PRODUCT RECOMMENDATIONS

The Formula 418 Stent is used in conjunction with equipment required for a conventional biliary stent procedure, including, but not limited to, a biliary access set, sheath introducer, guiding catheter, wire guide and inflation device.

Predilations of the duct must take into consideration proximal strictures which may hinder advancement of the stent in addition to strictures beyond the placement area which may prevent advancement of the stent across the primary stricture.

Wire Guide Use and Selection

The Formula 418 Stent delivery system is compatible with a .018 inch wire guide.

Guiding Catheter/Introducer Selection

Correct guiding catheter/introducer selection and technique are necessary for use of the stent. Ensure that the inside lumen of the guiding catheter/introducer is of sufficient size to allow unobstructed passage of the stent/balloon delivery catheter.

Stent Size Selection

The stent selected should have an expanded diameter approximately equal to, or slightly larger than, the estimated reference diameter of the bile duct (proximal and distal to the stricture) to be stented. The length of the stent should be chosen to adequately cover the length of the stricture.

INSTRUCTIONS FOR USE

Preparation of Balloon Catheter

1. Remove the stent/balloon catheter from the package, remove protective sleeve from distal tip of catheter and inspect the stent to ensure it has not been damaged.
2. Prepare balloon lumen with a standard 1:1 contrast-saline mixture as follows:
WARNING: Do not attempt pre-inflation technique to purge balloon lumen. Do not use air or any gaseous medium to inflate the balloon.
 - a. Using a 20 ml syringe containing 5 ml of contrast-saline mixture, apply negative pressure for 20-30 seconds.
 - b. Release pressure, allowing negative pressure to draw mixture into balloon lumen.
 - c. Detach syringe, leaving a meniscus of mixture on the hub of the balloon lumen.
 - d. Prepare inflation device in standard manner and purge to remove all air from syringe and tubing.
 - e. Attach inflation device to balloon lumen directly, ensuring no air bubbles remain at connection.
 - f. Pull negative pressure on inflation device.

CAUTION: Significant amounts of air in the balloon may cause uneven expansion of the stent and difficulty in stent deployment.

3. Moisten the stent with heparinized saline.
CAUTION: Do not use gauze sponges, as fibers may disrupt stent.
4. Flush the wire guide lumen of the balloon catheter in standard fashion to purge air.
5. Advance a .018 inch wire guide of appropriate length across target stricture.
6. Insert the appropriate guiding catheter/introducer.
7. Advance the premounted stent/balloon catheter over the wire into either the introducer valve or Tuohy-Borst "Y" adapter.
 - a. If using an introducer with a valve, make sure the flared end of the insertion tool (provided in the package) is loaded over the premounted stent on the balloon catheter. Pass the insertion tool loaded with the premounted stent through the introducer valve. Push the stent/balloon catheter into the body of the introducer. Slide the insertion tool proximally up the catheter shaft away from the guiding catheter/introducer. A slight contact of the stent with the introducer may be felt, but there must be no resistance.
 - b. If using a Tuohy-Borst "Y" adapter, advance the premounted stent/balloon catheter over the wire and into the fully opened large-bore Tuohy-Borst "Y" adapter. Gently advance the stent/balloon catheter completely through the Tuohy-Borst "Y" adapter and into the guiding catheter/introducer. A slight contact of the stent with the guiding catheter/introducer may be felt, but there must be no resistance.

CAUTION: If resistance is encountered, do not force passage. Resistance may indicate damage to stent.

Positioning the Stent

1. Ensure guiding catheter/introducer stability before advancing the stent delivery balloon into the biliary duct.
CAUTION: If initial guiding catheter/introducer position is lost, avoid pulling or pushing the guiding catheter/introducer over the stent. If this is done, the distal end of the guiding catheter/introducer may damage the stent.
2. Position the stent across the stricture, using both the distal and proximal balloon markers as reference points. Optimal placement requires the ends of the stent to extend beyond the margin of the duct segment to be stented.
CAUTION: If the stent/balloon catheter does not readily advance through the stricture, do not force. If the stent will not advance in spite of good guiding catheter/introducer support, consider dilating the stricture, or changing the wire guide or the guiding catheter/introducer. (Refer to instructions for Removal of Unexpanded Stent)

Balloon Expansion/Stent Deployment

1. Prior to stent expansion, use high-resolution fluoroscopy to verify that the stent has not been dislodged during positioning.
2. To expand the stent, inflate the balloon to the recommended expansion pressure indicated on product label. Do not move stent system during deployment.
3. Complete expansion and apposition of the stent against the bile duct wall is necessary for clinical success. Do not exceed rated burst pressure of the balloon as indicated on the product label.
4. Once the stent has been deployed, post-deployment inflation is at the discretion of the operator to achieve optimum fluoroscopic appearance. The inflated diameter of the post-dilatation balloon should not exceed the reference duct diameter by more than 10%.

Balloon Deflation and Removal

1. Completely deflate the balloon by pulling negative pressure on the inflation device or a 20 ml syringe.
CAUTION: Allow enough time for the balloon to fully deflate prior to removal.
2. Slowly withdraw the balloon catheter from the stent while maintaining negative pressure on the balloon. Maintain position of the guiding catheter/introducer to prevent stent from being drawn out of intended position. Observe under fluoroscopy to ensure that the balloon disengages from the stent.

Removal of Unexpanded Stent

Do not attempt to pull an unexpanded stent back into the guiding catheter/introducer. The stent/balloon catheter (stent delivery system) should be withdrawn until the proximal end of the stent is aligned with the distal tip of the guiding catheter/introducer. Withdraw the guiding catheter/introducer and stent delivery system as a single unit, leaving the wire guide in place.

WARNING: If stent is removed, do not attempt to reuse the device. Damage to the stent may occur upon removal.

Stent-Assisted, High-Pressure Balloon Dilatation

A deployed stent may be further expanded using a low-profile, high-pressure balloon catheter. The stent should have been well embedded into the bile duct wall by the deployment balloon before the stricture is recrossed. If recrossing is necessary, the position of the stent should be monitored and care should be taken to avoid displacing the stent with subsequent manipulation of wire guides and balloons. An appropriate balloon catheter should be used with a very flexible wire guide.

If a Formula 418 Stent is dilated with a high-pressure balloon, the balloon should be appropriately sized for the duct diameter. The inflated diameter of the post-dilatation balloon should not exceed the reference duct diameter by more than 10%. Use of a new balloon for post-deployment dilatation is suggested. Avoid balloon catheters with stiff tips that might displace the stent during advancement. Avoid balloons that “wing out” and/or re-wrap inconsistently, which might displace the stent during withdrawal of the balloon from the stent.

WARNING: If performing post-dilatation of the stent, use a balloon catheter indicated for biliary stent deployment.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.