**EN** **Formula 418®** Vascular Balloon-Expandable Stent System

3 Instructions for Use

**DA** **Formula 418®** vaskulært system med ballon-eksponderbar stent

9 Brugsanvisning

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**SV** **Formula 418®** vaskulärt ballongexpanderbart stentsystem

64 Bruksanvisning

## FORMULA 418® VASCULAR BALLOON-EXPANDABLE STENT SYSTEM

### DEVICE DESCRIPTION

The Formula 418 Vascular Stent is a balloon-expandable stent made of 316L stainless steel with a slotted tube configuration. It is pre-mounted on an over-the-wire balloon catheter, which serves as the delivery system (Formula 418 Vascular Balloon-Expandable Stent System). The stent is positioned between two radiopaque marker bands, which are located inside the balloon at the proximal and distal tapers of the balloon. The cannula design of the stent provides a low outside diameter profile, which permits use with a 5.0 French sheath or 6.0 French guiding catheter for 3-6 mm diameter sizes and a 6.0 French sheath and 7 French guiding catheter for 7-8 mm diameter sizes. The Formula 418 Vascular Balloon-Expandable Stent System is pre-mounted on 80 cm and 135 cm over-the-wire balloon catheter delivery systems. The stent is available in nominal expanded diameters of 3, 4, 5, 6, 7, and 8 mm with lengths of 12, 16, 20, 24 and 30 mm each. (Table 1)

**Table 1: Device Specifications (diameters and lengths)**

Order Number	Catheter Length (cm)	Stent Diameter (mm)	Stent Length (mm)	Minimum Guiding Catheter/Sheath* Inside Diameter (inches)	Guiding Catheter/Sheath (French size)
FOV418-18-80-3-12	80	3	12	0.068	6.0/5.0
FOV418-18-80-3-16	80	3	16	0.068	6.0/5.0
FOV418-18-80-3-20	80	3	20	0.068	6.0/5.0
FOV418-18-80-3-24	80	3	24	0.068	6.0/5.0
FOV418-18-80-3-30	80	3	30	0.068	6.0/5.0
FOV418-18-80-4-12	80	4	12	0.068	6.0/5.0
FOV418-18-80-4-16	80	4	16	0.068	6.0/5.0
FOV418-18-80-4-20	80	4	20	0.068	6.0/5.0
FOV418-18-80-4-24	80	4	24	0.068	6.0/5.0
FOV418-18-80-4-30	80	4	30	0.068	6.0/5.0
FOV418-18-80-5-12	80	5	12	0.068	6.0/5.0
FOV418-18-80-5-16	80	5	16	0.068	6.0/5.0
FOV418-18-80-5-20	80	5	20	0.068	6.0/5.0
FOV418-18-80-5-24	80	5	24	0.068	6.0/5.0
FOV418-18-80-5-30	80	5	30	0.068	6.0/5.0
FOV418-18-80-6-12	80	6	12	0.068	6.0/5.0
FOV418-18-80-6-16	80	6	16	0.068	6.0/5.0
FOV418-18-80-6-20	80	6	20	0.068	6.0/5.0
FOV418-18-80-6-24	80	6	24	0.068	6.0/5.0
FOV418-18-80-6-30	80	6	30	0.068	6.0/5.0
FOV418-18-80-7-12	80	7	12	0.078	7.0/6.0
FOV418-18-80-7-16	80	7	16	0.078	7.0/6.0
FOV418-18-80-7-20	80	7	20	0.078	7.0/6.0
FOV418-18-80-7-24	80	7	24	0.078	7.0/6.0
FOV418-18-80-7-30	80	7	30	0.078	7.0/6.0
FOV418-18-80-8-12	80	8	12	0.078	7.0/6.0
FOV418-18-80-8-16	80	8	16	0.078	7.0/6.0
FOV418-18-80-8-20	80	8	20	0.078	7.0/6.0
FOV418-18-80-8-24	80	8	24	0.078	7.0/6.0
FOV418-18-80-8-30	80	8	30	0.078	7.0/6.0
FOV418-18-135-3-12	135	3	12	0.068	6.0/5.0
FOV418-18-135-3-16	135	3	16	0.068	6.0/5.0

Order Number	Catheter Length (cm)	Stent Diameter (mm)	Stent Length (mm)	Minimum Guiding Catheter/Sheath* Inside Diameter (inches)	Guiding Catheter/Sheath (French size)
FOV418-18-135-3-20	135	3	20	0.068	6.0/5.0
FOV418-18-135-3-24	135	3	24	0.068	6.0/5.0
FOV418-18-135-3-30	135	3	30	0.068	6.0/5.0
FOV418-18-135-4-12	135	4	12	0.068	6.0/5.0
FOV418-18-135-4-16	135	4	16	0.068	6.0/5.0
FOV418-18-135-4-20	135	4	20	0.068	6.0/5.0
FOV418-18-135-4-24	135	4	24	0.068	6.0/5.0
FOV418-18-135-4-30	135	4	30	0.068	6.0/5.0
FOV418-18-135-5-12	135	5	12	0.068	6.0/5.0
FOV418-18-135-5-16	135	5	16	0.068	6.0/5.0
FOV418-18-135-5-20	135	5	20	0.068	6.0/5.0
FOV418-18-135-5-24	135	5	24	0.068	6.0/5.0
FOV418-18-135-5-30	135	5	30	0.068	6.0/5.0
FOV418-18-135-6-12	135	6	12	0.068	6.0/5.0
FOV418-18-135-6-16	135	6	16	0.068	6.0/5.0
FOV418-18-135-6-20	135	6	20	0.068	6.0/5.0
FOV418-18-135-6-24	135	6	24	0.068	6.0/5.0
FOV418-18-135-6-30	135	6	30	0.068	6.0/5.0
FOV418-18-135-7-12	135	7	12	0.078	7.0/6.0
FOV418-18-135-7-16	135	7	16	0.078	7.0/6.0
FOV418-18-135-7-20	135	7	20	0.078	7.0/6.0
FOV418-18-135-7-24	135	7	24	0.078	7.0/6.0
FOV418-18-135-7-30	135	7	30	0.078	7.0/6.0
FOV418-18-135-8-12	135	8	12	0.078	7.0/6.0
FOV418-18-135-8-16	135	8	16	0.078	7.0/6.0
FOV418-18-135-8-20	135	8	20	0.078	7.0/6.0
FOV418-18-135-8-24	135	8	24	0.078	7.0/6.0
FOV418-18-135-8-30	135	8	30	0.078	7.0/6.0

\*See manufacturer specifications for French (Fr) equivalent.

## INTENDED USE

The Formula 418 Vascular Balloon-Expandable Stent System is intended for stenting of peripheral arteries.

## CONTRAINDICATIONS

The Formula 418 Vascular Balloon-Expandable Stent System is contraindicated for use in:

- Patients for whom antiplatelet and/or anticoagulant therapy is contraindicated.
- Patients who have a lesion that cannot be crossed with a wire or a balloon angioplasty catheter.
- Patients who have stenoses that cannot be dilated to permit passage of the stent.
- Patients with bleeding disorders.
- Stenting of an arterial vessel where leakage from the artery could be exacerbated by placement of a stent.
- Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.

## WARNINGS

- The use of this device carries the associated risks of subacute thrombosis, vascular complications and/or bleeding events. Judicious selection of patients is necessary.
- Persons allergic to 316L stainless steel may suffer an allergic reaction to this implant.
- The stent should not be placed in a location with significant external compression which could crush the stent.
- The stent should not be deployed in areas subject to significant flexure or bending to reduce the risk of stent deformations or fracture.

## PRECAUTIONS

### General

- Implantation of the stent should be performed only by physicians who have received adequate training. Adequate training consists of specific training in vascular stenting techniques and procedures under the supervision of an interventional physician experienced in vascular stent implantation.
- Stent placement should only be performed at hospitals where emergency artery bypass graft surgery can be readily performed.
- Subsequent restenosis may require redilatation of the arterial segment containing the stent. The long-term outcome following repeat dilatation of an endothelialized Formula 418 Vascular Stent is unknown at present.

### Stent Placement

- Manipulation of the Formula 418 Vascular Balloon-Expandable Stent System requires fluoroscopic control.
- The stent is preloaded on its delivery balloon and should not be removed. Do not attempt to place stent on another balloon catheter for deployment.
- Deploy the stent over a wire guide not exceeding 0.018 inch diameter.
- DO NOT wipe or clean stent or catheter with organic solvent (i.e. isopropyl alcohol).
- DO NOT expose catheter to temperatures above 54°C (130°F).
- DO NOT use power injection systems with the delivery system.
- DO NOT rotate any part of the system during deployment.
- Special care must be taken not to handle or in any way disrupt the stent on the balloon. This is particularly important during removal of the catheter from packaging, placement over the wire guide and advancement through the large-bore Tuohy-Borst "Y" adapter and guiding catheter hub.
- Do not pre-inflate balloon prior to stent deployment (See Preparation of Balloon Catheter section in these Instructions for Use).
- Use only the appropriate balloon inflation medium. The standard inflation medium is a 1:1 mixture of contrast medium and normal saline. Do not use air or any gaseous substance as a balloon inflation medium.
- Expansion of the stent should not be undertaken if the stent is not properly positioned in the vessel (See Removal of Unexpanded Stent section in these Instructions for Use).
- Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure as indicated on the product label. Use of pressures higher than specified on the product label may result in a ruptured balloon, with possible intimal damage and dissection.
- Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stent and may cause acute closure of the vessel, requiring additional intervention (further dilatation or placement of additional stents).
- Use prior to the expiration date on the package label.
- Intended for one-time use (Do not re-sterilize).
- The use of mechanical atherectomy devices within an implanted stent may cause damage to the stent.
- Repositioning of the device after deployment is not possible.

### MAGNETIC RESONANCE IMAGING

Non-clinical testing has demonstrated that the Formula 418 Vascular Stent is **MR Conditional** according to ASTM F2503, *Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment*. A patient with this stent can be scanned safely anytime after placement under the following conditions.

- Static magnetic field of 3-Tesla or 1.5-Tesla
- Maximum spatial magnetic gradient of 720 Gauss/cm or less
- Product of the spatial gradient and the static magnetic field of 21.6 T<sup>2</sup>/m or less
- MR system reported whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg or less for 15 minutes of scanning

Non-clinical evaluation was conducted in an MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual). In non-clinical testing, a single and two overlapped Formula 7 mm x 20 mm stents produced a maximum temperature rise of 2.6 °C and 3.3 °C, respectively, during 15 minutes of MRI (i.e., for one scanning sequence) performed in a 3 Tesla MR system (3 Tesla/128 MHz, Excite, Software G3.0-052B, General Electric Healthcare) at an MR system reported whole-body-averaged SAR of 3.0 W/kg (associated with a calorimetry measured whole-body-averaged value of 2.8 W/kg).

The effect of heating in the MRI environment for stents with fractured struts is unknown.

## Image Artifact

MR image quality may be compromised if the area of interest is within the lumen or within approximately 15 mm of the Formula 7 mm x 20 mm stent. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic stent.

## MedicAlert Foundation

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted at:

- Mail: MedicAlert Foundation International  
2323 Colorado Avenue  
Turlock, CA 95382
- Phone: +1-888-633-4298 (toll free)  
+1-209-668-3333 (from outside the US)
- Fax: +1-209-669-2450
- Web: [www.medicalert.org](http://www.medicalert.org)

## POTENTIAL ADVERSE EVENTS

Adverse events that may be associated with the use of a vascular stent are presented in the following alphabetical list and include, but are not limited to:

- Abscess
- Allergic reaction to stainless steel or contrast agents
- Arrhythmias (ventricular fibrillation, ventricular tachycardia, other)
- Arteriovenous fistula
- Bowel infarct
- Death
- Dialysis
- Dissection
- Drug reaction to antiplatelet agents
- Drug reaction, allergic reaction to contrast media
- Emboli (air, tissue, or thrombotic emboli) resulting in tissue ischemia/infarction
- Emergency surgery to correct vascular complications
- Emergent renal artery bypass surgery
- Extremity ischemia/amputation
- Fever
- Gastrointestinal symptoms from anticoagulation/antiplatelet medication
- Hematoma at vascular access site
- Hemorrhage requiring transfusion
- Hypersensitivity reactions
- Hypotension/hypertension
- Infection and pain at vascular access site
- Intimal tear
- Kidney infarct
- Myocardial infarction
- Myocardial ischemia
- Nephrectomy
- Peripheral neuropathy
- Pseudoaneurysm at vascular access site
- Pseudoaneurysm formation
- Renal artery thrombosis, aneurysm, rupture, perforation, occlusion, spasm, or restenosis
- Renal insufficiency or failure
- Stent migration or embolization
- Stent misplacement
- Stroke/cerebral vascular accident
- Tissue necrosis or ulceration

## PRODUCT RECOMMENDATIONS

### Wire Guide Use and Selection

The Formula 418 Vascular Balloon-Expandable Stent System is compatible with 0.018 inch wire guides.

### Guiding Catheter/Introducer Selection

Correct guiding catheter/sheath selection and technique are necessary for use of the stent. Ensure that the inside lumen of the guiding catheter/sheath is of sufficient size to allow unobstructed passage of the Formula 418 Vascular Balloon-Expandable Stent System.

**NOTE:** Minimum inside diameter requirements for the guiding catheter/sheath are shown in **Table 1**.

Standard vascular catheter curves must be selected to provide adequate guiding catheter/sheath "back-up support" to achieve successful stent placement.

### Stent Size Selection

The stent selected should have an expanded diameter approximately equal to, or slightly larger than, the estimated reference diameter of the artery to be stented. The length of the stent should be chosen to adequately cover the length of the lesion. A compliance card is provided with the product to facilitate the selection of an appropriate size of stent.

## INSTRUCTIONS FOR USE

The Formula 418 Vascular Balloon-Expandable Stent System is used in conjunction with equipment required for a conventional PTA procedure including, but not limited to, a vascular access set, arterial sheath, guiding catheter, wire guide and inflation device.

### Lesion/Vessel Preparation

Stent placement must take into account proximal atherosclerotic plaque, which may inhibit advancement of the stent, as well as atherosclerotic plaque beyond the lesion which may prevent advancement of the device across the primary lesion. In preparing the vessel, optimal vasodilation, anticoagulant and antiplatelet therapy are essential. This device is intended for treatment following suboptimal angioplasty.

### Preparation of Balloon Catheter

1. Remove the Formula 418 Vascular Stent System from the package, remove protective sleeve from the distal tip of the catheter and inspect the stent to ensure it has not been damaged.
2. Prepare the balloon lumen with a standard 1:1 contrast-saline mixture as follows:

**WARNING: Do not attempt pre-inflation technique to purge the balloon lumen. Do not use air or any gaseous medium to inflate the balloon.**

- a. Using a 20 ml syringe containing 5 ml of contrast-saline mixture, apply negative pressure for 20-30 seconds.
- b. Release pressure, allowing negative pressure to draw mixture into balloon lumen.
- c. Detach syringe, leaving a meniscus of mixture on the hub of the balloon lumen.
- d. Prepare inflation device in standard manner and purge to remove all air from syringe and tubing.
- e. Attach inflation device to balloon lumen directly, ensuring no air bubbles remain at connection.
- f. Pull negative pressure on inflation device.

**CAUTION: Significant amounts of air in the balloon may cause uneven expansion of the stent and difficulty in deployment of the stent.**

3. Moisten the stent with heparinized saline.  
**CAUTION: Do not use gauze sponges, as fibers may disrupt stent.**
4. Flush the wire guide lumen of the balloon catheter in standard fashion to purge air.
5. Before inserting the Formula 418 Vascular Balloon-Expandable Stent System into the introducer or guiding catheter, be sure that a large-bore Tuohy-Borst "Y" adapter or introducer valve is on the guiding catheter/sheath.
6. Advance a 0.018 inch wire guide of appropriate length across the target lesion.
7. Insert the appropriate guiding catheter/sheath (see **Table 1**).
8. Advance the pre-mounted Formula 418 Vascular Balloon-Expandable Stent System over the wire into either the introducer valve or Tuohy-Borst "Y" adapter, making sure that the wire guide exits the wire guide port proximal to the balloon.
  - a. If using an introducer with a valve, make sure the flared end of the insertion tool (provided in the package) is loaded over the pre-mounted stent on the balloon catheter. Pass the insertion tool loaded with the pre-mounted stent through the introducer valve. Push the Formula 418 Vascular Balloon-Expandable Stent System into the body of the introducer. Slide the insertion tool proximally up the catheter shaft away from the

guiding catheter/sheath. A slight contact of the stent with the guiding catheter/sheath may be felt, but there must be no resistance.

- b. If using a Tuohy-Borst "Y" adapter, advance the pre-mounted Formula 418 Vascular Balloon-Expandable Stent System over the wire and into the fully opened large-bore Tuohy-Borst "Y" adapter. Gently advance the Formula 418 Vascular Balloon-Expandable Stent System completely through the Tuohy-Borst "Y" adapter, taking care to maintain wire guide position. Then advance the Formula 418 Vascular Balloon-Expandable Stent System into the guiding catheter/sheath. A slight contact of the stent with the guiding catheter/sheath may be felt, but there must be no resistance.

**WARNING: If resistance is encountered, do not force passage. Resistance may indicate damage to stent.**

### Positioning the Stent

1. Ensure guiding catheter/sheath stability before advancing the stent delivery balloon into the vascular artery.

**WARNING: If initial guiding catheter/sheath position is lost, avoid pulling or pushing the guiding catheter/sheath over the stent. If this is done, the distal end of the guiding catheter/sheath may damage the stent.**

2. Position the Formula Vascular Stent across the lesion, using both the distal and proximal balloon markers as reference points. Optimal placement requires the ends of the stent to extend beyond the margin of the lesion to be stented.

**WARNING: If the Formula 418 Vascular Balloon-Expandable Stent System does not readily advance through the vessel, do not force. If the stent will not advance in spite of good guiding catheter/sheath support, consider dilating proximal obstructing plaque, or changing the wire guide or the guiding catheter/sheath. (Refer to instructions for Removal of Unexpanded Stent.)**

**NOTE:** When stenting long lesions, be sure to cover the distal portion of the lesion first. It is very important not to leave distal lesions uncovered. Frequent injections of contrast around the uninflated Formula 418 Vascular Balloon-Expandable Stent System will allow visualization of the extent of the lesion and facilitate accurate placement of the stent.

### Balloon Expansion/Stent Deployment

Prior to stent expansion, utilize high-resolution fluoroscopy to verify that the stent has not been dislodged during positioning.

1. To expand the stent, inflate the balloon with a 1:1 contrast-saline mixture to the recommended expansion pressure indicated on product label. Do not move stent system during deployment.
2. Complete expansion and apposition of the stent against the vessel wall is necessary for clinical success. Do not exceed rated burst pressure of the balloon as indicated on the product label.
3. Once the stent has been deployed, post-deployment inflation is at the discretion of the operator to achieve optimum angiographic appearance. The stent may be post-dilated up to 1 mm beyond nominal labeled stent diameter.

### Balloon Deflation and Removal

1. Completely deflate the balloon, by pulling negative pressure with the inflation device or a 20 ml syringe. This usually requires 10 seconds or less.

**WARNING: Allow enough time for the balloon to fully deflate prior to removal.**

2. Slowly withdraw the balloon catheter from the stent while maintaining negative pressure on the balloon. Maintain position of the guiding catheter/sheath to prevent it from being drawn into the vessel. Observe under fluoroscopy to ensure that the balloon disengages from the stent.

### Removal of Unexpanded Stent

Do not attempt to pull an unexpanded stent back into the guiding catheter/sheath. The Formula 418 Vascular Balloon-Expandable Stent System should be withdrawn until the proximal end of the stent is aligned with the distal tip of the guiding catheter/sheath. Withdraw the guiding catheter/sheath and stent delivery system as a single unit, leaving the wire guide in place.

**WARNING: If stent is removed, do not attempt to reuse the device. Damage to the stent may occur upon removal.**

### Stent-Assisted, High-Pressure Balloon Angioplasty

A deployed Formula 418 Vascular Balloon-Expandable Stent may be further expanded using a low-profile, high-pressure balloon catheter. The stent should have been well embedded into the vessel wall by the deployment balloon before the lesion is recrossed. If recrossing is necessary, the position of the stent should be monitored and care should be taken to avoid displacing the stent with subsequent manipulation of wire guides and balloons. An appropriate balloon catheter should be used with a very flexible wire guide.

If a Formula 418 Vascular Balloon-Expandable Stent is dilated with a high-pressure balloon, the balloon should be appropriately sized for the vessel diameter. The stent may be post-dilated up to 1 mm beyond nominal labeled stent diameter. Use of a new balloon for post-deployment dilatation is suggested. Avoid balloon catheters with stiff tips that might displace the stent during advancement. Avoid balloons that “wing out” and/or re-wrap inconsistently, which might displace the stent during withdrawal of the balloon from the stent

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
**Cook Medical Europe Ltd**