EN

2

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Brugsanvisning

T _ F T C W G _ R E V O

FALLOPIAN TUBE CATHETERIZATION WIRE GUIDE WITH PLATINUM TIP

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Fallopian Tube Catheterization Wire Guide with Platinum Tip is nitinol with a stainless steel platinum tip coil.

INTENDED USE

The Fallopian Tube Catheterization Wire Guide with platinum tip is to be used through a dimensionally compatible (> 0.018 inch inner diameter) uterine access device or a selective salpingography catheter to assist in the evaluation of tubal patency.

CONTRAINDICATIONS

This device should not be used in the presence of:

- Pregnancy or suspected pregnancy
- Active or recent pelvic infection
- Unexplained or severe uterine bleeding
- Gynecologic malignancy (suspected or known)
- Known allergy to dye or contrast medium

WARNINGS

This wire guide is intended only to facilitate placement of an inner catheter; it is not intended for tubal recanalization and should not be advanced beyond the tubal isthmus.

PRECAUTIONS

- This product is intended for use by a physician or properly licensed practitioner trained and experienced in diagnostic and interventional techniques. Standard techniques for HSG and SSG procedures should be employed.
- This product is not intended for complete catheterization of the fallopian tube or tubal recanalization. There is no data available showing any clinical benefit for use in this manner.

POTENTIAL ADVERSE EVENTS

- Tubal perforation/dissection
- Damage to normal tubes
- Extravasation
- Infertility
- Ectopic pregnancy
- Discomfort, pain/cramping
- Bleeding
- Infection

INSTRUCTIONS FOR USE

1. Follow the instructions for use provided with the uterine access catheter and the selective salpingography (SSG) catheter with a working channel that accepts at least a 0.018 wire guide.
2. While scanning with fluoroscopy, insert the wire guide through the SSG catheter.
3. The wire guide will be visible under radiographic guidance and is designed to facilitate proximal fallopian tube access; do not advance the wire guide beyond the fallopian tube isthmus.
4. Advance the wire guide up to the distal fallopian tube isthmus unless resistance is met; if resistance is met, visualize the fallopian tube anatomy under fluoroscopy.
5. After the SSG catheter is properly positioned, remove the wire guide and proceed by following the SSG catheter's instructions for use.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2017
Bloomington, IN 47404 U.S.A.