COOK

MEDICAL

EN
3

Fallopian Tube Catheterization Set

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Bruksanvisning

T _ F T C _ R E V 6

FALLOPIAN TUBE CATHETERIZATION SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Fallopian Tube Catheterization Set consists of catheters and wire guides necessary to catheterize the fallopian tubes. A uterine access catheter with a working channel greater than 6 French or 9 French may be used to facilitate placement.

INTENDED USE / INDICATIONS FOR USE

The Fallopian Tube Catheterization Set is intended to be used through a uterine access device for fluoroscopic selective catheterization of the proximal fallopian tubes and injection of contrast medium in the evaluation of tubal patency.

CONTRAINDICATIONS

- Pregnancy or suspected pregnancy
- Active or recent pelvic infection
- Unexplained or severe uterine bleeding
- Gynecologic malignancy (suspected or known)
- Known allergy to dye or contrast medium

WARNINGS

- Do not use power injection systems for injection of contrast medium through the catheter lumens. Fallopian tube dissection may occur.
- The wire guides in the Fallopian Tube Catheterization Set are intended only to facilitate placement of the inner catheters. They are not intended for tubal recanalization and should not be advanced beyond the tubal isthmus.
- These catheter systems are only designed to access the fallopian tube from the uterine cavity. They should only be advanced as far as necessary to achieve a meaningful diagnostic study and no further.
- These catheter systems are not designed nor intended for embryo transfer or other assisted reproduction procedures.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in fallopian tube catheterization techniques. Standard techniques for placement of catheters and wire guides should be employed.
- Fallopian Tube Catheterization Set is not intended for complete catheterization of the fallopian tube or tubal recanalization. There are no data available showing any clinical benefit for use in this manner.

POTENTIAL ADVERSE EVENTS

The following adverse events have been reported or are of concern with the use of this type of device. Patients must be advised of the following risks:

- Ectopic pregnancy. Patients who have previously blocked fallopian tubes may be at increased risk for ectopic pregnancy after the tubes are opened. The magnitude of the risk is unknown.
- Damage to normal tubes. It is unknown whether this procedure may damage normal fallopian tubes or put the patient who has normal tubes at any increased risk for subsequent ectopic pregnancy.
- Tubal dissection. It is possible that this catheter may dissect into the wall of the tube, under the epithelial lining, and cause blockage or narrowing of the lumen of the fallopian tube. This may result in sterility, infertility or increased risk of ectopic pregnancy.
- Pain or discomfort. Some patients may experience cramping or other abdominal pain.
- Extravasation. Dye or contrast medium that exceeds normal anatomic constraints of the tube may be indicative of extravasation and consequent tubal damage.
- Tubal perforation, bleeding, or infection. There is a small risk of tubal perforation, bleeding and infection associated with the use of this type of device. The magnitude of the risk has not been quantified.

INSTRUCTIONS FOR USE

1. If using a uterine access device greater than 9 French, follow the Instructions for Use provided. Introduce the uterine access device into the cervix and affix it in position either on the cervix or within the uterine cavity. Alternatively, the 9 French introducing catheter may be disregarded and a uterine access catheter with greater than 6 French working channel may be used. If so, place the uterine access catheter per instructions, then advance the 5.5 French torque control catheter until the distal tip is positioned within the lower portion of the uterine cavity, then skip to step 6.
2. Advance the 5.5 French torque-control delivery catheter through the 9.0 French introducing catheter until both catheter tips are even.
3. If needed, connect the Tuohy-Borst adapter to the access/injection port of the uterine access catheter.

4. Insert the introducing catheter/delivery catheter assembly through the working channel of the uterine access catheter, advancing until the distal tips of the catheters are positioned within the lower portion of the uterine cavity.
5. Tighten the Tuohy-Borst adapter on the access catheter, if applicable, to prevent movement of the catheter and leakage of contrast medium.
6. Attach a syringe filled with the appropriate contrast medium to the Luer lock fitting on the delivery catheter. Under fluoroscopic guidance, inject contrast medium through the delivery catheter to visualize the uterine cavity and tubal ostium.
7. Detach the syringe and introduce the .035 inch Safe-T-J® wire guide through the delivery catheter into the uterine cavity.
8. With the wire guide leading and under fluoroscopic guidance, advance the catheter and wire guide together toward the tubal ostium and wedge the catheter tip into the ostium.
NOTE: If difficulty is encountered placing the delivery catheter at the tubal ostium, remove the Safe-T-J wire guide and replace it with the .035 inch straight wire guide.
9. Leaving the delivery catheter tip wedged in the tubal ostium, remove the wire guide.
10. Attach a syringe filled with the appropriate contrast medium. Inject 2-10 mL of contrast medium into the tubal ostium while observing under fluoroscopy for spillage from the fimbrial portion of the tube. **If such spillage is observed, tubal patency has been confirmed, and the procedure is complete.**
NOTE: If spillage of contrast medium from the fimbrial portion of the tube into the peritoneal cavity is not observed, it will be necessary to introduce the 3.0 French delivery catheter into the proximal fallopian tube to confirm tubal blockage or rule out a false positive result due to temporary tubal spasm.
11. Detach the syringe.
12. Introduce the 3.0 French delivery catheter and .015 inch mandril wire guide into the 5.5 French delivery catheter.
13. Advance the catheter and mandril wire guide together until the distal tip of the 3.0 French delivery catheter is even with the tip of the 5.5 French delivery catheter.
14. Begin to cannulate the intramural portion of the fallopian tube by first advancing the mandril wire guide. **NOTE:** The tip of the mandril wire guide will become more flexible as the length exiting the catheter tip increases.
15. Under fluoroscopic guidance, advance the mandril wire guide beyond the utero-tubal junction and into the tubal isthmus unless significant resistance is met.
16. Coaxially, advance the 3.0 French delivery catheter over the positioned mandril wire guide for a short distance or until resistance is felt.
WARNING: If significant resistance is felt, do not attempt to advance the catheter.
WARNING: Do not attempt to advance the catheter or mandril wire guide beyond the tubal isthmus.
17. Remove the mandril wire guide and attach a syringe filled with the appropriate contrast medium.
18. Inject 2-10 mL of contrast medium into the tubal ostium while observing fluoroscopically for spillage from the fimbrial portion of the tube. If spillage is observed, tubal patency is confirmed.
19. Using these same instructions, confirm tubal patency in the other fallopian tube.
20. Once tubal evaluation has been performed for both fallopian tubes, remove the catheters.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.