COOK**MEDICAL**

- EN**
3
Cope Pediatric Gastrointestinal Suture Anchor Sets
Instructions for Use
- CS**
5
Soupravy pediatrického gastrointestinálního kotvičního vlákna Cope
Návod k použití
- DA**
8
Cope pædiatrisk gastrointestinale suturankersæt
Brugsanvisning
- DE**
10
Cope pædiatrische gastrointestinale Fadenanker-Sets
Gebrauchsanweisung
- EL**
13
Σετ παιδιατρικού γαστρεντερικού αγκίστρου ράμματος Cope
Οδηγίες χρήσης
- ES**
16
Equipos de anclajes de sutura gastrointestinal pediátricos Cope
Instrucciones de uso
- FR**
19
Set d'ancre de suture gastro-intestinale pédiatrique Cope
Mode d'emploi
- HU**
22
Cope gyermekgyógyászati gastrointestinalis varratrögzőtő készletek
Használati utasítás
- IT**
25
Set di suture pediatrico con ancora gastrointestinale Cope
Istruzioni per l'uso
- NL**
28
Cope pediatische gastro-intestinale hechtankersets
Gebruiksaanwijzing
- NO**
31
Cope pediatrike gastrointestinale suturankersett
Bruksanvisning
- PL**
33
Pediatryczne zestawy kotwicy szwu do przewodu pokarmowego Cope
Instrukcja użycia
- PT**
36
Conjuntos pediátricos de ancoragem de sutura gastrointestinal Cope
Instruções de utilização
- SV**
39
Cope pediatrika gastrointestinala suturankersett
Bruksanvisning

T _ G I A P _ R E V 1

COPE PEDIATRIC GASTROINTESTINAL SUTURE ANCHOR SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Cope Pediatric Gastrointestinal Suture Anchor Set consists of an introducer needle with a preloaded suture anchor, a .025 inch wire guide with a spring coil tip, and a second suture anchor in a loading cartridge.

INTENDED USE

The Cope Pediatric Gastrointestinal Suture Anchor Set is intended for anchoring the anterior wall of the stomach to the abdominal wall prior to the introduction of interventional catheters and can stay in place for up to 14 days in infant, child, and adolescent populations.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Cope Pediatric Gastrointestinal Suture Anchor is **MR Conditional**. A patient with this device may be scanned safely in a MR system meeting the following conditions.

- Static magnetic field of 1.5 Tesla or 3.0 Tesla
- Maximum spatial magnetic gradient of 1900 Gauss/cm or less (19 T/m)
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) for 15 minutes of scanning or less (i.e., per scanning sequence)

Under the scan conditions defined above, the Cope Pediatric Gastrointestinal Suture Anchor device is expected to produce a maximum temperature rise of less than 2°C after 15 minutes of continuous scanning.

In nonclinical testing, the image artifact caused by the device extends approximately 25 mm from this device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

CONTRAINDICATIONS

None known

WARNINGS

Suture anchors may migrate and be retained in the gastric wall. In rare circumstances the suture anchors have exited through the skin.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of gastrointestinal suture anchors should be employed.
- Manipulation of the product requires fluoroscopic control.
- Withdrawal or manipulation of the distal spring coil portion of the mandril wire guide through the needle tip may result in breakage.

INSTRUCTIONS FOR USE

1. Insufflate the stomach with approximately 500 cc of air through a nasogastric tube.
2. Localize the puncture site and administer local anesthetic into the subcutaneous tissue down to the peritoneum. **NOTE:** The preferred puncture site into the gastric lumen is at the junction of the left anterior rib cage and the lateral border of the rectus abdominis muscle, provided there is no interposition of the liver or colon.
3. Remove the plug from the introducer needle and attach a 10 mL syringe, half-filled with contrast medium, to the introducer needle that has been preloaded with the suture anchor.
4. Advance the introducer needle through the anterior abdominal wall and then rapidly thrust the needle another 2-5 cm into the air-filled stomach. Care should be taken to not puncture through the back wall of the stomach. **Air should be freely aspirated. Injection of contrast medium under fluoroscopy will aid in confirming proper intragastric needle position.**
5. Remove the 10 mL syringe from the needle hub. While maintaining slight tension on the trailing suture, introduce the distal spring coil of the provided wire guide into the introducer needle and use it to push the suture anchor out of the needle into the stomach cavity. If no wire guide is provided, use a 0.025 inch diameter wire guide. **Maintaining slight tension on the suture during anchor introduction will prevent it from bunching and lodging within the needle lumen.**
6. Remove the introducer needle over the wire guide and trailing suture.
7. With the wire guide still in position, apply traction to the suture to pull the anterior wall of the stomach against the abdominal wall.
8. While maintaining traction on the suture anchor, secure the suture to the skin under slight tension.
9. Remove the wire guide.
10. To place the second suture anchor, reintroduce the introducer needle as described in step 4.
11. If one exists, advance the loading cartridge containing the second suture anchor into the hub of the introducer needle. If no loading cartridge exists, insert the remaining suture anchor into the introducer needle.

12. If applicable, remove the plug from the introducer cartridge. Advance the wire guide through the cannula of the loading cartridge into the introducer needle, pushing the second suture anchor into the abdominal cavity.
13. Follow steps 6 through 8 with the second suture anchor.
14. The wire guide may now be used to facilitate passage of fascial dilators and introduction of the gastrostomy catheter. Take care not to advance the entire distal spring coil segment of the wire guide into the abdominal cavity.
15. While maintaining traction on the suture anchors, wrap the ends of the suture around 2x2 gauze rolls and secure in place. Apply a standard wound dressing around the insertion site.

NOTE: The sutures may be left in place for up to two weeks while tract formation occurs, or cut earlier when deemed appropriate by the physician. Cutting the sutures releases the anchors into the stomach, allowing their passage via the gastrointestinal system.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2023
Bloomington, IN 47404 U.S.A.