Enterostomy Suture Anchor Set

Instructions for Use

EN
5

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Brugsanvisning

DA
6

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
Fig. 4

ENTEROSTOMY SUTURE ANCHOR SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Enterostomy Suture Anchor Set consists of an introducer needle, Amplatz Ultra Stiff Wire Guide, loading cartridge with two preloaded suture anchors and Coons Dilator.

INTENDED USE

The Enterostomy Suture Anchor Set is intended for anchoring the wall of a hollow viscus to the abdominal wall prior to introduction of interventional catheters.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for percutaneous placement should be employed.

INSTRUCTIONS FOR USE

Illustrations

1. Insufflate the hollow viscus with air through an enema balloon catheter. Use of glucagon to control bowel peristalsis is recommended.
2. **NOTE:** Ultrasound is useful in determining anatomic structures before needle introduction. For cecostomy, identify the gallbladder and urinary bladder before initial puncture.
3. Localize the puncture site and administer local anesthesia with a 27 gage needle (not included in the set) into the subcutaneous tissue and down to the peritoneum using fluoroscopic control.
4. Attach a Luer slip 10 ml syringe, half filled with (full-strength) 300 mg/ml concentration contrast medium, to the 18 gage introducer needle. Make a small incision with a #11 blade at the entry site.
5. Advance the introducer needle through the anterior abdominal wall and then, under fluoroscopic guidance, rapidly thrust the needle another 2 to 3 cm into the air-filled hollow viscus.
NOTE: Air should be freely aspirated. Injection of contrast medium under fluoroscopy must be performed to ensure proper needle position. (Fig. 1)
6. Remove the 10 ml syringe from the needle hub. Introduce the metal end of the loading cartridge containing the two suture anchors into the base of the introducer needle and Luer lock into place. (Fig. 2)
NOTE: Make sure the cannula tip of the loading cartridge is seated within the hub of the introducer needle to prevent suture anchors from jamming.
7. Advance the distal tip of the Amplatz Ultra Stiff Wire Guide into the loading cartridge. Push the two anchors out of the loading cartridge and through the introducer needle into the hollow viscus cavity. (Fig. 3) Confirm the suture anchors are within the hollow viscus.
8. Disconnect the loading cartridge and withdraw over the wire guide to expose the suture ends. Discard the cartridge. Remove the introducer needle over the wire guide and expose the suture ends, maintaining slight tension on the suture to prevent bunching. (Fig. 4)
NOTE: Normally, the sutures are left in place for two weeks to enable the tract to mature.
9. With the wire guide still in position, apply light traction to the sutures using hemostats to secure the anterior wall of the hollow viscus against the abdominal wall.
10. The wire guide may now be used to facilitate passage of fascial dilators and introduction of a procedural catheter.
11. The wire guide may be removed after introduction of dilators and/or procedural catheters is complete.
12. While maintaining traction on the suture anchors, wrap the ends of the suture around 2x2 gauze rolls and secure in place. Apply a standard wound dressing around the insertion site. **NOTE:** The suture anchors may be left in place for two weeks while tract formation occurs. Cutting the suture end at the skin level releases the anchors into the organ, allowing their passage via the gastrointestinal system.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.