<b>EN</b> 5	<b>Enterostomy Suture Anchor Set</b> Instructions for Use
<b>CS</b> 6	<b>Souprava kotvicího vlákna pro enterostomii</b> Návod k použití
<b>DA</b> 8	<b>Enterostomi-suturankersæt</b> Brugsanvisning
<b>DE</b> 9	<b>Fadenanker-Set für die Enterostomie</b> Gebrauchsanweisung
<b>EL</b> 11	<b>ΣΕΤ αγκίστρου ράμματος εντεροστομίας</b> Οδηγίες χρήσης
<b>ES</b> 13	<b>Equipo de anclajes de sutura para enterostomía</b> Instrucciones de uso
<b>FR</b> 15	<b>Set d'ancre de suture pour entérostomie</b> Mode d'emploi
<b>HU</b> 16	<b>Enterosztómias varratrögző készlet</b> Használati utasítás
<b>IT</b> 18	<b>Set di suture con ancora per enterostomia</b> Istruzioni per l'uso
<b>NL</b> 20	<b>Hechtankerset voor enterostomie</b> Gebruiksaanwijzing
<b>NO</b> 21	<b>Enterostomisuturankersett</b> Bruksanvisning
<b>PL</b> 23	<b>Zestaw kotwicy szwu do enterostomii</b> Instrukcja użycia
<b>PT</b> 25	<b>Conjunto de ancoragem de sutura de enterostomia</b> Instruções de utilização
<b>SV</b> 26	<b>Suturankerset för enterostomi</b> Bruksanvisning


Fig. 3


Fig. 4

## ENTEROSTOMY SUTURE ANCHOR SET

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### DEVICE DESCRIPTION

The Enterostomy Suture Anchor Set consists of an introducer needle, Amplatz Ultra Stiff Wire Guide, loading cartridge with two preloaded suture anchors, plug to secure the suture anchors, and Coons Dilator.

### INTENDED USE

The Enterostomy Suture Anchor Set is intended for anchoring the wall of a hollow viscus to the abdominal wall prior to the introduction of interventional catheters and can stay in place for up to 14 days in child, adolescent, and adult populations.

### MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Enterostomy Suture Anchor is **MR Conditional**. A patient with this device may be scanned safely in a MR system meeting the following conditions.

- Static magnetic field of 1.5 Tesla or 3.0 Tesla
- Maximum spatial magnetic gradient of 1900 Gauss/cm or less (19 T/m)
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) for 15 minutes of scanning or less (i.e., per scanning sequence).

Under the scan conditions defined above, the Enterostomy Suture Anchor device is expected to produce a maximum temperature rise of less than 2°C after 15 minutes of continuous scanning.

In nonclinical testing, the image artifact caused by the device extends approximately 25 mm from this device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

### For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International  
2323 Colorado Avenue  
Turlock, CA 95382

Phone: 888-633-4298 (toll free)  
209-668-3333 from outside the US

Fax: 209-669-2450

Web: [www.medicalert.org](http://www.medicalert.org)

### CONTRAINDICATIONS

None known

### WARNINGS

Suture anchors may migrate and be retained in the gastric wall. In rare circumstances the suture anchors have exited through the skin.

### PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of gastrointestinal suture anchors should be employed.
- Manipulation of the product requires fluoroscopic control.
- Withdrawal or manipulation of the distal spring coil portion of the mandril wire guide through the needle tip may result in breakage.

### INSTRUCTIONS FOR USE

1. Insufflate the hollow viscus with air through an enema balloon catheter. Use of glucagon to control bowel peristalsis is recommended.
2. **NOTE:** Ultrasound is useful in determining anatomic structures before needle introduction. For cecostomy, identify the gallbladder and urinary bladder before initial puncture.
3. Localize the puncture site and administer local anesthesia with a 27 gage needle (not included in the set) into the subcutaneous tissue and down to the peritoneum using fluoroscopic control.
4. Attach a Luer slip 10 mL syringe, half-filled with (full-strength) contrast medium, to the introducer needle. Make a small incision with a #11 blade at the entry site.
5. Advance the introducer needle through the anterior abdominal wall, and then under fluoroscopic guidance rapidly thrust the needle another 2 to 3 cm into the air-filled hollow viscus.

**NOTE:** Air should be freely aspirated. Injection of contrast medium under fluoroscopy must be performed to ensure proper needle position. (Fig. 1)

6. Remove the 10 mL syringe from the needle hub. Introduce the metal end of the loading cartridge containing the two suture anchors into the base of the introducer needle and Luer lock into place. (Fig. 2)

**NOTE:** Make sure the cannula tip of the loading cartridge is seated within the hub of the introducer needle to prevent suture anchors from jamming.

7. Remove the plug from the proximal end of the loading cartridge and advance the distal tip of the Amplatz Ultra Stiff Wire Guide into the loading cartridge. Push the two anchors out of the loading cartridge and through the introducer needle into the hollow viscus cavity. (Fig. 3) Confirm the suture anchors are within the hollow viscus.
8. Disconnect the loading cartridge and withdraw over the wire guide. Discard the cartridge. Remove the introducer needle over the wire guide, maintaining slight tension on the sutures to prevent bunching. (Fig. 4)  
**NOTE:** Normally, the sutures are left in place for two weeks to enable the tract to mature.
9. With the wire guide still in position, apply light traction to the sutures using hemostats to secure the anterior wall of the hollow viscus against the abdominal wall.
10. The wire guide may now be used to facilitate passage of fascial dilators and introduction of a procedural catheter.
11. The wire guide may be removed after introduction of dilators and/or procedural catheters is complete.
12. While maintaining traction on the suture anchors, wrap the ends of the suture around 2x2 gauze rolls and secure in place. Apply a standard wound dressing around the insertion site.  
**NOTE:** The suture anchors may be left in place for two weeks while tract formation occurs. Cutting the suture end at the skin level releases the anchors into the organ, allowing their passage via the gastrointestinal system.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd