COOK**MEDICAL****CE**
0123

EN 3	Cope Gastrointestinal Suture Anchor Sets Instructions for Use
CS 5	Soupravy gastrointestinálního kotvicího vlákna Cope Návod k použití
DA 8	Cope gastrointestinale suturankersæt Brugsanvisning
DE 11	Gastrointestinale Nahtanker-Sets nach Cope Gebrauchsanweisung
EL 14	Σετ γαστρεντερικού αγκίστρου ράμματος Cope Οδηγίες χρήσης
ES 17	Equipos de anclajes de sutura gastrointestinal Cope Instrucciones de uso
FR 20	Sets d'ancre de suture gastro-intestinale Cope Mode d'emploi
HU 23	Cope gasztrointesztinális varratrögzítő készletek Használati utasítás
IT 26	Set di suture con ancora gastrointestinale Cope Istruzioni per l'uso
NL 29	Cope gastro-intestinale hechtankersets Gebruiksaanwijzing
NO 32	Cope gastrointestinale suturankersett Bruksanvisning
PL 34	Zestawy kotwicy szwu Cope do przewodu pokarmowego Instrukcja użycia
PT 38	Conjuntos de ancoragem em T para sutura gastrointestinal Cope Instruções de utilização
SV 41	Cope gastrointestinale suturankersæt Bruksanvisning

T _ G I A _ R E V 6

COPE GASTROINTESTINAL SUTURE ANCHOR SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Cope Gastrointestinal Suture Anchor Set consists of an introducer needle with a preloaded suture anchor, a special .018/.035 inch wire guide with a spring coil tip, and a second suture anchor in a loading cannula.

INTENDED USE

The Cope Gastrointestinal Suture Anchor Set is intended for anchoring the anterior wall of the stomach to the abdominal wall prior to introduction of interventional catheters.

MRI INFORMATION

Non-clinical testing has demonstrated that the Cope Gastrointestinal Suture Anchor is **MR Conditional** according to ASTM F2503, *Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment*.

A patient with this metallic anchor may be scanned safely anytime after placement under the following conditions.

- Static magnetic field of 3.0 Tesla or less
- Maximum spatial magnetic gradient of 1600 Gauss/cm or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) for 15 minutes of scanning or less (i.e., per scanning sequence)

Static Magnetic Field

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

In non-clinical testing, the Cope Gastrointestinal Suture Anchor produced a maximum temperature rise of 1.6 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3.0 Tesla System (General Electric Excite, HDx, Software 14X.M5) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole-body-averaged value of 2.7 W/kg).

Image Artifact

MR image quality may be compromised if the area of interest is within approximately 25 mm of the position of the Cope Gastrointestinal Suture Anchor as found during non-clinical testing using T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic anchor.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of gastrointestinal suture anchors should be employed.
- Manipulation of the product requires fluoroscopic control.
- Withdrawal or manipulation of the distal spring coil portion of the mandril wire guide through the needle tip may result in breakage.

INSTRUCTIONS FOR USE

1. Insufflate the stomach with approximately 500 cc of air through a nasogastric tube.
2. Localize the puncture site and administer local anesthetic to the peritoneum. **NOTE:** The preferred puncture site into the gastric lumen is at the junction of the left anterior rib cage and the lateral border of the rectus abdominis muscle, provided there is no interposition of the liver or colon.
3. Attach a 10 ml syringe, half filled with contrast medium, to the introducer needle that has been preloaded with the suture anchor.
4. Advance the introducer needle through the anterior abdominal wall and then rapidly thrust the needle another 2-5 cm into the air-filled stomach. **Air should be freely aspirated. Injection of contrast medium under fluoroscopy will aid in confirming proper intragastric needle position.**
5. Remove the 10 ml syringe from the needle hub. While maintaining slight tension on the trailing suture, introduce the distal spring coil .035 inch portion of the wire guide into the needle and use it to push the suture anchor out of the needle into the stomach cavity. **Maintaining slight tension on the suture during anchor introduction will prevent it from bunching and lodging within the needle lumen.**

6. Remove the introducer needle over the wire guide and trailing suture needle. **The .018 inch part of the wire guide must be positioned adjacent to the suture needle to allow clearance and removal of the introducer needle.**
7. With the wire guide still in position, apply traction to the suture to pull the anterior wall of the stomach against the abdominal wall.
8. While maintaining traction on the suture anchor, secure the suture to the skin under slight tension.
9. Remove the wire guide.
10. To place the second suture anchor, reintroduce the introducer needle as described in step 4.
11. Advance the loading cannula containing the second suture anchor into the hub of the introducer needle.
12. Advance the wire guide through the cannula into the introducer needle, pushing the second suture anchor into the abdominal cavity.
13. Follow steps 6 through 8 with the second suture anchor.
14. The wire guide may now be used to facilitate passage of fascial dilators and introduction of the gastrostomy catheter. Take care not to advance the entire distal spring coil segment of the wire guide into the abdominal cavity.

NOTE: The sutures may be left in place for up to two weeks while tract formation occurs, or cut earlier when deemed appropriate by the physician. Cutting the sutures releases the anchors into the stomach, allowing their passage via the gastrointestinal system.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd