EN
3**Single-Lumen and Double-Lumen Carey-Alzate-Coons Gastrojejunostomy Catheters and Sets**

Instructions for Use

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Bruksanvisning

T _ G J S _ R E V 4

SINGLE-LUMEN AND DOUBLE-LUMEN CAREY-ALZATE-COONS GASTROJEJUNOSTOMY CATHETERS AND SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Carey-Alzate-Coons Gastrojejunostomy Sets consist of a single- or double-lumen Friction-Lock™ Malecot Gastrojejunostomy Catheter, and Peel-Away® introducer, as well as an appropriately sized access needle, wire guide, dilators and other accessories for use in percutaneous placement of a gastrojejunostomy tube.

Sets also include a Cope Gastrointestinal Suture Anchor Set, consisting of an introducer needle with a preloaded suture anchor, a special .018/.035 inch wire guide with a spring coil tip, and a second suture anchor in a loading cannula.

INTENDED USE

The Single-Lumen Carey-Alzate-Coons Gastrojejunostomy Catheters and Sets are intended for percutaneous placement of a catheter for gastroenteric feeding. The Double-Lumen Carey-Alzate-Coons Gastrojejunostomy Catheters and Sets are intended for placement of a double-lumen catheter for simultaneous feeding in the jejunum and aspiration of gastric contents in the stomach.

The Cope Gastrointestinal Suture Anchor Set is intended for anchoring the anterior wall of the stomach to the abdominal wall prior to introduction of interventional catheters.

MRI INFORMATION

Cope Gastrointestinal Suture Anchor

Non-clinical testing has demonstrated that the Cope Gastrointestinal Suture Anchor is **MR Conditional** according to ASTM F2503, *Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment*. A patient with this metallic anchor may be scanned safely anytime after placement under the following conditions.

- Static magnetic field of 3.0 Tesla or less
- Maximum spatial magnetic gradient of 1600 Gauss/cm or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) for 15 minutes of scanning or less (i.e., per scanning sequence)

Static Magnetic Field

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

In non-clinical testing, the Cope Gastrointestinal Suture Anchor produced a maximum temperature rise of 1.6 °C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3.0 Tesla System (General Electric Excite, HDx, Software 14X.M5) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole-body-averaged value of 2.7 W/kg).

Image Artifact

MR image quality may be compromised if the area of interest is within approximately 25 mm of the position of the Cope Gastrointestinal Suture Anchor as found during non-clinical testing using T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic anchor.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedAlert Foundation. The MedAlert Foundation can be contacted in the following manners:

Mail: MedAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in gastrojejunostomy techniques. Standard techniques for percutaneous placement of gastrojejunostomy tubes should be employed.
- Single-lumen catheter sets include a Silhouette™ Transitionless Peel-Away® introducer, which has a uniquely fabricated transition between dilator and sheath. Do not turn or pull back dilator from sheath prior to introduction. Manipulation of dilator will result in misalignment of sheath and dilator. Do not attempt to flush sheath prior to insertion.
- A TFE-coated wire guide must be used with this catheter.
- Manipulation of the product requires fluoroscopic control.
- Withdrawal or manipulation of the distal spring coil portion of the mandril wire guide through the needle tip may result in breakage.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

Placement of Cope Gastrointestinal Suture Anchor Set

1. Insufflate the stomach with approximately 500 cc of air through a nasogastric tube.
2. Localize the puncture site and administer local anesthetic to the peritoneum. **NOTE:** The preferred puncture site into the gastric lumen is at the junction of the left anterior rib cage and the lateral border of the rectus abdominis muscle, provided there is no interposition of the liver or colon.
3. Attach a 10 ml syringe, half filled with contrast medium, to the introducer needle that has been preloaded with the suture anchor.
4. Advance the introducer needle through the anterior abdominal wall and then rapidly thrust the needle another 2-5 cm into the air-filled stomach. **Air should be freely aspirated. Injection of contrast medium under fluoroscopy will aid in confirming proper intragastric needle position.**
5. Remove the 10 ml syringe from the needle hub. While maintaining slight tension on the trailing suture, introduce the distal spring coil .035 inch portion of the wire guide into the needle and use it to push the suture anchor out of the needle into the stomach cavity. **Maintaining slight tension on the suture during anchor introduction will prevent it from bunching and lodging within the needle lumen.**
6. Remove the introducer needle over the wire guide and trailing suture needle. **The .018 inch part of the wire guide must be positioned adjacent to the suture needle to allow clearance and removal of the introducer needle.**
7. With the wire guide still in position, apply traction to the suture to pull the anterior wall of the stomach against the abdominal wall.
8. While maintaining traction on the suture anchor, secure the suture to the skin under slight tension.
9. Remove the wire guide.
10. To place the second suture anchor, reintroduce the introducer needle as described in step 4.
11. Advance the loading cannula containing the second suture anchor into the hub of the introducer needle.
12. Advance the wire guide through the cannula into the introducer needle, pushing the second suture anchor into the abdominal cavity.
13. Follow steps 6 through 8 with the second suture anchor.
14. The wire guide may now be used to facilitate passage of fascial dilators and introduction of the gastrostomy catheter. Take care not to advance the entire distal spring coil segment of the wire guide into the abdominal cavity.

NOTE: The suture may be left in place for two weeks while tract formation occurs, or it may be cut after gastrostomy catheter placement. This releases the anchor into the stomach, allowing its passage via the gastrointestinal system.

Placement of Wire Guide

1. Using the supplied introducer needle, make a puncture adjacent to the site used to introduce the suture anchors, advancing the needle into the stomach. **NOTE:** This puncture should be performed at an approximately 45-degree angle cephalad.
2. Insert the supplied heavy-duty straight safety wire through the introducer needle. Under fluoroscopic guidance, manipulate the wire guide into the antrum of the stomach. Most commonly, the wire guide coils in the fundus and then exits to the antrum.
3. Remove the introducer needle and advance the seeking catheter over the wire guide.

4. Advance the wire guide and seeking catheter until the catheter tip is at the pylorus, then negotiate through the pylorus and advance the wire guide into the duodenum. **NOTE:** Take care not to coil too much seeking catheter in the fundus of the stomach. Counter forces, which make it difficult to advance the seeking catheter, may develop. Frequently, a rotational motion on the seeking catheter allows easier passage over the wire guide.
5. Continue advancing the wire guide and seeking catheter beyond the ligament of Treitz, then remove the seeking catheter.

Dilation and Placement of Peel-Away® Introducer

For placement of a single-lumen gastrojejunostomy tube introducer:

1. Dilate the tract through the abdominal wall with the 10.0 French and 14.0 French dilators in sequence.
NOTE: Because the stomach wall is relatively pliable, it is often difficult to be certain that the dilator is within the stomach. Ensure that each dilator enters the stomach by visualizing, fluoroscopically, the dilator curving over the wire guide positioned in the fundus.
2. Remove the dilator, leaving the wire guide in place.
3. Advance the Peel-Away introducer over the wire guide.
4. Ensure that the introducer is within the stomach by visualizing, fluoroscopically, the tip of the introducer beginning to curve in the fundus of the stomach.
5. Turn the dilator one-half turn to break the seal between the sheath and dilator.
6. Pull back the dilator to remove, leaving the sheath portion of introducer and wire guide in place.

For placement of a double-lumen gastrojejunostomy tube introducer:

1. Dilate the tract through the abdominal wall using a coaxial dilator technique. Begin by first inserting the supplied guiding catheter over the wire guide into the fundus of the stomach. Then serially dilate over the guiding catheter and wire guide, beginning with the 12.0 French dilator and continuing up in size. Do not advance the distal tip of the dilator beyond the distal tip of the guiding catheter during dilation.
NOTE: Because the stomach wall is relatively pliable, it is often difficult to be certain that the dilator is within the stomach. Ensure that each dilator enters the stomach by visualizing, fluoroscopically, the dilator curving over the wire guide positioned in the fundus.
2. With the guiding catheter still in place, insert the Peel-Away introducer.
3. Ensure that the introducer is within the stomach by visualizing, fluoroscopically, the tip of the introducer beginning to curve in the fundus of the stomach. Remove the dilator and guiding catheter, leaving the sheath portion of introducer and wire guide in place.

Placement of Gastrojejunostomy Catheter

1. Advance the gastrojejunostomy catheter over the wire guide and through the sheath. **NOTE:** Rotating the catheter with a pulsing motion will help facilitate passage of the catheter through the pylorus and duodenum and into the jejunum.
2. To position the Friction-Lock Malecot portion of the catheter totally within the stomach, it may be necessary to peel away the external portion of the sheath. **It is essential that the Friction-Lock Malecot be entirely within the stomach prior to removal of the Peel-Away sheath.**
3. After the catheter is fully inserted, peel away the remainder of the sheath.
4. While holding the external portion of the catheter with one hand, just proximal to the white suture wing, grasp the suture wing and slide it forward until the black indicator mark on the external portion of the catheter is visible. This will indicate proper expansion of the Friction-Lock Malecot. Gentle traction on the catheter should demonstrate the Friction-Lock Malecot to be within the stomach and against the anterior gastric wall.
5. Remove the wire guide and inject contrast to identify the catheter tip position in the jejunum.
6. To prevent inadvertent straightening of the catheter, tie a suture completely around the catheter and through the medial holes on the butterfly wing. **NOTE:** In most cases, it is not necessary to suture the butterfly wing to the anterior abdominal wall, but if this is necessary, it can be accomplished via the lateral holes.
7. Dress the site in a routine manner.
8. Place a male Luer lock cap on the hub(s) of all catheter lumens or extensions. **NOTE:** If using a double-lumen gastrojejunostomy catheter, close the clamps at this time.

Unlocking Friction-Lock Malecot

1. Release suture, if present, from the suture wing.
2. Holding the suture wing, advance the inner coaxial catheter 2-3 cm, if possible, to relax the Malecot.

Post-Placement Care

- It is not uncommon for a small amount of intraperitoneal air to be present post-procedure.
- Feedings may be instituted immediately, although it is recommended to begin with dilute feedings and advance gradually. **NOTE:** If signs of peritonitis develop, feedings should be discontinued and appropriate antibiotics instituted. Feedings should not be resumed until all parameters are normal.
- Suture anchor(s) may be left in place for up to two weeks while tract formation occurs. Cutting the suture end(s) at the skin level releases the anchor(s) into the organ, allowing passage via the gastrointestinal system.
- If using a Luer tip syringe for flushing purposes, remove feeding adapter and connect Luer converter to the hub of the jejunal catheter. Remove Luer converter and connect the feeding adapter after flushing the catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open pouches. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.