COOK

MEDICAL

CE
0088

EN
3

Cope Gastrojejunostomy Cannula

Instructions for Use

DA
3

Cope gastrojejunostomikanyle

Brugsanvisning

DE
4

Cope Gastrojejunostomie-Kanüle

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EL
5

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T _ G R M _ R E V O

COPE GASTROJEJUNOSTOMY CANNULA

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Cope Gastrojejunostomy Cannula is a 13 gage diameter, 30 cm long stainless steel cannula with a malleable distal tip. Its inside diameter accepts wire guides and catheters up to 5.5 French in diameter.

INTENDED USE

The Cope Gastrojejunostomy Cannula is intended for directional control and support of wire guides and catheters through the pylorus during percutaneous gastrojejunostomy procedures.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- Fluoroscopic control is required for use of this product.
- Use caution during cannula advancement to avoid damaging the gastric mucosa or perforating the stomach wall.
- Always have the catheter tip leading the cannula tip during cannula advancement.

INSTRUCTIONS FOR USE

NOTE: Prior to use, the distal tip of the cannula may be shaped into a slight curve to aid in positioning the cannula tip at the pyloric canal.

1. After gaining percutaneous access into the gastric lumen, insert a Peel-Away introducer sheath appropriately sized to the jejunal feeding tube to be placed.
2. Insert a 5 French multipurpose catheter through the lumen of the cannula until the catheter tip protrudes slightly beyond the cannula tip.
3. After opacifying the stomach with dilute contrast, advance the cannula and 5 French catheter as a unit through the introducer sheath. Verify the catheter tip extends beyond the cannula tip throughout the procedure.
4. Under fluoroscopic monitoring, use a gentle, circular probing motion to direct the cannula tip into the pyloric canal.
5. Holding the cannula in position, insert a J-tipped wire guide through the catheter, and advance the catheter and wire guide through the pylorus and duodenum to the jejunum.
6. Remove the J-tipped wire guide and insert an exchange length wire guide.
7. Maintaining the wire guide in position, remove the cannula and catheter.
8. Insert the jejunal feeding tube over the wire guide and advance it to its desired location.
9. Remove the wire guide, flush the jejunal catheter, and cap the tube or connect to an enteric feeding set.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2012
Bloomington, IN 47404 U.S.A.