- EN**
3 **Greene Renal Transplant Stent Set**
Instructions for Use
- CS**
5 **Souprava stentu Greene pro transplantaci ledvin**
Návod k použití
- DA**
7 **Greene stentsæt til nyretransplantation**
Brugsanvisning
- DE**
9 **Greene Nierentransplantat-Stentset**
Gebrauchsanweisung
- EL**
11 **Σετ ενδοπρόσθεσης νεφρικής μεταμόσχευσης Greene**
Οδηγίες χρήσης
- ES**
13 **Equipo de stent para trasplante renal Greene**
Instrucciones de uso
- FR**
15 **Set d'endoprothèse de Greene pour transplantation rénale**
Mode d'emploi
- HU**
18 **Greene vesetranszplantációs szentkészlet**
Használati utasítás
- IT**
20 **Set standard con stent per trapianto renale Greene**
Istruzioni per l'uso
- NL**
22 **Greene renale transplantatie-stentset**
Gebruiksaanwijzing
- NO**
24 **Greene stentsett for nyretransplantasjon**
Bruksanvisning
- PL**
26 **Zestaw stentu Greene do transplantacji nerki**
Instrukcja użycia
- PT**
28 **Conjunto de stent para transplante renal Greene**
Instruções de utilização
- SV**
31 **Greene njurtransplantationsstentset**
Bruksanvisning

GREENE RENAL TRANSPLANT STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

Set Contains:

- Greene Renal Transplant Stent
- Stent positioner
- Appropriate diameter hydrophilic wire guide
- Open End Ureteral Catheter

INTENDED USE

Greene Renal Transplant Stent Set is used to establish temporary internal drainage from the ureteropelvic junction to the bladder. It is indicated for use following renal transplant.

CONTRAINDICATIONS

None known

WARNINGS

Formation of knots in multi-length stents may occur. This may result in injury to the ureter during removal and/or the need for additional surgical intervention. The presence of a knot should be considered if significant resistance is encountered during attempts at removal.

PRECAUTIONS

- Maximum indwelling time is less than or equal to 3 months. If the patient's status permits, the stent may be replaced with a new stent.
- Complications of ureteral stent placement are documented in literature. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures.
- Do not forcefully advance any component during placement or removal of the catheter. If resistance is encountered, stop. Determine the cause of the resistance before proceeding.
- These stents are not intended as permanent indwelling devices.
- A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the wire guide or stent should be avoided.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for the reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

- Extravasation
- Occlusion
- Migration
- Hemorrhage
- Sepsis
- Perforation of the urinary tract

- Peritonitis
- Encrustation
- Urinary tract infection
- Loss of renal function

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Greene Renal Transplant Stent is **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 Tesla
- Maximum magnetic field spatial gradient of 1900 gauss/cm (19 T/m) or less

Under the scan conditions provided above, the Greene Renal Transplant Stent is not expected to experience RF induced heating above background heating.

The image artifact extends approximately 2.5 mm from the Greene Renal Transplant Stent as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

INSTRUCTIONS FOR USE

Wire Guide Preparation

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. To activate the hydrophilic coating on the wire guide, fill a syringe with sterile water or sterile saline solution and attach it to the flushing port on the wire guide holder. Inject enough solution to wet the wire guide surface entirely. **NOTE:** For optimal performance, rehydrate the hydrophilically coated wire guide after exposure to ambient environment.
3. Remove the wire guide from the holder.

Stent Placement

1. Pass the flexible wire guide tip to the renal pelvis. Tortuosity in the obstructed ureter often can be resolved using the wire guide and the open-end ureteral catheter in combination. Remove catheter before attempting to place stent.
2. Pass the stent over the wire guide through the cystoscope. Under direct vision, advance the stent into the ureter with the stent positioner. Have an assistant hold the wire guide in position to prevent advancement of the wire guide into the renal parenchyma.
3. Watch for the distal end of the stent at the ureterovesical junction. The non-radiopaque ink bands serve as a visual indicator regarding where the coils will form as the stent is placed. At that point, halt advancement of the stent.
4. As an assistant removes the wire guide, hold the stent in position with the positioner. The stent pigtail will form spontaneously. Carefully remove the positioner from the cystoscope.

NOTE: If necessary, final adjustment can be made with endoscopic forceps. The stent may be removed by gentle withdrawal traction using endoscopic forceps.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd