**Cook-Z® Stents: Gianturco-Rösch
Tracheobronchial Design**

Instructions for Use

**EN
6**

**Cook-Z® stents: Gianturco-Rösch
trakeobronkial design**

Brugsanvisning

**DA
8**

**Cook-Z® Stents: Tracheobronchiales
Gianturco-Rösch-Design**

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**DE
10**

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
Fig. 7

COOK-Z® STENTS: GIANTURCO-RÖSCH TRACHEOBRONCHIAL DESIGN

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Cook-Z Stent consists of a stainless steel self-expanding double-body z-stent and a Peel-Away® stent loader assembly. The stent is introduced through a Cook-Z Stent Introducer Set, which consists of a .038 inch Amplatz Extra Stiff Wire Guide, radiopaque introducer, and radiopaque stent positioner.

NOTE: Enclosed with this product is a tracking form for device implant registration. Complete this form and return it in the postage-paid envelope provided. Complete a form for each product opened, whether it is used or discarded.

INDICATIONS FOR USE

The stent is intended for use in patients with extrinsic narrowing of the trachea and/or bronchus as a result of tumor encasement or compression in which other interventional techniques did not give satisfactory results. The stent has been designed to maintain patency of stenosed trachea or bronchus due to extrinsic compression caused by a tumor. **It is indicated as a palliative measure, particularly for patients with end-stage malignant airway obstruction. It is not intended for intravascular use.**

The device is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of tracheobronchial stents should be employed.

CONTRAINDICATIONS

- Patients who will require repeated endotracheal laser treatments
- Patients taking daily or alternate daily doses of steroid medications

WARNINGS

- This device is not intended for removal.
- **The safety and effectiveness of this device for use in the vascular system has not been established and can result in serious harm and/or death.**
- **The use of metallic tracheal stents may preclude the success of subsequent surgical procedures and should be considered only for patients after all alternative therapies have been exhausted. Please consider the use of tracheal surgical procedures or placement of silicone stents before using metallic tracheal stents.**
- Placement of this device requires strict attention to technical details. Review the following information carefully prior to proceeding with stent placement.

PRECAUTIONS

- The stent is packaged pre-loaded in the Peel-Away stent loader assembly. Do not remove stent from assembly; reloading will be difficult.
- The stent is not designed to inhibit tumor growth. Malignant tissue may grow through stent and occlude tracheobronchial passages.
- Upon successful deployment of the stent, use a soft, rubber catheter for routine endotracheal suctioning. **Use of a stiff cannula for endotracheal suctioning may lead to the stent becoming caught in the suctioning apparatus.**
- **Do not attempt to balloon-expand the stent after deployment.** Rupture of the balloon may lead to migration of the device during removal.
- The correct size Cook-Z Stent should be selected according to individual tracheal and/or bronchial lumen diameters. Use of a too-small stent may lead to stent migration.

PRODUCT RECOMMENDATIONS

Stent Sizing

- Obtain a radiograph of the chest.
- Measure the diameter of the normal tracheobronchial passage immediately adjacent to the area intended for stenting.
- The labeled stent diameter should be approximately 15-25% larger than the estimated normal lumen diameter of the tracheobronchial tree.

Multiple Stents

If multiple stent placements are required for a patient, the following recommendations should be considered:

- The most distal area of narrowing should be stented first, followed by more proximal locations. **Positioning a stent distal to a stent already in place can be difficult. Reintroduction of the introducer sheath through the stent may displace it.**
- Stents placed in line should touch one another or slightly overlap.

- Stents placed within multiple bronchi should touch one another at the junction point.
- When stenting multiple bronchi that come together at the same point, a wire guide should be positioned across each narrowed bronchus to be stented prior to placing the first stent. This will prevent complete closure of other bronchi as each stent is deployed.

INSTRUCTIONS FOR USE

Illustrations

NOTE: Whenever possible and appropriate, general anesthesia is recommended before insertion of tracheobronchial stents.

1. Localize the stricture via chest x-ray, tomography, CT, etc.
2. Delineate the stricture by placing the end of a flexible bronchoscope at the exact level of the stricture. **NOTE:** The procedure should always be performed under fluoroscopic control. However, the use of a flexible bronchoscope will assist in accurate delineation and placement.
3. Insert the Amplatz Extra Stiff Wire Guide through the bronchoscope's operating channel and advance the wire guide beyond the stricture.
4. Remove the bronchoscope, leaving the wire guide in place.
5. Advance a dilatation balloon over the wire guide, position it across the area to be stented, and dilate the stricture. Mark and measure the area to determine the diameter and length of stent(s) required.
6. Advance the introducer sheath assembly over the wire guide and position it so that the distal marker of the sheath lies beyond the stricture. (Fig. 1)

CAUTION: If an endotracheal tube is in place, ensure that its inner diameter (ID) is at least 8 mm prior to introduction of stent introducer.

7. Remove the dilator, leaving the wire guide and stent introducer sheath in place.
8. Advance the metal cannula and stent, pre-loaded within the Peel-Away stent loader, over the wire guide (non-cuffed end first). (Fig. 2)
9. Once the metal cannula and Peel-Away stent loader have been advanced over the wire guide, remove the metal cannula. **NOTE:** Do not remove the metal cannula until the wire guide has been inserted through the entire Peel-Away stent loader assembly. Otherwise, the wire guide may exit through the struts of the stent.
10. Advance the Peel-Away stent loader assembly **firmly** into the proximal fitting of the stent introducer sheath until it stops. (Fig. 3)
11. Advance the stent positioner over the wire guide and up to the Peel-Away stent loader. Maintaining the Peel-Away stent loader assembly **firmly** against the end of the introducer sheath, advance the stent positioner through the Peel-Away stent loader, pushing the stent inside the stent introducer sheath. (Fig. 4)
12. Advance the stent positioner until the hub of the positioner meets the cuffed end of the Peel-Away stent loader. This will position the stent at the distal tip of the stent introducer.
13. Peel away the stent loader by grasping the knobs and pulling outward. (Fig. 5)

CAUTION: Do not advance the stent positioner. Doing so will result in inadvertent deployment of the stent.

14. Carefully verify that the position of the stent within the introducer sheath is centered across the area to be stented prior to stent deployment.
15. Maintaining position of the stent positioner, under fluoroscopic control, **slowly pull back** the stent introducer sheath with radiopaque marker to deploy the stent across the strictured area. (Fig. 6)

CAUTION: Do not push stent through the introducer sheath to deploy.

16. Perform a post-placement chest radiograph to confirm stent position and tracheobronchial patency. (Fig. 7) **Long-term patency with the stent has not been established; therefore, periodic evaluation is advised.**

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.