- EN**  
3 **H-30 Holmium Laser Fiber (Multi-Use)**  
Instructions for Use
- CS**  
8 **Holmiový laserový optický kabel H-30 (pro více použití)**  
Návod k použití
- DA**  
13 **H-30 Holmium laserfiber (flergangsbrug)**  
Brugsanvisning
- DE**  
18 **H-30 Holmium-Laserfaser (Mehrfachgebrauch)**  
Gebrauchsanweisung
- EL**  
23 **Ίνα λείζερ ολμίου H-30 (πολλών χρήσεων)**  
Οδηγίες χρήσης
- ES**  
29 **Fibra de láser de holmio H-30 (multiuso)**  
Instrucciones de uso
- FR**  
35 **Fibre laser holmium H-30 (usage multiple)**  
Mode d'emploi
- HU**  
40 **H-30 holmiumlézerszál (többször használatos)**  
Használati utasítás
- IT**  
45 **Fibra laser a olmio H-30 (riutilizzabile)**  
Istruzioni per l'uso
- NL**  
51 **H-30 holmiumlaserfiber (meermalig gebruik)**  
Gebruiksaanwijzing
- NO**  
56 **H-30 holmiumlaserfiber (flerbruks)**  
Bruksanvisning
- PL**  
61 **Światłowód lasera holmowego H-30 (do wielorazowego użytku)**  
Instrukcja użycia
- PT**  
67 **Fibra de laser de hólmio H-30 (várias utilizações)**  
Instruções de utilização
- SV**  
72 **H-30 holmiumlaserfiber (för flergångsbruk)**  
Bruksanvisning



Fig. 1

From laser  
Od laseru  
Fra laser  
Vom Laser  
Από το λέιζερ  
Del láser  
Vient du laser  
A lézer felől  
Dal laser  
Van laser  
Fra laser  
Od lasera  
Do laser  
Från lasern

→ To surgical site  
K místu chirurgického zákroku  
Til operationssted  
Zur Operationsstelle  
Προς το σημείο της χειρουργικής  
επέμβασης  
Al sitio quirúrgico  
Vers le site chirurgical  
A műtéti hely felé  
Al sito chirurgico  
Naar operatieplaats  
Til kirurgsted  
Do pola operacyjnego  
Até ao local da cirurgia  
Till operationsstället

Fig. 2

## H-30 HOLMIUM LASER FIBER (MULTI-USE)

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### INTENDED USE

Intended for stone fragmentation, incision, excision, ablation, and coagulation (hemostasis) when attached to the H-30 Holmium Laser System or the Odyssey Laser System for the indications for which the lasers have been cleared.

### CONTRAINDICATIONS

This device is contraindicated for treatment of patients for whom endoscopic procedures are not recommended. Refer to the laser system manual for further contraindications.

### WARNINGS

- All personnel present in the laser hazard zone must wear all suggested protective devices. Wear the laser safety eyewear that the laser manufacturer specifies.
- The holmium laser fiber is validated only for the cleaning procedure, steam sterilization parameters, and STERRAD® cycles that are listed below and must be reprocessed by use of one of the methods listed below.
- Baskets, wire guides and other ureteroscopic accessories may be damaged by direct contact with the laser treatment beam.
- Do not bend fiber at sharp angles. If visible light (aiming beam) can be seen leaking from the fiber, fiber failure may result when therapeutic energy is applied as the fiber is deflected beyond the optical limits of total internal reflection. Discard fiber if visible light leakage is observed.
- Fiber should not be clamped with forceps or other securing instruments as it could result in fiber damage or breakage.
- The device must be thoroughly cleaned and sterilized per the instructions provided in this document before re-use. Ensure the device and all its crevices are completely dry before use.
- Any deviation(s) to the reprocessing instructions prescribed in this IFU may damage the device and result in a device which has compromised sterility or is not fully reprocessed to relevant standards and guidelines.
- The reprocessing instructions provided below have been validated as being capable of preparing the laser fiber for re-use. It remains the responsibility of the processor to ensure that the processing as actually performed using equipment, materials and personnel in the processing facility achieves the desired result. This requires validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided below should be properly evaluated for effectiveness and potential adverse consequences.

### LIMITATIONS ON REPROCESSING

The laser fiber can be reprocessed up to 20 times, which may not reflect the actual number of uses. The device has been subjected to 20 autoclave and Sterrad cycles and has shown very limited wear with no degradation in power transmission performance or mechanical integrity.

### PRECAUTIONS

- To avoid reducing the life of the fiber, follow these precautions:
  - Do not improperly handle the fiber.
  - Do not lase at high power for extended periods of time.
  - Do not lase continuously with the fiber tip in contact with the tissue.
  - Do not lase with a fiber that has a contaminated or damaged proximal end.
  - Do not operate a laser with poor beam alignment or focus.
  - Do not subject the fiber to sharp bends during handling, use, or storage.
- Always keep the proximal connector end dry and free of contaminants.
- Discard any laser fiber that is cracked or broken or does not meet minimum power transmission standards.
- Do not use this laser fiber in the presence of flammable anesthetics or combustible materials.
- Ferrule dust cap should stay attached when fiber is not connected to the laser system.
- Do not exceed recommended power limits.

Fiber Core Size	Power Limit (Single Cavity Lasers)
273 $\mu\text{m}$	15 W
365 $\mu\text{m}$	20 W
550 $\mu\text{m}$	30 W
940 $\mu\text{m}$	30 W

**NOTE:** In addition to the above recommended fiber-handling protocol and precautions, all other applicable hospital protocols that are required to ensure patient and hospital personnel safety should be followed.

### PRODUCT RECOMMENDATIONS

Always wear safety eyewear and medical grade gloves when reprocessing laser fibers.

### POTENTIAL ADVERSE EVENTS

The potential adverse events associated with holmium laser fibers may include but are not limited to:

- Infection
- Thermal damage
- Hematoma
- Perforation
- Discomfort
- Post-procedure fever and leukocytosis (associated with tissue destruction)
- Bleeding
- Hypertension
- Vasovagal response
- Edema
- Delay in healing

### INSTRUCTIONS FOR USE

1. Open and remove the tray from the sterile pouch.
2. Remove the fiber clip from the tray.
3. Peel and remove the backing from the fiber clip. (**Fig. 1**)
4. Adhere the fiber clip to the desired location.
5. Remove the proximal connector from the tray.

**NOTE:** Keep the output of tubing in line with the fiber to minimize drag and prevent kinking as the fiber exits the tubing.

6. Hand the proximal connector to a nonsterile attendant while maintaining sterility of distal fiber.

**NOTE:** The laser system is NOT sterile. Anyone who touches the laser system to insert the fiber will no longer be sterile.

7. Take off the ferrule dust cap and insert the proximal connector into the laser system. Screw in the proximal connector until it is finger-tight.
8. Slide the fiber into the fiber clip's lower clip to hold the fiber in place.

**NOTE:** Use the upper clip to hold the fiber tip when the fiber is not in use. (**Fig. 2**)

9. Check the integrity of laser fiber by using visible aiming beam from laser system. Observe visible light profile out of laser fiber tip. The profile should be a well-defined round or oval pattern. If output beam is not round/oval pattern, or visible light leaks from any spot along the length of the fiber shaft, then the fiber core is broken and the fiber should be discarded as hazardous waste.
10. Follow the laser manufacturer's directions for use.

11. After use, gently unscrew the proximal connector from the laser system and immediately replace the ferrule dust cap.

**NOTE:** Take care to avoid contamination of the proximal connector.

12. Discard the fiber clip.

## **FIBER REPROCESSING INSTRUCTIONS**

### **Initial Treatment at Point of Use**

1. Apply two sprays of 70/30 isopropyl alcohol (IPA) onto a new lint-free wipe and wipe the length of the laser fiber from the proximal to the distal end to remove residual organic matter, blood, and irrigation solution.

**NOTE:** Any debris should not be allowed to dry on the fiber or stay on the fiber for more than 1 hour after use. If debris has dried, then cut the soiled portion of the fiber with a cutting device.

2. Repeat the previous step, using a new lint-free wipe with 70/30 IPA to wipe the laser fiber a second time.
3. Coil the fiber to a diameter that is no less than 8 inches (20 cm) before transporting it between various departments.
4. Ensure the ferrule dust cap is securely fastened to the proximal connector.
5. Follow facility protocols on handling and transport of contaminated devices.
6. It is recommended that instruments be reprocessed within 1 hour following use. If the device is allowed to dry for more than an hour and the contaminants have dried, then cut the soiled portion of the fiber with a cutting device.

### **Preparation Before Cleaning: Cutting, Stripping, and Cleaving**

#### **Proximal Face Inspection**

1. Having removed the ferrule dust cap, use a microscope with 20X magnification to inspect the proximal face of the fiber for damage. If visible damage such as debris, scratches, chips, or pits is observed, then the fiber may damage the laser system and should be discarded as hazardous waste.
2. Put the ferrule dust cap on the connector end immediately after removing the fiber from the microscope. This ferrule dust cap should stay attached throughout the rest of the fiber reprocessing.

#### **Cutting:**

3. Using a cutting device, cut 2.5 cm of the distal tip of the fiber. Discard the cut portion of the fiber as hazardous waste.

**NOTE:** Do not use a cleaver when cutting is specified. Scissors are acceptable.

**NOTE:** Be sure to cut any length of fiber that has been flexed within a flexible endoscope and might have received significant stress from bending.

**NOTE:** Be sure to use a cutting device to remove any length of fiber upon which residual organic material has adhered and from which it has not been removed.

#### **Stripping:**

4. Hold the appropriately sized Cook Fiber Stripping Tool in one hand and the fiber in the other.
5. Insert the distal end of the fiber into the fitting hole of the Fiber Stripping Tool until the fiber tip reaches the first indicated marking on the side body of the Fiber Stripping Tool. The marking length is equivalent to approximately 1 cm of fiber length.
6. Squeeze both buttons on the Fiber Stripping Tool to close the blades, and keep the tool in this closed position while simultaneously pulling the fiber away from the tool. Discard the stripped portion of the fiber appropriately, as hazardous waste.

#### **Cleaving:**

7. Using the cleaver, score the stripped portion of the fiber as close to the blue jacket as possible. Make a single score line that is as perpendicular to the fiber axis as possible. Push down gently on the fiber with the cleaver, but do not attempt to cut all the way through the fiber with the cleaver, as this will result in a poor fiber surface.
8. Grasp the fiber on opposite sides of the score mark. Pull the fiber on opposite sides of the score mark in a straight, opposing direction until the fiber separates. Do not bend to separate the fiber, because bending will result in a poor surface. Discard the cleaved distal end of the fiber.
9. Using the same method as above in steps 4-6, strip the buffer to expose 4-6 mm (0.4-0.6 cm) of the inner glass fiber.

### **Fiber Integrity and Length Inspection:**

10. Coil the fiber to a diameter that is no less than 8 inches (20 cm), and check the integrity of the length of the fiber shaft as well as the distal surface by launching a visible light source down the fiber ferrule (the green aiming beam from the laser system can be used as the visible light source). If visible light leaks from any spot along the length of the fiber shaft, then the fiber core is broken and the fiber should be discarded as hazardous waste.
11. Check the visible beam exit profile out of fiber tip by projecting it onto a flat, white surface.
  - a. A well-cleaved surface is indicated by a well-defined circular or slightly oval pattern.
  - b. A poorly-cleaved surface is indicated by a pattern with one or more arrays or "comet tails" extending from the central circle. Repeat the cutting, stripping, and cleaving process from the **Preparation Before Cleaning: Cutting, Stripping, and Cleaving** section.

**NOTE:** This step will also help verify that the fiber's integrity is still intact and that no cracks or other damage exist.
12. To ensure that the fiber length is adequate for clinical use, verify that the fiber's length is no less than 165 cm. Discard the fiber as hazardous waste if the fiber length is less than 165 cm.

### **Cleaning Instructions**

#### **CAUTION: Automated cleaning methods are not compatible with this device.**

1. Apply two sprays of 70/30 IPA onto a new lint-free wipe and wipe the length of the laser fiber from the proximal to the distal end to remove residual organic matter, blood, and irrigation solution.
2. Repeat the previous step three more times, using a new lint-free wipe with 70/30 IPA each time.
3. Visually inspect the device in adequate lighting to ensure that it is free of residual organic matter, blood, and other debris.

**NOTE:** If the device still has residual material on its shaft, then repeat the above cleaning procedures or cut the affected length of the fiber and repeat the steps in the **Preparation Before Cleaning: Cutting, Stripping, and Cleaving** and **Cleaning Instructions** sections. If no residual materials are noticeable on the device, then proceed to the **Sterilization** section.

### **Drying**

1. No post-cleaning drying of the device is required.

### **Inspection and Maintenance**

1. Visually inspect the device in adequate lighting to ensure that it is free of residual organic matter, blood, and other debris.

**NOTE:** If the device still has residual material on its shaft, then repeat the above **Cleaning Instructions** or cut the affected length of the fiber and repeat the steps in the **Preparation Before Cleaning: Cutting, Stripping, and Cleaving** and **Cleaning** sections. If no residual materials are noticeable on the device, then proceed to the **Sterilization** section.
2. The laser fiber does not require any maintenance.

### **Packaging**

1. For steam sterilization, coil the fiber to a diameter that is no less than 8 inches (20 cm) and either wrap the device in a sterilization wrap or place it in an appropriately sized sterilization pouch. Package the device following the sterilization wrap or pouch manufacturer's instructions. Do not exceed more than 2 layers of wraps or pouches.
2. For sterilization using Sterrad systems, coil the fiber to a diameter that is no less than 8 inches (20 cm) and place the device in an appropriately sized Sterrad compatible pouch following the pouch manufacturer's instructions.
3. The sterile barrier system should comply with AAMI/ISO 11607 and AAMI/ISO TIR16775.

### **Sterilization**

**NOTE:** Ensure that the device is completely dry before sterilization by looking for moisture deposition on the surfaces of the device.

1. Sterilization should be performed in accordance with industry standards and with a validated sterilization cycle.
2. Store and handle reprocessed fiber assemblies with the same care as similar new devices, following all previously mentioned warnings, notes, and precautions.
3. Place the appropriately packaged device (see **Packaging** section above) into an ISO/AAMI 17665 compliant moist heat sterilization system (steam sterilizer) or a validated Sterrad system.
4. Select any one of the sterilization cycles indicated below:

## Steam Sterilization

Steam Sterilization (Pre-Vacuum)		
Temperature	132 °C	134 °C
Time	4 minutes	3 - 18 minutes
Drying Cycle	20 minutes	20 minutes

**Caution:** The holmium laser fiber is validated only for the above-listed steam sterilization parameters and must be resterilized according to one of the listed parameters. A maximum of 20 steam sterilization cycles are validated. Elevated temperatures and time will degrade product integrity.

## STERRAD® Sterilization

ASP STERRAD® cycles	
STERRAD® 100S	Short cycle
STERRAD® 100NX	Standard cycle
STERRAD® 100NX	Express cycle
STERRAD® 100NX	Flex cycle
STERRAD® NX	Standard cycle
STERRAD® NX	Advanced cycle

**Caution:** The holmium laser fiber is validated only for the above-mentioned Sterrad system cycles and must be resterilized with one of the Sterrad system cycles listed above. A maximum of 20 Sterrad cycles are validated.

## Storage

The device can be stored in the sterilization pouch/wrap without losing sterility for a duration as prescribed by the pouch/wrap manufacturer. Refer to the sterilization pouch/wrap specifications for further information.

## Additional Information

All the reprocessing instructions detailed in this IFU have been validated for effectiveness by testing each process under worst-case settings.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Reusable. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2022
Bloomington, IN 47404 U.S.A.