EN 3 **Cook Turbo-Flo® HD Acute Hemodialysis Catheter**

Instructions for Use

DA 7 **Cook Turbo-Flo® HD akut hæmodialysekateter**

Brugsanvisning

DE 11 **Cook Turbo-Flo® HD Akut-Hämodialyse-Katheter**

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

1

2

COOK TURBO-FLO® HD ACUTE HEMODIALYSIS CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Cook Turbo-Flo HD Acute Hemodialysis Catheter is a double lumen catheter constructed of polyurethane and incorporating two parallel D-shaped lumens. To minimize blood recirculation, the sideports of the arterial lumen are proximal to the sideport and endhole of the venous lumen. The catheter also contains a swivel fitting, which allows free catheter rotation without removal of securement sutures. This device does not contain natural rubber latex.

The proximal lumen is indicated in red and is used for arterial withdrawal. The distal lumen is indicated in blue and is used for providing venous return.

INTENDED USE

The Cook Turbo-Flo HD Acute Hemodialysis Catheter is intended for acute hemodialysis, apheresis and hemofiltration via percutaneous insertion into the subclavian, jugular, or femoral veins. The device is a short-term use catheter (less than 30 days).

CONTRAINDICATIONS

- Infection at the access site
- Venous stenosis severe enough to prohibit percutaneous entry into the affected vessel

WARNINGS

- Do not use if the packaging or product has been compromised by physical damage.
- Failure to remove the inner obturator after catheter introduction may cause vascular injury and reduced performance.
- Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system. Tip position should be verified by x-ray and monitored on a routine basis. Periodic lateral view x-ray is suggested to assess tip location in relation to the vessel wall. Tip position should appear to be parallel to the vessel wall.
- To avoid vascular injury, do not use excessive force when advancing dilators. Use the smallest size dilator catheter placement will allow. The wire guide must always lead the dilator by several centimeters. Do not advance the dilator more than a few centimeters into the vessel.
- To distend great vessels and to prevent inadvertent air aspiration during subclavian and jugular catheter insertion, the patient should be placed in Trendelenburg position.

PRECAUTIONS

- These products are intended for use by physicians trained and experienced in hemodialysis. Standard techniques for placement of vascular access catheters and wire guides should be employed.
- Strict adherence to aseptic technique is required during placement of this device.
- Any signs of infection, redness or swelling at the access site should be treated with appropriate therapy immediately.
- The catheter should not be used for long-term indwelling applications.
- Do not re-sterilize the catheter.
- Do not cut, trim or modify the catheter or components prior to placement or intraoperatively.
- Patient movement can cause catheter tip displacement. Use should be limited to controlled hospital situations. Catheters placed via either a jugular or subclavian vein have demonstrated forward tip movement of 1-3 cm with neck and shoulder motion.
- If lumen flow is impeded, do not force injection or withdrawal of fluids. Notify attending physician immediately.
- Select appropriate puncture site and length of catheter by assessing patient anatomy and condition.
- Use of electrocardiography (ECG), ultrasound and/or fluoroscopy is suggested for accurate catheter placement.
- Trays that include lidocaine also include a lidocaine label, which should be used to label the syringe for lidocaine injection. Any unused lidocaine should be disposed of at the conclusion of the procedure.
- Prior to insertion, the catheter shaft should not be wiped with or immersed in ethyl alcohol, isopropyl alcohol, or other alcohols, acetone or other non-polar solvents.

POTENTIAL ADVERSE EVENTS

- Infection
- Hematoma
- Hemorrhage
- Catheter thrombosis
- Central venous stenosis

- Central venous thrombosis
- Hemothorax
- Vascular thrombosis
- Bleeding
- Pneumothorax
- Fever
- Perforation of vessel
- Stenosis
- Air embolism
- Bacteremia
- Cardiac arrhythmia
- Cardiac tamponade
- Endocarditis
- Vessel trauma

PRODUCT RECOMMENDATIONS

Recommended Lumen Volumes

Catheter Size	Catheter Length	Lumen Volumes	
		Proximal	Distal
12.0 Fr Double Lumen	15 cm	1.3 mL	1.4 mL
	20 cm	1.4 mL	1.6 mL
	25 cm	1.6 mL	1.7 mL

Flow Rate vs Pressure

Size	Lumen	Flow Rate (mL/min)	Average Maximum Pressure (mmHg) FORWARD FLOW	Average Maximum Pressure (mmHg) REVERSE FLOW
12 Fr 25 cm	Arterial	100 ± 5	-43.05	35.2
		200 ± 5	-90.8	73.6
		300 ± 5	-151.8	121.0
		400 ± 5	-224.1	180.1
	Venous	100 ± 5	41.2	-49.05
		200 ± 5	87.3	-103.23
		300 ± 5	141.7	-171.8
		400 ± 5	211.1	-252.5

Recirculation Rates

Size	Mean Recirculation (%)
12 Fr 25 cm	1.41

The following variables should be considered when selecting appropriate catheter and length:

- Patient history
- Size and age of patient
- Access site available
- Unusual anatomical variables
- Proposed use and duration of treatment plan

Puncture Site

The angle of the catheter tip to the vessel wall should be checked carefully. Blackshear reviewed the medical literature of catheter perforations, which have confirming x-rays, and found that an incident angle of the catheter to the vessel wall greater than 40 degrees was more likely to perforate.¹

Another critical factor that can cause a catastrophic event is the choice of puncture site. Findings by Tocino and Watanabe indicate that the left subclavian and left jugular veins should be avoided when practical. Eighty percent of the perforations or erosions were found when these vessels were used. In addition, they have observed that the tip curve of a wedged catheter can be detected with lateral view x-ray.²

The above discussion is meant to be a guide for puncture site. As more data become available, other causal factors may become evident, but present information suggests that the left subclavian veins, left jugular veins, and femoral veins should be used only when other sites are not available.

Suggested Catheter Lumen Utilization

#1 Distal exit port (endhole) -- indicated in blue and used for providing venous return


#2 Proximal exit port -- indicated in red and used for arterial withdrawal

Suggested Catheter Maintenance

- Catheter entry site must be prepared and maintained in a manner consistent with standard procedure for hemodialysis catheterization.
- Prior to catheter insertion both distal and proximal lumens should be filled with normal saline solution or heparinized saline solution, depending on institutional protocol.
- After catheter is placed and prior to use, tip position and lumen patency should be confirmed by free aspiration of venous blood. If blood is not freely aspirated, the physician should immediately reevaluate catheter tip position.
- Before using any lumen already locked with heparin, the lumen should be aspirated. After use, the lumen should again be flushed with normal saline to prevent systemic heparinization before reestablishing heparin lock.
- Strict aseptic technique must be adhered to while using and maintaining catheter.

Access Sites of Choice

1. High Internal Jugular
2. External Jugular
3. Low Internal Jugular
4. Supraclavicular
5. Infraclavicular
6. Femoral (not shown)

INSTRUCTIONS FOR USE

1. **For subclavian and jugular insertions:** Place the patient in the Trendelenburg position and turn the patient's head slightly to the side opposite the intended insertion site to expose the insertion site. The Trendelenburg position may prevent air embolism. Ultrasound guidance is recommended if it is available.
For femoral insertions: Place the patient on his/her back and palpate both femoral veins for site selection and assessment. On the same side of the insertion site, flex the knee and abduct the thigh, then place the foot across the opposite leg. Ultrasound guidance may be used.
2. Create a strict, sterile operative field prior to procedure, utilizing maximum barrier protection such as, but not limited to, sterile gown, drapes, cap, gloves, and mask.
3. If desired, perform an ultrasound to check patency and confirm position of the vein to be cannulated.

4. If applicable, remove the Luer-lock cap from each extension.
5. Prepare the catheter for insertion by flushing each of the lumens with normal saline, and clamp the proximal extension.
6. Attach the injection cap to the proximal extension. Leave the distal extension uncapped for wire guide passage.
NOTE: To prevent systemic heparinization, heparin concentrations should be chosen and monitored carefully. Catheter heparinization should be determined by institutional protocol and clinical judgment.
7. Infiltrate the subcutaneous tissues at the catheter insertion site with local anesthetic.
8. Introduce a thinwall percutaneous entry needle into the vessel. **Venous blood should be easily aspirated to confirm position of the needle tip within the vessel.**
9. Position the distal end of the Safe-T-J wire guide into the J-tip straightener of the thumb advancing holder.
10. Pass the straightened wire guide through the needle and advance it 5-10 cm into the vessel.
If resistance is encountered during wire guide insertion, do not force the wire guide. Withdrawal of the wire guide through the needle should be avoided; breakage may result. Remove the wire guide and needle as a unit, then begin again with a new needle and wire guide.
11. While maintaining wire guide position, withdraw the needle.
12. Enlarge the puncture site using a number 11 scalpel blade, if required. If dilation is required, a dilator can be advanced over the wire guide and removed prior to insertion of the hemodialysis catheter.
CAUTION: To prevent vascular injury, do not use excessive force when advancing the dilator. Use the smallest size dilator catheter placement will allow. The wire guide must always lead the dilator by several centimeters. Do not advance the dilator more than a few centimeters into the vessel.
13. Measure the catheter to be used against the patient to determine the approximate length of catheter needed from puncture site to central venous tip position.
14. Introduce the hemodialysis catheter over the wire guide. While maintaining wire guide position, advance the catheter into the vessel with a gentle twisting motion. (**Fig. 1**)
NOTE: Do not advance the catheter tip beyond the distal tip of the wire guide. Always have the wire guide leading during catheter placement. Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, the catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC. Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system.
15. After the catheter is in position, remove the wire guide and inner obturator. (**Fig. 2**) Venous blood should be easily aspirated, and saline should flush easily into the hub. The winged hub can now be sutured into place. If the catheter is not introduced to its full length, additional suture should be carefully placed around the catheter and affixed to the skin at the entry site. This will help prevent backward or forward catheter movement.
NOTE: If resistance is met when attempting to remove the wire guide and inner obturator after the catheter has been placed, the wire guide may be kinked at the catheter tip within the vessel. If this occurs, withdraw the catheter approximately 2-3 centimeters and attempt to remove the wire guide. If resistance is met again, remove the wire guide and the catheter together as a unit. Do not pull the wire guide back forcibly through any component.
16. The lumens should now be flushed with 5-10 cc normal saline prior to use or establishment of heparin lock.
17. Each lumen should be capped and locked with solution based on specific hospital policy or protocol.
Recommended lumen volumes are indicated on the catheter extensions and in the Product Recommendations section of this booklet. If a high concentration of heparin is used, it is important to fill with the correct lumen volume to avoid systemic heparinization.
NOTE: A wire guide that is at least twice as long as the catheter is recommended for catheter exchange.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use only. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store the product in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

1. Blackshear RH, Gravenstein N. "Critical Angle of Incidence for Delayed Vessel Perforation by Central Venous Catheter: A Study of *In Vivo* Data." *Annals of Emergency Medicine*, 1992; 21:242.
2. Tocino IM, Watanabe A. "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition." *American Journal of Roentgenology*, 1986; 146:487-490.

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1. Blakeshear RH, Gravenstein N. "Critical Angle of Incidence for Delayed Vessel Perforation by Central Venous Catheter: A Study of *In Vivo* Data." *Annals of Emergency Medicine*, 1992; 21:242.
2. Tocino IM, Watanabe A. "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition." *American Journal of Roentgenology*, 1986; 146:487-490.

MANUFACTURER
COOK INCORPORATED
750 Daniels Way
Bloomington, IN 47404 U.S.A.

EC REPRESENTATIVE
WILLIAM COOK EUROPE ApS
Sandet 6, DK-4632
Bjæverskov, DENMARK