**COOK**

**MEDICAL**

**CE**  
0088

**EN**  
3

## Heavy Duty Double Flexible Tipped Wire Guide

Instructions for Use

**DA**  
4

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T \_ H D O C \_ R E V O

## HEAVY DUTY DOUBLE FLEXIBLE TIPPED WIRE GUIDE

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### DEVICE DESCRIPTION

The Heavy Duty Double Flexible Tipped Wire Guide is constructed of stainless steel.

Refer to the product labeling for specific model information regarding wire guide diameter and length.

### INTENDED USE

Intended to facilitate placement of devices used during diagnostic and interventional procedures.

### CONTRAINDICATIONS

Not intended for use in the coronary or neurovascular arteries.

### WARNINGS

- Possible allergic reactions should always be considered.
- Wire guide stiffness should be considered.
- If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take appropriate steps to relieve the resistance. Advancement and withdrawal of the wire guide should be performed slowly and carefully.

### PRECAUTIONS

- Manipulation of wire guides requires fluoroscopic visualization.
- End hole size and length of the catheter must be taken into consideration to ensure proper fit between wire guide and catheter.
- Altering tip configuration of the wire guide may result in wire or vessel damage.
- Withdrawal or manipulation of the wire guide through a needle tip may result in breakage.
- The products are intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

### POTENTIAL ADVERSE EVENTS

Using this product may cause the following adverse events.

- Arrhythmias
- Death
- Thrombosis and embolization
- Vessel dissection or perforation

## INSTRUCTIONS FOR USE

1. Inspect the wire guide prior to use for tip shape, bends, kinks or coil separation. If the wire guide is damaged, DO NOT USE.
2. Using sterile technique, localize and puncture the target vessel with a needle cannula.
3. Insert the distal end of the wire guide through the needle cannula and into the vessel.
4. Remove the needle cannula, leaving the wire guide within the vessel.
5. If applicable, dilate the vessel and make an appropriate skin incision.
6. Pass the catheter over the wire guide, within the lumen of the vessel, to the desired position. **NOTE:** Distance markers on the wire guide can help identify the depth of the wire in the vessel.
7. Carefully remove the wire guide from the catheter.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2015
Bloomington, IN 47404 U.S.A.