- EN**
5 Hilal Embolization Microcoils
Instructions for Use
- CS**
6 Embolizační mikrospirálky Hilal Embolization Microcoils
Návod k použití
- DA**
8 Hilal Embolization Microcoils
Brugsanvisning
- DE**
10 Hilal Embolisationsspiralen Microcoils
Gebrauchsanweisung
- EL**
11 Microcoils εμβολισμού Hilal
Οδηγίες χρήσης
- ES**
13 Microespirales de embolización Hilal Microcoils
Instrucciones de uso
- FR**
15 Microcoils d'embolisation Hilal
Mode d'emploi
- HU**
17 Hilal Microcoils embolizációs mikrospirálok
Használati utasítás
- IT**
19 Microspirali per embolizzazione Hilal Microcoils
Istruzioni per l'uso
- NL**
21 Hilal Microcoils voor embolisatie
Gebruiksaanwijzing
- NO**
22 Hilal Embolization Microcoils
Bruksanvisning
- PL**
24 Mikrospirale embolizacyjne Hilal Microcoils
Instrukcja użycia
- PT**
26 Microcoils Hilal de embolização
Instruções de utilização
- SV**
28 Hilal Embolization Microcoils
Bruksanvisning

.018 INCH MICROCOILS **0,018 INCH MICROCOILS** **MICROCOILS DE 0,018 INCH** **MICROCOILS VAN 0,018 INCH** **MICROCOILS DE 0,018 INCH**

MIKROSPIRÁLKY 0,018 INCH **MICROCOILS 0,018 INCH** **0,018 INCH MÉRETŰ MIKROSPIRÁLÓK** **MIKROSPIRALER PÅ 0,018 INCH** **0,018 INCH MICROCOILS**

0,018 INCH MICROCOILS **MICROESPIRALES DE 0,018 INCH** **MICROSPIRALI DA 0,018 INCH** **MIKROSPIRALE O ŠREDNICY 0,018 INCH**



1

3

2

4

5

1 Weaving

2 Sheath

3 Catheter

1 Pohyb sem a tam

2 Sheath

3 Katetr

1 Snoning

2 Sheath

3 Kateter

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2 Schleuse

3 Katheter

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3 Kateter

HILAL EMBOLIZATION MICROCOILS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Hilal Embolization Microcoils are made of platinum with spaced synthetic fibers and supplied preloaded in a loading cartridge.

INTENDED USE

Hilal Embolization Microcoils are intended for arterial and venous embolization in the peripheral vasculature.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that single and multiple Hilal Embolization Microcoils are **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned after placement under the following conditions:

- Static magnetic field of 3.0 Tesla or 1.5 Tesla only
- Maximum magnetic field spatial gradient of 1700 gauss/cm (17 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode)

Under the scan conditions provided above, the Hilal Embolization Microcoils are not expected to result in a temperature rise of more than 3.0 °C after 15 minutes of continuous scanning.

The image artifact extends approximately 5.7 mm from the Hilal Embolization Microcoils as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners.

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

CONTRAINDICATIONS

None known

WARNINGS

- Positioning of embolization coils should be done with particular care. Coils should not be left too close to the inlets of arteries and should be intermeshed with previously placed coils if possible. A minimal but sufficient arterial blood flow should remain to hold the coils against the previously placed coils until a solid clot ensures permanent fixation. The purpose of these suggestions is to minimize the possibility of loose coils becoming dislodged and obstructing a normal and essential arterial channel.
- Hilal Embolization Microcoils are not recommended for use with polyurethane catheters or catheters with sideports. If a catheter with sideports is used, the embolization coil may lodge in the sideport or pass inadvertently through it. Use of a polyurethane catheter may also result in lodging of the embolization coil within the catheter.
- If difficulties occur when deploying the embolization coil, withdraw the wire guide, coil and angiographic catheter simultaneously as a unit.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in embolization techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Perform an angiogram prior to embolization to determine correct catheter position.
- Prior to introduction of the embolization coil, flush the angiographic catheter with saline.
- **If using a Hilal Embolization Microcoil, ensure that the delivery catheter has an internal diameter (ID) of .018 to .025 inch.**

PRODUCT RECOMMENDATIONS

The following table offers specific recommendations.

Coil Size Diameter	Catheter Type & Size	Wire Guide Type & Size
.018 inch	MCS-2.5 MCS-2.8	STF-18 TSF-18 TSFB-18

Coil Delivery Technique and Coil Size Selection

Long-term occlusion depends on achieving cross-sectional occlusion of the blood vessel, and coaxial catheters provide the ability to control placement of coils and permanent occlusion.

- Coaxial technique: The use of an outer guiding sheath/catheter is the most important step in preventing coil elongation and uncertain long-term occlusion. The outer guiding sheath/catheter provides support, and the inner catheter provides finer selective maneuvers. (Fig. 5)

In general, the first coil selected should have a diameter that is 20% larger, or at least 2 mm oversized, than the vessel that is being occluded.

INSTRUCTIONS FOR USE

1. Perform an angiogram prior to embolization to determine optimal catheter position.
2. Push the loading cartridge completely into the delivery catheter. (Fig. 1)
3. Advance the Luer lock connector fitting toward the catheter hub and lock into place. (Fig. 2)
4. Use the pusher stylet to load the Hilal Embolization Microcoil into the delivery catheter. **NOTE:** The stylet must be pushed as far as possible into the loading cartridge to ensure proper loading. Remove the pusher stylet and loading cartridge. (Fig. 3)
5. To obtain secure placement of the Hilal Embolization Microcoil, we recommend the use of the pusher stylet to push the Microcoil further into the delivery catheter. (Fig. 4)
6. Once coil is positioned into delivery catheter, verify position of catheter tip prior to deployment. Coil design permits delivery into the target vessel by saline/contrast flush or by push technique using an appropriately sized wire guide or pusher.
7. Perform final angiogram to confirm coil position within target vessel.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2017
Bloomington, IN 47404 U.S.A.