EN
3

Holmium Laser Fiber (Multi-Use)

Instructions for Use

CS
6

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
Fig. 2

HOLMIUM LASER FIBER (MULTI-USE)

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

INTENDED USE

Used with holmium laser systems for the treatment of kidney, ureter, and bladder stones, as well as soft tissue ablation.

CONTRAINDICATIONS

None known

WARNINGS

All personnel present in the laser hazard zone must wear all suggested protective devices. Wear laser safety eyewear per the laser manufacturer's specifications.

PRECAUTIONS

- A number of different factors affect the life of any particular fiber, including:
 - Extended lasing at high power
 - Continuous lasing with the fiber tip in contact with tissue
 - Poor reprocessing (minimum length of fiber removed in reprocessing should be between 2.5 and 3 inches)
 - Lasing with a contaminated or damaged proximal end
 - Improper handling
 - Poor laser beam alignment or focus
- Never subject fiberoptics to sharp bends in handling, use, or storage.
- Always keep connector-end dry and free from contaminants.
- Discard any fiberoptic assembly that is cracked or broken, or does not meet minimum transmission standards.
- Do not exceed recommended power limits.

Fiber Core Size	Power Limit (Single Cavity Lasers)
150 μm	6 W
200 μm	10 W
273 μm	15 W
365 μm	20 W

NOTE: In addition to the above recommended fiber handling protocol and precautions, all other applicable hospital protocols required to ensure patient and hospital personnel safety should be followed.

NOTE: Due to the extremely small fiber core size, the output power from 150 or 200 μm core fibers may be more than 20% lower than set power.

INSTRUCTIONS FOR USE

1. Open and remove the packaging tray from the sterile pouch.
2. Remove the Fiber Clip from the packaging tray.

3. Peel and remove the backing from the Fiber Clip. (**Fig. 1**)
4. Adhere the clip to the desired location.
5. Remove the proximal connector **and tubing** from the packaging tray.
6. Hand the proximal connector to a non-sterile attendant.
NOTE: The laser system is NOT sterile. Anyone who touches the laser system to insert the fiber will no longer be "scrubbed in."
7. Open any laser aperture cover or door and insert the proximal connector into the laser system; screw in finger-tight.
8. Pull the coiled tubing into the sterile field and allow the fiber to slide out of the tubing.
NOTE: Keep output of tubing in line with fiber to minimize drag as the fiber is exiting the coil.
9. Slide the fiber into the lower clip to hold the fiber in place.
NOTE: Use the upper clip to hold the fiber tip when not in use. (**Fig. 2**)
10. Follow the laser manufacturer's directions for use.
11. Remove the fiber from the Fiber Clip and discard.

SUGGESTED FIBER REPROCESSING INSTRUCTIONS

All cleaning/disinfection/sterilization methods listed herein have been validated.

Step One: Cleaning

1. Carefully disconnect the fiber optic delivery system from the laser and immediately replace the ferrule dust cap.
2. Use lint-free wipes to remove all residual organic matter, blood and irrigation solution to allow the sterilization medium to contact the surface of the device.
NOTE: Ethylene oxide gas does not effectively penetrate residual organic matter.
NOTE: Any debris should not be allowed to dry on the fiber. If debris has dried, use a lint-free wipe (pre)saturated with an enzymatic cleaner/detergent for the first wipe(s) or remove the portion of soiled fiber.
3. Using a lint-free towel or gauze and sterile 70/30 IPA, wipe down the length of fiber, focusing on the contaminated distal end and any cracks or crevices where organic soil may collect.
4. Prior to sterilization and reuse, this product should be inspected to verify that it functions properly. All edges and surfaces should be smooth and free of abrasions.
WARNING: Do not immerse the fiber optic assembly.
5. Coil the fiber to a diameter of no less than 6 inches (15 cm).
WARNING: High-level disinfection is not recommended for this device.
WARNING: Automated cleaning methods are not compatible with this device.

Step Two: Inspection

1. Using a microscope with ~20X magnification, inspect the proximal fiber face for damage. If visible damage (such as debris, scratches, chips, or pits) is observed, the fiber may damage the laser system and should be discarded.
2. Replace the dust cap on the connector end immediately after removing from microscope. This dust cap should stay attached throughout the rest of the fiber reprocessing.

Step Three: Cleaving and Stripping

1. Using scissors or other cutting device, remove approximately 2.5-3 inches (6 to 8 cm) of the distal tip of the fiber.
NOTE: Be sure to remove any fiber length that has been flexed within a flexible endoscope and might have received significant bending stresses.
NOTE: Be sure to remove any fiber length upon which residual organic material has adhered and has not been removed with the previous alcohol wipe-down.
2. Hold the appropriate fiber stripping tool in one hand and the fiber in the other.
3. Insert the fiber into the stripping tool until the fiber tip reaches the 1 inch (2.5 cm) rule marking.
4. Squeeze the fiber stripping tool to close the blades and maintain this closed position while simultaneously pulling the fiber from the tool.
5. Using the cleaver, score the fiber as close to the buffer as possible. Make a single score line as perpendicular to

- the fiber axis as possible. Do not attempt to cut the fiber with the cleaver as this will result in a poor surface.
- Grasp the fiber on opposite sides of the score mark. Pull the sides of the score mark in a straight, opposing direction until it separates. Do not bend the fiber while cleaving as this will result in a poor surface. Discard the distal end of fiber.
 - Using the same method as above, strip the buffer to expose 4-6 mm of the bare fiber.
 - Check the integrity of the fiber length as well as the distal surface by launching a visible laser down the fiber while it is not flexed to a diameter of less than 6 inches (15 cm). (The aiming beam from the laser system works well as the visible laser.) Any spot along the length of fiber where the visible laser light is leaking out and readily visible is a break in the fiber core and indicates that the fiber should be discarded.
 - Check the visible laser's exit profile by projecting it onto a flat, white surface.
 - A good cleaved surface is indicated by a well-defined circular or slightly oval pattern.
 - A poor cleaved surface is indicated by a pattern with one or more rays or "comet tails" extending from the central circle. Repeat the Cleaving and Stripping process from step 2.
- NOTE:** This step will also help verify that the fiber integrity is still intact and there are no cracks or other damage. To ensure adequate usability, verify that the fiber length is no less than 2 m.

Step Four: Sterilization

- Sterilization should be performed in accordance with industry standards, with a validated sterilization cycle.
- Store and handle reprocessed fiber assemblies with the same care as similar new devices, following all previously mentioned warnings, notes, and precautions.
- Coil the fiber no tighter than an 8 inch (20 cm) diameter and insert it into a suitable sterilization tray or pouch, or wrap it in sterilization paper in preparation for sterilization.
- Acceptable sterilization methods are: ETO, H₂O₂ gas via Sterrad® Sterilizers, and Steam (Autoclave). Recommended parameters are below.

NOTE: Steam sterilization should be used least often (when time is critical) as it may reduce the number of times the fiber assembly can be reused.

Steam Sterilization (gravity displacement or pre-vacuum)			Sterrad® Sterilization (Models: 100S, 100NX, NX)	
	Min. (pre-vacuum)	Max. (gravity displacement)		
Temperature	132°C	135°C	Hydrogen Peroxide Concentration	7 – 18 mg/L
Time	4 min	10 min	Exposure time	20-30 minutes
Drying Cycle	20 min	30 min	Chamber Temperature	40 – 50°C
			Plasma Power	400 Watts

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Reusable. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404, U.S.A
COOK IRELAND LTD.