**EN** **Roadrunner® UniGlide™ Hydrophilic Wire Guide**

Instructions for Use

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Bruksanvisning

## ROADRUNNER® UNI GLIDE™ HYDROPHILIC WIRE GUIDE

**CAUTION: Federal law (USA) restricts this device to the sale and use by or on the order of a physician.**

### DEVICE DESCRIPTION

Utilizing proprietary processes, the Roadrunner UniGlide Hydrophilic Wire Guide is constructed from a steerable, metallic core with a polymer coating. A hydrophilic coating is applied over the radiopaque polymer jacket. Wire guide Length, Diameter, Coating, Tip, and Core configurations are indicated on the product label.

**STANDARD  
SHAFT**

Standard Shaft Configuration

**STIFF  
SHAFT**

Stiff Shaft Configuration

**ANGLED**

Angled Tip Configuration

**STRAIGHT**

Straight Tip Configuration

### INTENDED USE

The Roadrunner UniGlide Hydrophilic Wire Guide is intended for use in facilitating the delivery of percutaneous catheters into the peripheral vasculature.

### CONTRAINDICATIONS

None known

### WARNINGS

- The hydrophilic wire guide may slide entirely into the catheter, sheath introducer, vessel dilator or other device because of its low sliding friction. To prevent this, keep at least 5 cm of the wire protruding from the device fitting at all times.
- To prevent possible tissue damage, care should be taken when manipulating a device over a wire guide during the device's placement and withdrawal. If resistance is felt during device placement, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of resistance cannot be determined, remove the wire guide and device as a unit to prevent possible damage and/or complications.
- When using a wire guide, potential exists for thrombus formation or emboli, arterial or venous wall damage and/or plaque dislodgment. The physician should be familiar with the literature concerning the complications of angiography.
- Cook does not recommend a particular technique for the use of this wire guide. The steps contained in the following directions are for information purposes only. Each physician should evaluate their appropriateness according to individual patient condition and his or her medical training and experience.


### PRECAUTIONS

- Do not resterilize and/or reuse this device.
- Carefully read all instructions prior to use. Observe all warnings and cautions. Failure to do so may result in complications.
- **Avoid manipulating or withdrawing the hydrophilic wire guide back through a metal needle or cannula. A sharp edge may scrape the coating or shear the wire guide. A catheter, introducer sheath or vessel dilator should replace the needle as soon as the wire guide has been inserted in the vessel.**
- It is recommended that a plastic torque device be used to handle the hydrophilic wire guide. Use of a metal torque device may damage the wire.

### INSTRUCTIONS FOR USE

#### Preparation

1. Before removing the hydrophilic wire guide from its dispenser, inject sterile heparinized saline solution into the luer lock hub end of the dispenser.
2. Inject enough solution to fill the dispenser coil. This will completely cover the wire guide surface and activate the hydrophilic coating.
3. Remove the hydrophilic wire guide from its dispenser by gently withdrawing the wire's tip.
4. If the hydrophilic wire guide cannot easily be removed from its dispenser, inject more heparinized saline solution into the dispenser and then try again.

#### Wire Guide Use

1. Fill catheter or other device with heparinized saline solution before and during use to ensure smooth movement of the hydrophilic wire guide within the device.
2. When wet with saline solution or blood, the hydrophilic wire guide will become lubricious. Use of sterilized gauze moistened with heparinized saline solution facilitates handling the wire.
3. Insert the wire guide into the device and advance to desired position.

**NOTE:** If movement of the wire within the device becomes diminished, remove the wire guide and reactivate the hydrophilic coating by wetting its entire surface with heparinized saline solution.

#### After Use

1. Excess blood can be removed from the wire guide by wiping the surface once or twice with gauze moistened with heparinized saline solution. **DO NOT use dry gauze** as this may damage the wire surface resulting in increased resistance when the wire is reinserted into the device.
2. Use of alcohol, antiseptic solutions or other solvents must be avoided, because they may adversely affect the surface of the hydrophilic wire guide.

3. After cleaning the wire, replace it by its proximal end into its dispenser filled with heparinized saline solution. This wire guide can only be used during the same procedure on the same patient.
4. Upon reinsertion into the device, if the wire guide does not move with initial ease, exchange it for a new UniGlide™ Hydrophilic Wire Guide.

#### HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

#### REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd