EN
3**HSMA Holmium Laser Fiber (Single-Use)**

Instructions for Use

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Fig. 2

HSMA HOLMIUM LASER FIBER (SINGLE-USE)

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

INTENDED USE

This device is intended for incision/excision, ablation, and coagulation (hemostasis) when attached to a cleared Ho:YAG laser system comprised of any standard SMA-type connector and an output power of up to 100 watts, for the indications for which the laser system has been cleared.

CONTRAINDICATIONS

None known

WARNINGS

- All personnel present in the laser hazard zone must wear all suggested protective devices. Wear laser safety eyewear per the laser manufacturer's specifications.
- Do not attempt to reprocess.
- Baskets, wire guides and other ureteroscopic accessories may be damaged by direct contact with the laser treatment beam.
- Do not bend fiber at sharp angles. If visible light (aiming beam) can be seen leaking from the fiber, fiber failure may result when therapeutic energy is applied as the fiber is deflected beyond the optical limits of total internal reflection.
- Fiber should not be clamped with forceps or other securing instruments as it could result in fiber damage or breakage.

PRECAUTIONS

- A number of different factors affect the life of any particular fiber, including:
 - Extended lasing at high power
 - Continuous lasing with the fiber tip in contact with tissue
 - Lasing with a contaminated or damaged proximal end
 - Improper handling
 - Poor laser beam alignment or focus
- Never subject fiberoptics to sharp bends in handling, use, or storage.
- Always keep connector-end dry and free from contaminants.
- Discard any fiberoptic assembly that is cracked or broken, or does not meet minimum transmission standards.
- Ferrule dust cap should stay attached when the fiber is not connected to the laser system.
- Do not use this laser fiber in the presence of flammable anesthetics or combustible materials.
- Do not exceed recommended power limits.

Fiber Core Size	Power Limit
150 μm	6 W
200 μm	10 W
273 μm	15 W

Fiber Core Size	Power Limit
365 μm	20 W
550 μm	100 W
940 μm	100 W

NOTE: In addition to above recommended fiber handling protocol and precautions, all other applicable hospital protocols required to ensure patient and hospital personnel safety should be followed.

NOTE: Due to the extremely small fiber core size, the output power from 150 or 200 μm core fibers may be more than 20% lower than set power.

INSTRUCTIONS FOR USE

1. Open and remove the packaging tray from the sterile pouch.
2. Remove the Fiber Clip from the packaging tray.
3. Peel and remove the backing from the Fiber Clip. (**Fig. 1**)
4. Adhere the clip to the desired location.
5. Remove the proximal connector from the packaging tray.

NOTE: Keep output of tubing in line with fiber to minimize drag as the fiber is exiting the tubing.

6. Hand the proximal connector to a non-sterile attendant while maintaining sterility of distal fiber.

NOTE: The laser system is NOT sterile. Anyone who touches the laser system to insert the fiber will no longer be "scrubbed in."

7. Open any laser aperture cover or door and insert the proximal connector into the laser system; screw in finger-tight.
8. Slide the fiber into the lower clip to hold the fiber in place.

NOTE: Use the upper clip to hold the fiber tip when not in use. (**Fig. 2**)

9. Follow the laser manufacturer's directions for use.
10. Remove the fiber from the Fiber Clip and discard.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd