COOK

MEDICAL

CE
0088

EN
2

Haskal Transjugular Intrahepatic Portal Access Set

Instructions for Use

DA
3

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T _ H T I P A _ R E V 4

HASKAL TRANSJUGULAR INTRAHEPATIC PORTAL ACCESS SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Haskal Transjugular Intrahepatic Portal Access Set consists of:

- TFE transjugular catheter and transjugular needle
- Check-Flo® II Introducer set
- Sheath with radiopaque band
- Aurous Centimeter Sizing Catheter
- Torcon NB Advantage angiographic catheter
- TFE Van Andel angiographic catheter
- The Firm LT Roadrunner Wire Guide
- TFE wire guide, Amplatz Extra Stiff
- TFE curved Newton LLT Wire Guide
- TFE Newton LLT Wire Guide
- 12 Fr dilator
- Micropuncture Introducer Set:
 - coaxial catheter
 - entry needle
 - stainless steel wire guide

INTENDED USE

The Haskal Transjugular Intrahepatic Portal Access Set is intended for transjugular liver access in diagnostic and interventional procedures.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Manipulation of products requires fluoroscopic guidance.
- The patient should be monitored (EKG, oximetry, blood pressure, pulse) during the procedure.
- After introduction and periodically during the procedure, the transjugular introducer sheath should be flushed with saline.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

- Intraperitoneal hemorrhage
- Puncture site hematoma
- Cardiac arrhythmia
- Arteriovenous fistula

INSTRUCTIONS FOR USE

1. Utilizing micropuncture or standard 18 gage access techniques, introduce an appropriate .035 inch wire guide into the inferior vena cava via jugular venous access.
2. Using the selective hepatic catheter and a wire guide of choice, catheterize the right hepatic vein (RVH) or the most adequate hepatic vein branch. Leave the wire guide in a safe, distal position, and remove the catheter.
3. Advance the transjugular introducer sheath assembly over the wire guide as distal as possible into the hepatic vein. Remove the wire and dilator.
4. Introduce the TFE transjugular catheter and needle assembly through the sheath. The needle should not protrude past the catheter end hole.
5. Advance the catheter/needle assembly into the hepatic vein. Orient the TFE transjugular catheter inferiorly and rotate anteriorly (arrow baseplate on needle indicates direction of needle curve).
6. Wedging the catheter/needle assembly against the vein wall, thrust the needle in one movement forward through the hepatic parenchyma and toward the portal system.
7. Connect a syringe with contrast medium to the needle, apply suction and withdraw the needle until blood is seen. Inject a small amount of contrast medium in order to confirm position in the portal system.

8. Advance the TFE transjugular catheter over the transjugular needle further into the portal vein. Remove the transjugular needle.
9. Introduce a Newton LLT wire guide of choice into the portal branch and select the main portal vein (MPV).
10. Advance the TFE transjugular catheter and sheath into position across the parenchymal tract.
11. Advance the TFE Van Andel catheter through the TFE transjugular catheter into the portal vein.
12. Proceed with any indicated interventional procedures. (For critical measurement of the tract length or to ascertain magnification factor, the Auros Sizing Catheter may be placed across the tract and used for injection.)

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2014
Bloomington, IN 47404 U.S.A.