EN 3	HiWire™ Hydrophilic Wire Guide Instructions for Use
CS 4	Hydrofilní vodící drát HiWire™ Návod k použití
DA 5	HiWire™ hydrofil kateterleder Brugsanvisning
DE 6	HiWire™ hydrophiler Führungsdraht Gebrauchsanweisung
EL 7	Υδρόφιλος συρμάτινος οδηγός HiWire™ Οδηγίες χρήσης
ES 8	Guía hidrofílica HiWire™ Instrucciones de uso
FR 9	Guide hydrophile HiWire™ Mode d'emploi
HU 11	HiWire™ hidrofíl vezetődrót Használati utasítás
IT 12	Guida idrofila HiWire™ Istruzioni per l'uso
NL 13	HiWire™ hydrofiele voerdraad Gebruiksaanwijzing
NO 14	HiWire™ hydrofil ledevaier Bruksanvisning
PL 15	Hydrofilny przewodnik HiWire™ Instrukcja użycia
PT 16	Fio guia hidrófilo HiWire™ Instruções de utilização
SV 17	HiWire™ hydrofil ledare Bruksanvisning

HIWIRE™ HYDROPHILIC WIRE GUIDE

Federal law (USA) restricts this device to the sale and use by or on the order of a physician.

PACKAGE CONTENTS:

1 - HiWire™ Hydrophilic Wire Guide

NOTE: Wire guide Length, Diameter, Coating, Tip, and Core configurations are indicated on the product label.

DESCRIPTION:

Utilizing proprietary processes, the HiWire™ Hydrophilic Wire Guide is constructed from a steerable, metallic core with a polymer coating. A hydrophilic coating is applied over the radiopaque polymer jacket.

INTENDED USE:

To facilitate the placement of devices during diagnostic and interventional procedures.

STORAGE:

Store in a cool, dry place.

Standard Shaft Configuration

Stiff Shaft Configuration


Angled Tip Configuration

Straight Tip Configuration

PREPARATION FOR USE:

1. Before removing the hydrophilic wire guide from its dispenser, inject sterile heparinized saline solution into the luer lock hub end of the dispenser.
2. Inject enough solution to fill the dispenser coil. This will completely cover the wire guide surface and activate the hydrophilic coating.
3. Remove the hydrophilic wire guide from its dispenser by gently withdrawing the wire's tip.
4. If the hydrophilic wire guide cannot easily be removed from its dispenser, inject more heparinized saline solution into the dispenser and then try again.

DIRECTIONS FOR USE:

1. Fill catheter or other device with heparinized saline solution before and during use to ensure smooth movement of the hydrophilic wire guide within the device.
2. When wet with saline solution or blood, the hydrophilic wire guide will become lubricious. Use of sterilized gauze moistened with heparinized saline solution facilitates handling the wire.
3. Insert the wire guide into the device and advance to desired position.

NOTE: If movement of the wire within the device becomes diminished, remove the wire guide and reactivate the hydrophilic coating by wetting its entire surface with heparinized saline solution.

AFTER USE:

1. Excess blood can be removed from the wire guide by wiping the surface once or twice with gauze moistened with heparinized saline solution. **DO NOT use dry gauze** as this may damage the wire surface resulting in increased resistance when the wire is reinserted into the device.
2. Use of alcohol, antiseptic solutions or other solvents must be avoided, because they may adversely affect the surface of the hydrophilic wire guide.
3. After cleaning the wire, replace it by its proximal end into its dispenser filled with heparinized saline solution. This wire guide can only be used during the same procedure on the same patient.
4. Upon reinsertion into the device, if the wire guide does not move with initial ease, exchange it for a new HiWire™ Hydrophilic Wire Guide.

WARNINGS:

1. The hydrophilic wire guide may slide entirely into the catheter, sheath introducer, vessel dilator or other device because of its low sliding friction. To prevent this, keep at least 5 cm of the wire protruding from the device fitting at all times.
2. To prevent possible tissue damage, care should be taken when manipulating a device over a wire guide during the device's placement and withdrawal. If resistance is felt during device placement, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of resistance cannot be determined, remove the wire guide and device as a unit to prevent possible damage and/or complications.
3. When using a wire guide, potential exists for thrombus formation or emboli, arterial or venous wall damage and/or plaque dislodgment. The physician should be familiar with the literature concerning the complications of angiography.
4. Cook does not recommend a particular technique for the use of this wire guide. The steps contained in the preceding directions are for information purposes only. Each physician should evaluate their appropriateness according to individual patient condition and his or her medical training and experience.

CAUTIONS:

- Prior to use, inspect for damage. If damaged, DO NOT USE.
- This device is sterilized by ethylene oxide gas and intended for one-time use only. Do not use any unit if its package is opened or damaged. Do not resterilize and/or reuse.
- Carefully read all instructions prior to use. Observe all warnings and cautions. Failure to do so may result in complications.

- Avoid manipulating or withdrawing the hydrophilic wire guide back through a metal needle or cannula. A sharp edge may scrape the coating or shear the wire guide. A catheter, introducer sheath or vessel dilator should replace the needle as soon as the wire guide has been inserted in the vessel.
- It is recommended that a plastic torque device be used to handle the hydrophilic wire guide. Use of a metal torque device may damage the wire.

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© COOK 2015
Bloomington, IN 47404 U.S.A.