COOK

MEDICAL

CE
0088

EN
3

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Bruksanvisning

T _ I D K B L _ R E V O

MREYE® KOPANS BREAST LESION LOCALIZATION NEEDLES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The MReye Kopans Breast Lesion Localization Needle consists of a bevel-tip needle with centimeter markings and a hookwire. The hookwire features a burnish mark that allows the operator to ensure that the hookwire is within the needle tip during needle manipulation. The needles are manufactured from a non-ferromagnetic material to prevent occurrence of significant image distortion during Magnetic Resonance Imaging (MRI).

INTENDED USE

The MReye Kopans Breast Lesion Localization Needle is intended for preoperative marking of nonpalpable breast lesions.

CONTRAINDICATIONS

None known

WARNINGS

- Under no circumstances should a hookwire engaged in tissue be pulled out without surgical intervention.
- **Do not exceed a magnet strength of 1.5 Tesla.**

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for localization of breast lesions should be employed.
- Following placement of the hookwire, the portion protruding outside of the breast should be bent and taped to the skin to prevent inadvertent movement.
- Final hookwire position should be confirmed by appropriate imaging modality.
- The hookwire should be used as a guide for the surgeon, not a retractor.
- Circumferential image distortion may occur around the body of the needle when the needle is oriented perpendicular to the magnetic field.
- RF heating of the needle may occur during MRI procedures. Heating effects are theoretically greatest when the needle is oriented along the length of the bore, near the bore wall.
- If using fast scan protocols, there is a theoretical possibility that peripheral nerve stimulation may occur under worst-case conditions (needle with inserted trocar oriented along the length of the bore, near the bore wall).
- Before using MReye Needles, physicians should have training and experience with interventional Magnetic Resonance Imaging.
- Different sizes of needle and hookwire on the image may be obtained. The following scanner parameters are recommended:
A spin-echo sequence (TR/TE 200ms/15ms) will provide an image depicting the needle and hookwire 1:1.
A gradient-echo sequence (TR/TE/flip angle 336ms/5ms/90°) will provide an image depicting the needle and hookwire enlarged approximately 2:1.
- Evaluation of image distortion was performed using a 1.5 Tesla scanner and SE (Spin Echo) and GRASS (Gradient-Recalled Acquisition in the Steady State) MRI sequences, with the device in two different orientations (0° and 90°) and the MRI parameters as shown in the following table. Results of this testing show that under worst-case conditions in terms of artifact size and needle position, circumferential image distortion of up to 2 cm may occur around the body of the needle when the needle is oriented perpendicular to the magnetic field using the GRASS sequence. Under these conditions, error in determining needle tip position may be as much as 1 cm.

MRI Parameters for Image Distortion Testing

MRI Pulse Sequence

Parameter	SE (Spin Echo)	GRASS (Gradient-Recalled Acquisition in Steady State)
Static Field Strength	1.5 Tesla	1.5 Tesla
Needle Orientation	0°, 90°	0°, 90°
TR	300 msec	100 msec
TE	16 msec	20 msec
Flip Angle	90°	30°
Bandwidth	16 kHz	16 kHz
Field of View	Axial, Sagittal: 24x12 Coronal: 40x30	Axial, Sagittal: 24x12 Coronal: 40x30
Frequency	Axial (0°): "x"	Axial (0°): "x"
Encoding Gradient	Sagittal: (90°): "z"	Sagittal: (90°): "z"
Direction	Coronal: (0°, 90°): "z"	Coronal: (0°, 90°): "z"
Slice Thickness	5 mm	5 mm
Data Acquisition Matrix Size	256x256	256x256
Maximum Gradient Field Strength	1 G/cm	1 G/cm
Phantom Filler T1	2 ¹ sec	2 ¹ sec
Phantom Filler T2	2 ¹ sec	2 ¹ sec

¹ Aqueous solution of 0.45% saline (phantom filler) is assumed to have T1 and T2 value of approximately 2 seconds (similar to water).

POTENTIAL ADVERSE EVENTS

- Infection
- Internal bleeding and/or bleeding at the puncture site
- Pneumothorax
- Puncture or injury to target organ or to a nearby organ that is traversed by the needle
- Hemoptysis
- Peritonitis
- Damage of tissue adjacent to the needle

INSTRUCTIONS FOR USE

NOTE: The burnish mark on the hookwire should be seen just outside the needle hub to ensure that the hookwire is within the needle tip during needle manipulation.

1. Utilizing imaging guidance, introduce the needle into the lesion.
2. During MR procedures, contrast medium may be injected for indication of the tumor.
3. Check the position of the needle.
4. Advance and release the hookwire.
5. Remove the needle.
6. Bend the hookwire protruding from the breast and tape flat to the skin.
7. Verify final hookwire position.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.