EN 3	Breast Lesion Needle with Müller-Schimpfle Localization Coil Instructions For Use
CS 4	Jehla s lokalizační spirálkou Müller-Schimpfle k lokalizaci léze v prsu Návod k použití
DA 7	Brystlæsionsnål med Müller-Schimpfle lokaliseringsspiral Brugsanvisning
DE 9	Brustläsionskanüle mit Müller-Schimpfle Lokalisationsspirale Gebrauchsanweisung
EL 11	Βελόνα βλάβης μαστού με σπείρα εντοπισμού Müller-Schimpfle Οδηγίες χρήσης
ES 13	Aguja para lesiones mamarias con espiral de localización Müller-Schimpfle Instrucciones de uso
FR 15	Aiguille pour lésion mammaire à spirale de localisation Müller-Schimpfle Mode d'emploi
HU 17	Emlő léziós tű Müller-Schimpfle lokalizáló spirállal Használati utasítás
IT 19	Ago per lesioni mammarie con spirale di localizzazione Müller-Schimpfle Istruzioni per l'uso
NL 21	Borstlaesienaald met Müller-Schimpfle-lokalisatiecoil Gebruiksaanwijzing
NO 23	Brystlesjonsnål med Müller-Schimpfle lokaliseringsspiral Bruksanvisning
PL 25	Igła do zmian w piersi ze zwojem lokalizującym Müller-Schimpfle Instrukcja użycia
PT 28	Agulha com espiral de localização Müller-Schimpfle para lesão mamária Instruções de utilização
SV 30	Bröstlesionsnål med Müller-Schimpfle lokaliseringsspiral Bruksanvisning
ZH 32	帶Müller-Schimpfle定位環的乳腺病灶定位針 使用说明

BREAST LESION NEEDLE WITH MÜLLER-SCHIMPFLE LOCALIZATION COIL

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Breast Lesion Needle with Localization Coil consists of a two-part positioning needle with a trocar tip and a pusher, and a coil preloaded on a stylet with a cannula stopper. The needle tip is etched for improved visibility during ultrasound imaging and the needle cannula is marked in 1 cm increments, with double markings at 5 cm increments for monitoring of needle depth during the procedure.

INTENDED USE

The Breast Lesion Needle with Localization Coil is intended for preoperative marking of non-palpable breast lesions, where the localization coil serves as a marker for the surgeon.

MR STATEMENTS

Non-clinical testing demonstrated that the Breast Lesion Needle with Localization Coil is MR conditional. A patient with this device can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less
- Maximum whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg for 15 minutes of scanning.

In non-clinical testing, the Breast Lesion Needle with Localization Coil produced a temperature rise of less than 2.6°C at a maximum whole-body-averaged specific absorption rate (SAR) of 2.9 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 1.5 Tesla/64 MHz, Magnetom, Siemens Medical Solutions MR scanner.

In non-clinical testing, the Breast Lesion Needle with Localization Coil produced a temperature rise of less than 3.6°C at a maximum whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 3 Tesla/128 MHz, Excite, Software G3.0-052B, General Electric Healthcare MR scanner.

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the following items for the Breast Lesion Needle with Localization Coil:

- (1) Stylet, needle, and coil in place
- (2) Coil only

Therefore, optimization of MR imaging parameters to compensate for the presence of these devices may be necessary.

(1) Stylet, needle cannula, and coil

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	2,003-mm ²	21-mm ²	1,388-mm ²	177-mm ²
Imaging Plane	parallel	perpendicular	parallel	perpendicular

(2) Coil only

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	181-mm ²	96-mm ²	332-mm ²	356-mm ²
Imaging Plane	parallel	perpendicular	parallel	perpendicular

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in breast lesion localization and MR techniques. Standard techniques for needle placement should be employed.
- Final coil position should be confirmed by appropriate imaging modality.

INSTRUCTIONS FOR USE

1. Introduce the two-part needle into the lesion.
2. Check the position of the needle.
3. Remove the inner trocar, leaving the needle cannula in place.
4. Introduce the stylet with the preloaded coil and stopper cannula.
5. Remove the stylet, leaving the coil and stopper cannula in place.
6. Remove the stopper cannula, leaving the coil in the needle cannula.
7. Check the position of the needle cannula.
8. Introduce the pusher into the needle cannula and advance it to release the coil in the lesion.
9. Verify final coil position.
10. Remove the needle cannula.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.