- EN**
3 **Double Pigtail Stent Set**
Instructions for Use
- CS**
5 **Souprava double pigtailového stentu**
Návod k použití
- DA**
8 **Stentsæt med dobbelt grisehale**
Brugsanvisning
- DE**
11 **Doppelpigtail-Stentset**
Gebrauchsanweisung
- EL**
14 **Σετ ενδοπρόσθεσης διπλού σπειροειδούς άκρου**
Οδηγίες χρήσης
- ES**
18 **Equipo de stent de doble pigtail**
Instrucciones de uso
- FR**
20 **Set d'endoprothèse double pigtail**
Mode d'emploi
- HU**
23 **Kettős pigtail sztentkészlet**
Használati utasítás
- IT**
26 **Set con stent a doppio pigtail**
Istruzioni per l'uso
- NL**
29 **Dubbele pigtail stentset**
Gebruiksaanwijzing
- NO**
32 **Stentsett med dobbel pigtail**
Bruksanvisning
- PL**
35 **Zestaw stentu z dwiema końcówkami pigtail**
Instrukcja użycia
- PT**
38 **Conjunto de stent de duplo pigtail**
Instruções de utilização
- SV**
41 **Stentset med dubbla pigtails**
Bruksanvisning

DOUBLE PIGTAIL STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

- Radiopaque double pigtail stent (Polyurethane, C-Flex®, C-Flex® TPE, Endo-Sof™, Sof-Flex®, or LSe Sof-Flex®)
- Radiopaque ureteral catheter
- Radiolucent stent positioner
- Wire guide (*size and configuration may vary*)

NOTE: Set components may vary.

NOTE: Some stents are available with an AQ® hydrophilic coating.

NOTE: If no wire guide is provided with this set, the following is recommended:

3.0 Fr. accepts .018"	7.0 Fr. accepts .038"
3.7 Fr. accepts .025"	8.0 Fr. accepts .038"
4.0 Fr. accepts .025"	8.2 Fr. accepts .038"
4.7 Fr. accepts .035"	10 Fr. accepts .038"
5.0 Fr. accepts .035"	12 Fr. accepts .038"
6.0 Fr. accepts .038"	

INTENDED USE

The Double Pigtail Stent is used for temporary internal drainage from the ureteropelvic junction to the bladder. Double pigtail ureteral stents have been employed to relieve obstruction in a variety of benign, malignant, and post-traumatic conditions. These stents may be placed using endoscopic, percutaneous, or open surgical techniques.

The wire guide is used to gain ureteral access, to establish a tract, and to assist in the placement, replacement, and exchange of medical devices during urological and gynecological procedures.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- Manipulation of the wire guide requires appropriate imaging control. Use caution not to force or over manipulate the wire guide when gaining access.
- When using a wire guide through a metal cannula/needle, use caution as damage may occur to the outer coating.
- When exchanging or withdrawing an instrument over the wire guide, secure and maintain the wire guide in place under fluoroscopy in order to avoid unexpected wire guide displacement.
- These wire guides are not intended for PTCA use.

- Complications of ureteral stent placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to each patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures
- Endo-Sof™ stents must not remain indwelling more than twelve months. Polyurethane, C-Flex®, Sof-Flex®, and LSe Sof-Flex®, stents must not remain indwelling more than six months. If the patient's status permits, the stent may be replaced with a new stent.
- These stents are not intended as permanent indwelling devices.
- Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered.
- A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the wire guide or stent should be avoided. Use of a 0-degree scope lens is recommended. Scopes larger than 21.0 French are suggested.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.
- Hydrophilically coated wires are very slippery when wet. Always maintain control of the wire guide when manipulating it through any device.
- If problems occur using this device, please call your Cook Urological sales representative or contact our Customer Service department at the address/phone listed at www.cookmedical.com.

INSTRUCTIONS FOR USE

Wire Guide Preparation

NOTE: If your device did not include a wire guide, then the instructions below are not applicable. If your device did include a wire guide, please refer to the instructions below for proper use.

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. If using a hydrophilically coated wire guide, fill a syringe with sterile water or sterile saline solution and attach it to the flushing port on the wire guide holder. Inject enough solution to wet the wire guide surface entirely. This will activate the hydrophilic coating. **NOTE:** For optimal performance, rehydrate the hydrophilically coated wire guide after exposure to ambient environment or, after extended use, replace it with a new hydrophilically coated wire guide.
3. Remove the wire guide from the holder.
4. If using a movable core wire guide, adjust mandril position by advancing or withdrawing the mandril within the wire guide to achieve the required degree of flexibility.

Endoscopic Stent Placement

NOTE: If using an AQ stent, immerse the stent in sterile water or isotonic saline prior to use to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

1. Pass the flexible wire guide tip beyond the obstruction to the renal pelvis. Tortuosity in the obstructed ureter often can be resolved using the wire guide and an open-ended ureteral catheter in combination.

NOTE: Endhole size and length of the device must be taken into consideration to ensure proper fit between the wire guide and the device.

2. Using a baseline pyelogram, estimate the length of stent needed by taking the ureteral length and adding 1 cm. Accurate measurement enhances drainage efficiency and patient comfort.

3. Pass the stent over the wire guide through the cystoscope. Under direct vision, advance the stent into the ureter with the stent positioner. Have an assistant hold the wire guide in position to prevent advancement of the wire guide into the renal parenchyma.
4. Watch for the distal end of the stent at the ureterovesical junction. At that point, halt advancement of the stent. As an assistant removes the wire guide, hold the stent in position with the positioner. The stent pigtail will form spontaneously. Carefully remove the positioner from the cystoscope.

NOTE: If necessary, final adjustment can be made with endoscopic forceps. The stent may be removed easily by gentle withdrawal traction using endoscopic forceps. Fluoroscopy facilitates stent placement; however, standard radiography may be used.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "Self Retained Internal Ureteral Stents: A New Approach," *The Journal of Urology*, 119 (1978), 731.

T. W. Hepperlen, H. K. Mardis, H. Kammandel: "The Pigtail Ureteral Stent in the Cancer Patient," *The Journal of Urology*, 121 (1979), 17.

M. F. Camacho, R. Pereiras, H. Carrion, M. Bondhus, V. A. Politano: "Double-Ended Pigtail Ureteral Stent: Useful Modification To Single End Ureteral Stent," *Urology*, 13 (1979), 516-520.

M. A. Koyle, A. D. Smith: "Ureteral Stents," Chapter 31 in A. D. Smith, W. R. Castañeda-Zuñiga, J. G. Bronson: *Endourology: Principles and Practice*. Thieme Inc., New York, New York, 1986.

C-FLEX® is a registered trademark of Concept Polymer Technologies, Inc.

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© COOK 2012
C-FLEX® är ett registrerat varumärke som tillhör Concept Polymer Technologies, Inc.
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.