COOK**MEDICAL****CE**
0123

EN 2	Integral Filiform Urethral Dilator Instructions for Use
BG 3	Интегрален филиформен уретрален дилататор Инструкции за употреба
CS 5	Integrovaný filiformní uretrální dilatátor Návod k použití
DA 6	Integral filiform dilatator til urethra Brugsanvisning
DE 8	Integrierter filiformer Urethradilatator Gebrauchsanweisung
EL 9	Ενσωματωμένος νηματοειδής διαστολέας ουρήθρας Οδηγίες χρήσης
ES 11	Dilatador uretral filiforme integral Instrucciones de uso
ET 12	Integreeritud kapillaarotsaga ureetra dilataator Kasutusjuhend
FI 14	Kiinteä filiforminen virtsaputken laajennin Käyttöohjeet
FR 15	Dilatateur filiforme urétral intégré Mode d'emploi
HR 17	Integrirani filiformni dilatator uretre Upute za uporabu
HU 18	Integral filiform urethrális dilatátor Használati utasítás
IT 20	Dilatatore uretrale con guida filiforme integrata Istruzioni per l'uso
LT 21	Vientisas siūlo formos šlaplės plėtiklis Naudojimo instrukcija
LV 23	Viendalīgs diegveida urīnizvadkanāla dilatators Lietošanas instrukcija
NL 24	Integrale filiforme urethrale dilatator Gebruiksaanwijzing
NO 26	Integrert filiform urinrørdilatator Bruksanvisning
PL 27	Zintegrowany nitkowany rozszerzacz do cewki moczowej Instrukcja użycia
PT 29	Dilatador uretral filiforme integral Instruções de utilização
RO 30	Dilatator uretral filiform integral Instrucţiuni de utilizare
SK 32	Integrálny nitkovitý uretrálny dilatátor Návod na použitie
SL 33	Integralni tanki dilatator za sečnico Navodila za uporabo
SV 35	Integrerad uretral dilatator för filiform Bruksanvisning
TR 36	Entegre Filiform Üretral Dilatör Kullanma Talimatı

T _ I F U D - M _ R E V O

INTEGRAL FILIFORM URETHRAL DILATOR

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

The Integral Filiform Urethral Dilator is manufactured from non-radiopaque tubing in a total length of 38 cm. Integral Filiform Urethral Dilator is manufactured with an 18 Fr shaft tapered and bonded to a 7 Fr non-radiopaque filiform tip. The filiform tip measures 6 cm in length. The device French size is indicated on the end of the dilator without the filiform tip.

Performance Characteristics

- Filiform tip aids in traversing and dilating a urethral stricture, eliminating the need for an additional filiform
- Bonded taper will allow a gradual dilatation of the stricture

Device Compatibility

This device does not have a specific intended device compatibility.

Patient Population

The target patient population for these devices is adult patients with urethral strictures who require dilation of the stricture areas of the urethra.

Intended User

These devices are intended for physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact With Body Tissue

These devices interact with the urethral tract including the urethra and bladder.

Operating Principle

These devices are intended to be introduced directly into the urethra. The filiform tip can traverse a urethral stricture, eliminating the need for an additional filiform. The filiform tip gradually transitions to a larger diameter dilator portion to dilate the stricture.

INTENDED USE

These devices are dilators intended for traversing and dilating urethral strictures by using a filiform tip that gradually transitions to a larger diameter.

INDICATIONS FOR USE

These devices are used for traversing and dilating urethral strictures in adult patients.

CLINICAL BENEFITS

Urethral dilation to widen the diameter of the urethra.

CONTRAINDICATIONS

None

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use

PRECAUTIONS

No Known Precautions

POTENTIAL ADVERSE EVENTS

- Bleeding
- False passage
- Fistula
- Hematoma
- Hematuria
- Infection
- Pain
- Perforation, urethral
- Rectal injury
- Urethral injury
- Lower urinary tract symptoms (poor/diminished urinary stream, dysuria, incomplete emptying of the bladder, (recurrent) urinary tract infections,

overflow incontinence, frequency, hesitancy, urgency, urinary retention, increasing post void residual, and/or pain)

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

Introduce the filiform tip of the dilator into the urethral meatus and slowly advance it beyond the obstruction.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions contraindications, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical. This shall also be reported to the competent authority the device was used in.

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© COOK 2024
Cook Medical Europe Ltd.
William A. Cook Australia Pty. Ltd.