COOK

MEDICAL

CE
0123

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ES Vista lateral del asa

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INDY OTW® VASCULAR RETRIEVER

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Indy OTW® Vascular Retriever consists of a 4-loop wire snare attached to an inner catheter (which is attached proximally to a metal cannula) and an outer Flexor sheath. The device is compatible with a 0.035 inch wire guide for over-the-wire introduction. The wire snare is composed of nitinol, for shape memory, and platinum, for visibility. The Flexor sheath has a radiopaque band incorporated within the sheath material to identify the location of the sheath's distal tip.

Performance Characteristics

- Catheter inner lumen for over-the-wire use.
- Snare with overlapping loops expands when extended from and collapses when retracted into sheath to facilitate foreign body capture.
- Radiopaque wire incorporated into snare material and radiopaque marker band incorporated into sheath material for positioning under imaging techniques.
- Tuohy-Borst valve aids in manipulation or retrieval of foreign body by maintaining tension on the snare catheter.
- Hemostasis valve minimizes leakage.
- Stainless steel coiling incorporated into the sheath material for resistance to kinking.
- Rounded catheter and sheath tip to facilitate atraumatic insertion into patient anatomy.

Device Compatibility

- Standard luer syringe.
- 0.035 inch wire guide.
- Guiding catheter or introducer sheath with a larger inner diameter than the outer diameter of the Indy sheath.
- Heparinized saline.

Patient Population

The target population for this device is adult patients with an intravascular foreign body who require and are suitable for endovascular intervention based on the assessment of their treating physician/healthcare provider. The underlying patient condition, the coagulation status of the patient, and the location of the foreign body should be considered.

Intended User

This device is intended for use by physicians trained and experienced in intravascular diagnostic and interventional techniques.

Contact with Body Tissue

This device contacts circulating blood for limited duration (≤ 24 hours).

Operating Principle

This device manipulates or retrieves a foreign body from the vasculature by capturing the body in wire loops and either repositions the object or retrieves the object through an *in situ* catheter or sheath. The shape memory properties of the wire loops and ability of the wire loops to be collapsed by the extension of the Flexor sheath facilitate the capture of the foreign body. Fluoroscopic imaging identifies and facilitates control of both the wire loops and distal tip of the Flexor sheath throughout the procedure.

INTENDED USE

This device is intended to manipulate or retrieve a foreign body within the central and peripheral circulatory systems by capturing the foreign body in wire loops.

INDICATIONS FOR USE

This device is indicated for use in adult patients with medical conditions related to the presence of a vascular foreign body and for whom manipulation or retrieval of the foreign body would provide a clinical benefit.

CLINICAL BENEFITS

This device provides a minimally-invasive means to manipulate or retrieve an intravascular foreign body.

CONTRAINDICATIONS

This device is not intended to be placed in the coronary or neuro vasculature.

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- This device contains nitinol, a nickel titanium alloy. Allergic reactions to nickel should be considered.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

PRECAUTIONS

- Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- This device should be used in patients who are suitable for endovascular interventions with particular attention paid to the patient's coagulation status, presence of excessive thrombus formation on the foreign body, and the risk of vessel injury posed by retrieval.
- This device should be used in patients with vessel diameters compatible with the associated device components.
- Excessive force should not be used to manipulate or retrieve foreign objects.
- Manipulation of this device requires fluoroscopic control.
- Visually inspect this device before use to ensure it is undamaged. For example, the shaft should be free of kinks.
- Always check fit through the intended guiding catheter or introducer sheath prior to use.
- Do not use after the expiration date printed on the label.

POTENTIAL ADVERSE EVENTS

- Arrhythmia
- Bleeding / Hemorrhage
- Device and foreign body entrapment
- Embolization
- Heart injury
- Infection
- Myocardial infarction
- Pulmonary embolism
- Stroke
- Thrombosis
- Vessel injury

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Upon removal from package, ensure the outer diameter (OD) of the retriever is appropriate for the inner diameter of the introducer sheath or guiding catheter.
2. Using the side flush port, flush the Flexor sheath with heparinized saline.
3. Using the cannula flushing port, flush the inner catheter with heparinized saline.
4. Insert the retriever over-the-wire through an *in situ* guiding catheter or introducer sheath and advance it to the desired position.
5. While holding the Flexor sheath in position, loosen Tuohy-Borst and advance the inner catheter through the Flexor sheath until the snare emerges from the distal tip. **NOTE:** The snare will expand upon emergence from the Flexor sheath.
6. Position the retriever so that the foreign body, such as a wire guide or catheter, is positioned within the snare. While maintaining the snare position, slide the Flexor sheath forward to capture the foreign body.
7. Tighten the Tuohy-Borst to maintain tension on the catheter controlling the foreign body, withdraw the retriever, including the snare and Flexor sheath assembly, to a peripheral location.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of. You can use the basic unique device identifier for this device – (0827002CINC050-211_A-1J5) – to find the Summary of Safety and Clinical Performance (SSCP) on the EUDAMED website. When the EUDAMED is available, use the following link: ec.europa.eu/tools/eudamed

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2022
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd