
COOK

MEDICAL

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Fig.1

INDY OTW™ VASCULAR RETRIEVER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Indy OTW Vascular Retriever consists of a 4-loop nitinol wire snare attached to an inner catheter (which is attached proximally to a metal cannula) and an outer Flexor sheath. The device is compatible with a .035 inch wire guide for over-the-wire introduction. The Flexor sheath has a radiopaque band incorporated within the sheath material to identify the location of the sheath's distal tip. (**Fig. 1**)

INTENDED USE

The Indy OTW Vascular Retriever is intended to snare a foreign body and withdraw it to a peripheral vascular location.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Manipulation of the product requires fluoroscopic control.
- Visually inspect the product before use to ensure it is undamaged.
- The outer diameter of the device's Flexor sheath is 8 French (2.63 mm) nominal. Always check fit through the intended guiding catheter or introducer sheath prior to use.
- Do not use after the expiration date printed on the label.

INSTRUCTIONS FOR USE

1. Upon removal from package, ensure the outer diameter (OD) of the retriever is appropriate for the inner diameter of the introducer sheath.
2. Using the side flush port, flush the Flexor sheath with heparinized saline.
3. Using the cannula flushing port, flush the inner catheter with heparinized saline.
4. Insert the retriever over-the-wire through an *in situ* guiding catheter or introducer sheath and advance it to the desired position.
5. While holding the Flexor sheath in position, loosen Tuohy-Borst and advance the inner catheter through the Flexor sheath until the snare emerges from the distal tip. **NOTE:** The snare will expand upon emergence from the Flexor sheath.
6. Position the retriever so that the foreign body, such as a wire guide or catheter, is caught within the snare. While maintaining the snare position, slide the Flexor sheath forward to capture the foreign body.
7. Tighten Tuohy-Borst to maintain tension on the catheter controlling the foreign body, withdraw the retriever, including the snare and Flexor sheath assembly, to a peripheral location and retrieve the foreign body.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.