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

1

2

COOK CERVICAL RIPENING BALLOON WITH ADJUSTABLE STYLET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Cook Cervical Ripening Balloon is a silicone double balloon catheter with an adjustable-length malleable stylet. The maximum balloon inflation is 80 mL/balloon.

INTENDED USE

The Cook Cervical Ripening Balloon is indicated for mechanical dilation of the cervical canal prior to labor induction at term when the cervix is unfavorable for induction.

CONTRAINDICATIONS

- Patient receiving or planning to undergo exogenous prostaglandin administration
- Placenta previa, vasa previa, or placenta percreta
- Transverse fetal orientation
- Prolapsed umbilical cord
- Prior hysterotomy, classic uterine incision, myomectomy or any other full-thickness uterine incision
- Pelvic structural abnormality
- Active genital herpes infection
- Invasive cervical cancer
- Abnormal fetal heart-rate patterns
- Breech presentation
- Maternal heart disease
- Multiple gestational pregnancy
- Polyhydramnios
- Presenting part above the pelvic inlet
- Severe maternal hypertension
- Any contraindication to labor induction
- Ruptured membranes

WARNINGS

- Concomitant use of the Cook Cervical Ripening Balloon with exogenous prostaglandins may increase the risk of adverse events associated with prostaglandin administration, including, but not limited to: uterine hyperstimulation, impaired utero-placental circulation, tachysystole, uterine rupture, placental abruption, amniotic fluid embolism, pelvic pain, retained placenta, severe genital bleeding, shock, fetal bradycardia, fetal death, and maternal death.
- The stylet should only be used to traverse the tip of the catheter through the cervix and should be removed as soon as the uterine balloon is above the level of the internal uterine opening (internal os) prior to full insertion of the catheter. Aggressive insertion may result in injury to the baby.
- The product should not be left indwelling for longer than 12 hours.
- The safety and effectiveness of the Cook Cervical Ripening Balloon has not been established among women with an obstetrical history of low transverse caesarean section.
- The safety and effectiveness of extra-amniotic saline infusion with the Cook Cervical Ripening Balloon has not been established.
- Always inflate the balloon with a sterile saline. Never inflate with air, carbon dioxide or any other gas.

- Do not overinflate. Using excessive pressure to inflate the balloon on this device can cause the balloon to rupture.
- If spontaneous rupture of membranes occurs while the Cook Cervical Ripening Balloon is in place, there is a risk that the uterine balloon could become entangled in the umbilical cord, necessitating emergent cesarean delivery.

PRECAUTIONS

If fetal membranes rupture spontaneously while this device is in place, it is recommended that both balloons be deflated and the device removed in preparation for spontaneous active labor contractions.

POTENTIAL ADVERSE EVENTS

Risks associated with use of the Cook Cervical Ripening Balloon and labor induction may include, but are not limited to:

- Placental abruption
- Uterine rupture
- Spontaneous rupture of membranes
- Spontaneous onset of labor
- Device expulsion
- Device entrapment and/or fragmentation
- Maternal discomfort during and after insertion
- Failed dilation or need for cesarean delivery
- Cervical laceration
- Bleeding
- Risk of pre-term labor and birth in subsequent pregnancy

INSTRUCTIONS FOR USE

Patient Preparation

1. Perform an abdominal ultrasound to confirm singleton, vertex presentation and to rule out partial or complete placenta previa and/or placenta percreta.
2. Place the patient in the lithotomy position.
3. Insert a large vaginal speculum to gain cervical access.
4. Clean the cervix with an appropriate cleaning solution to prepare for device insertion.

Using The Moldable Stylet

1. Loosen the fitting on the proximal hub of the stylet and adjust the wire so that the distal tip of the stylet is even with the distal tip of the Cervical Ripening Balloon. (**Fig. 1**)
2. Tighten the fitting so that the wire does not move during manipulation and seat the adjustable handle firmly into the blue port labeled "S." (**Fig. 2**)
3. Use the Cervical Ripening Balloon with stylet to traverse the cervix if necessary.

NOTE: Once the cervix has been traversed and the uterine balloon is above the level of the internal uterine opening (internal os), remove the stylet before further advancing the catheter.

4. Advance the Cervical Ripening Balloon until both balloons have entered the cervical canal.
5. Inflate the uterine balloon with 40 mL of sterile saline using a standard 20 mL Luer-lock syringe through the red Check-Flo valve marked "U".
6. Once the uterine balloon is inflated, pull the device back until the uterine balloon is against the internal cervical os.
7. The vaginal balloon should now be visible outside the external cervical os. Inflate the vaginal balloon with 20 mL of sterile saline using a standard 20 mL Luer-lock syringe through the green Check-Flo valve marked "V".

8. Once the balloons are situated on each side of the cervix and the device has been fixed in place, remove the speculum.
 9. Add more fluid to each balloon in turn, in 20 mL increments, until each balloon contains 80 mL (maximum) of fluid. **NOTE:** Do not over inflate the balloon.
 10. If desired, the proximal end of the catheter may be taped to the patient's thigh.
- NOTE:** The device is not intended to be in place for longer than 12 hours. Time the placement of the device 12 hours prior to the planned induction.

Device Removal

Deflate both balloons through the corresponding valves marked "U" and "V" and remove vaginally.

NOTE: If the membranes rupture spontaneously before removal of the device, it is recommended to deflate the balloons and remove the device to facilitate active labor management.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

Atad J, Hallak M, Ben-David Y, Auslender R, Abramovici H. *Ripening and dilation of the unfavorable cervix for induction of labor by a double balloon device: experience with 250 cases.* Br J Obstet Gynaecol. 1997;104 (pt 1): 29.

Atad J, Hallak M, Auslender R, Porat-Packer T, Zarfati D, Abramovici H. *A randomized comparison of prostaglandin E2, oxytocin and the double-balloon device in inducing labor.* Obstet Gynecol. 1996; 87:223-227.

Sherman D, Frenkel E, Tovbin J, Arieli S, Caspi E, Bukovsky I. *Ripening of the unfavorable cervix with extra-amniotic catheter balloon: clinical experience and review.* Obstet Gynecol Survey. 1996; 51 (10): 621-627.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.