- EN**
3 **Cook® Silicone Balloon HSG Catheter**
Instruction for Use
- CS**
4 **Katetr pro hysterosalpingografii Cook® se silikonovým balónkem**
Návod k použití
- DA**
6 **Cook® SHG-cervikalkateter med ballon**
Brugsanvisning
- DE**
8 **Cook® Silikon-Hysterosalpingographie-Ballonkatheter**
Gebrauchsanweisung
- EL**
10 **Καθετήρας HSG με μπαλόνι από σιλίκι Cook®**
Οδηγίες χρήσης
- ES**
13 **Catéter de HSG con balón de silicona Cook®**
Instrucciones de uso
- FR**
15 **Sonde à ballonnet silicone Cook® pour HSG**
Mode d'emploi
- HU**
17 **Cook® szilikon ballonos HSG katéter**
Használati utasítás
- IT**
19 **Catetere per isterosalpingografia con palloncino in silicone Cook®**
Istruzioni per l'uso
- NL**
21 **Cook® HSG-katheter met silicone ballon**
Gebruiksaanwijzing
- NO**
23 **Cook® HSG-kateter med silikonballong**
Bruksanvisning
- PL**
25 **Cewnik Cook® z silikonowym balonem do HSG**
Instrukcja użycia
- PT**
27 **Cateter de balão em silicone Cook® para histerossalpingografia**
Instruções de utilização
- SV**
29 **Cook® HSG-kateter med silikonballong**
Bruksanvisning

COOK® SILICONE BALLOON HSG CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Cook Silicone Balloon HSG Catheter is a polyurethane double lumen catheter with a silicone balloon and syringe.

INTENDED USE / INDICATIONS FOR USE

The Cook Silicone Balloon HSG Catheter is used for delivery of contrast medium or saline into the uterine cavity and fallopian tubes for the evaluation of tubal patency using hysterosalpingography, or to access the uterine cavity for saline infusion sonohysterography.

CONTRAINDICATIONS

This device should not be used in the presence of:

- Hemorrhage
- Active pelvic infection
- Sexually transmitted disease
- Pregnancy

WARNINGS

Always inflate the balloon with a sterile liquid. Never inflate with air, carbon dioxide or any other gas. To reduce the risk of perforation, if any resistance is felt while inserting the catheter, do not force the catheter against the resistance.

PRECAUTIONS

- This procedure should only be performed by a physician or properly licensed practitioner trained in diagnostic and interventional techniques. Standard techniques for SHG and HSG procedures should be employed.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects. The check valve used for the balloon inflation is manufactured from PVC (polyvinylchloride), which contains DINP (Diisononyl phthalate). This portion of the device does not have direct patient contact. The catheter is not made with phthalates.
- Do not overinflate. Using excessive volume to inflate the balloon on this device can cause the balloon to rupture. Refer to product inflation line on the balloon device for appropriate balloon volume.

POTENTIAL ADVERSE EVENTS

- Hypersensitivity to contrast media
- Discomfort, pain/cramping
- Infection

INSTRUCTIONS FOR USE

Hysterosalpingography (HSG)

1. Pretest the catheter balloon by first priming the enclosed syringe and then filling it with a labeled volume of sterile saline. **NOTE:** Ensure tip of syringe is contacting check valve (black plug).
2. Fill the instillation syringe with sterile saline and eliminate air from the syringe.
3. Prime the instillation lumen of the catheter with contrast medium.

4. Position a vaginal speculum within the vagina so the external cervical os can be visualized during catheter insertion.
5. Cleanse the cervix with an appropriate cleansing solution.
6. Insert the catheter tip through cervix and into the uterus.
7. Inflate the balloon by attaching tip of enclosed syringe into the check valve, inject sterile saline solution and remove syringe from the check valve.
8. Pull back gently on the catheter shaft to seat the inflated balloon against the internal cervical os.
9. Under fluoroscopy inject sufficient contrast medium through the catheter into the uterus.
10. Once desired image has been recorded, deflate the balloon by firmly pressing tip of syringe into the check valve, engaging into check valve will allow emptying the content of balloon.
11. Remove the catheter from the patient when the examination is complete.

Sonohysterography (SHG)

1. Pretest the catheter balloon by first priming the enclosed syringe and then filling it with a labeled volume of sterile saline. **NOTE:** Ensure tip of syringe is contacting check valve (black plug).
2. Fill the instillation syringe with sterile saline and eliminate air from the syringe.
3. Prime the instillation lumen of the catheter with sterile saline.
4. Position a vaginal speculum within the vagina so the external cervical os can be visualized during catheter insertion.
5. Cleanse the cervix with an appropriate cleansing solution.
6. Insert the catheter tip through cervix and into the uterus.
7. Inflate the balloon by attaching tip of enclosed syringe into the check valve, inject sterile saline solution and remove syringe from the check valve.
8. Pull back gently on the catheter shaft to seat the inflated balloon against the internal cervical os.
9. Carefully remove the speculum.
10. Insert the transvaginal probe.
11. Inject sufficient saline through the catheter into the uterus.
12. Perform the ultrasound procedure.
13. Once desired image has been recorded, deflate the balloon by firmly pressing tip of syringe into the check valve, engaging into check valve will allow emptying the content of balloon.
14. Remove the catheter from the patient when the examination is complete.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCE

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.