**Guardia™ Obturator**

**Obturátor Guardia™**

**Guardia™ obturator**

**Guardia™ Obturator**

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T \_ J \_ G O \_ R E V O

## Guardia™ Obturator

Used to supplement and assist uterine access of an approved embryo transfer device for placement of in vitro fertilized (IVF) embryos into the uterine cavity. Intended for one-time use.

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

**NOTE: These devices are intended to support the use of an approved embryo transfer catheter. These tools are not intended to directly transfer gametes/embryos into the uterine cavity; rather, they are to aid in traversing a difficult or stenotic cervix or to perform a trial/mock transfer to better identify the uterine environment.**

**CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is damaged.**

### CONTRAINDICATIONS

This catheter should not be used on a patient with an active vaginal or intrauterine infection, a sexually transmitted disease, a recent uterine perforation, a recent pregnancy (or is currently pregnant), or if the patient currently has an intrauterine device.

### PRECAUTIONS

All assisted reproduction catheters used to deliver gametes or embryos into the uterine cavity should be comprised of embryo compatible materials.

Infection may occur due to bacterial contamination of the device during vaginal manipulation and result in urinary tract infection (UTI), pelvic inflammatory disease (PID), or uterine infection. Recommendations to minimize occurrence include use of only embryo compatible materials, flushing catheters and other accessories with sterile, compatible culture media, and closely adhering to sterile techniques.

Bleeding may occur as a result of trauma due to insertion of the catheter through the cervix and has been reported to be associated with a lower pregnancy rate. A simple and atraumatic transfer method has been noted to offer the best conditions for success.

### SUGGESTED INSTRUCTIONS FOR USING EMBRYO TRANSFER CATHETER ASSIST DEVICE

1. Place patient in a lithotomy or knee/chest position.
2. Insert vaginal speculum to expose cervix.
3. Loosen the fitting on the proximal hub and adjust the wire (see figure A for reference) so that the distal tip of the obturator is even with the distal tip of the guide catheter. (see figure B for reference)
4. Tighten the fitting so that the wire does not move during manipulation. (see figure A for reference)


Figure A


Figure B

5. Use obturator with guide catheter to traverse cervix if necessary. Remove obturator prior to placement of trial/mock catheter. (May be performed under ultrasound guidance if desired.)
6. Introduce trial/mock catheter through guide sheath.
7. Remove trial/mock catheter from guide sheath.
8. If so desired, move on to embryo transfer using approved catheter and refer to instructions for use for that approved catheter.
9. Upon completion of trial/mock procedure, remove the catheter/catheters, and discard.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.