EN 3	Guardia™ Obturator Instruction for Use
CS 4	Obturátor Guardia™ Návod k použití
DA 6	Guardia™ obturator Brugsanvisning
DE 7	Guardia™ Obturator Gebrauchsanweisung
EL 9	Επιπωματικό Guardia™ Οδηγίες χρήσης
ES 11	Obturador Guardia™ Instrucciones de uso
FR 13	Obturateur Guardia™ Mode d'emploi
HU 14	Guardia™ obturátor Használati utasítás
IT 16	Otturatore Guardia™ Istruzioni per l'uso
NL 18	Guardia™ obturator Gebruiksaanwijzing
NO 19	Guardia™ obturator Bruksanvisning
PL 21	Mandryn Guardia™ Instrukcja użycia
PT 23	Obturador Guardia™ Instruções de utilização
SV 24	Guardia™ obturator Bruksanvisning

GUARDIA™ OBTURATOR

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

INTENDED USE / INDICATIONS FOR USE

Used to supplement and assist uterine access of a cleared, dimensionally compatible embryo transfer device for placement of in vitro fertilized (IVF) embryos into the uterine cavity.

CONTRAINDICATIONS

This obturator should not be used in the presence of:

- Active vaginal or intrauterine infection
- Sexually transmitted disease
- Recent uterine perforation
- Recent or current pregnancy
- Confirmed or suspected intrauterine device

PRECAUTIONS

- All assisted reproduction catheters used to deliver gametes or embryos into the uterine cavity should be comprised of embryo compatible materials.
- Infection may occur due to bacterial contamination of the device during vaginal manipulation and result in urinary tract infection (UTI), pelvic inflammatory disease (PID), or uterine infection. Recommendations to minimize occurrence include use of only embryo compatible materials, flushing catheters and other accessories with sterile, compatible culture media, and closely adhering to sterile techniques.
- Bleeding may occur as a result of trauma due to insertion of the catheter through the cervix and has been reported to be associated with a lower pregnancy rate. A simple and atraumatic transfer method has been noted to offer the best conditions for success.

PRODUCT RECOMMENDATIONS

NOTE: These devices are intended to support the use of an approved embryo transfer catheter. These tools are not intended to directly transfer gametes/embryos into the uterine cavity; rather, they are to aid in traversing a difficult or stenotic cervix or to perform a trial/mock transfer to better identify the uterine environment.

	Guardia™ Access and Guardia™ AccessET	Guardia™ Access Nano	Soft-Pass™	Pivet
J-UOB-2828		X		
J-UOB-4028	X		X	X
J-UOB-4028-ST	X		X	X

INSTRUCTIONS FOR USE

1. Place patient in a lithotomy or knee/chest position.
2. Insert vaginal speculum to expose cervix.

- Loosen the fitting on the proximal hub and adjust the wire (see figure A for reference) so that the distal tip of the obturator is even with the distal tip of the guide catheter. (see figure B for reference)
- Tighten the fitting so that the wire does not move during manipulation. (see figure A for reference)


Figure A


Figure B

- Use obturator with guide catheter to traverse cervix if necessary. Remove obturator prior to placement of trial/mock or transfer catheter. (May be performed under ultrasound guidance if desired.)
- If so desired, move on to trial/mock procedure or embryo transfer using approved catheter and refer to instructions for use for that approved catheter.
- Upon completion of procedure, remove the catheter/catheters, and discard.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd