EN
3

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Brugsanvisning

T _ J - S P A D _ R E V 4

EMBRYO TRANSFER CATHETER ASSIST DEVICE

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Guardia™ Access and Soft-Pass™ Embryo Transfer Catheter Assist Devices are comprised of obturators.

INTENDED USE / INDICATION FOR USE

Used to supplement and assist uterine access of a cleared, dimensionally compatible Embryo Transfer Device for placement of in vitro fertilized (IVF) embryos into the uterine cavity.

CONTRAINDICATIONS

This device should not be used in the presence of:

- Active pelvic infection
- Sexually transmitted disease
- Pregnancy
- Intrauterine device
- Uterine perforation

WARNINGS

None known

PRECAUTIONS

- All assisted reproduction catheters used to deliver gametes or embryos into the uterine cavity should be comprised of embryo compatible materials.
- Infection may occur due to bacterial contamination of the device during vaginal manipulation. This may result in urinary tract infection (UTI), pelvic inflammatory disease (PID), or uterine infection. Recommendations to minimize occurrence include use of only embryo compatible materials, flushing catheters and other accessories with sterile, compatible culture media, and closely adhering to sterile techniques.
- Bleeding may occur as a result of trauma due to insertion of the catheter through the cervix and has been reported to be associated with a lower pregnancy rate. A simple and atraumatic transfer method has been noted to offer the best conditions for success.
- Embryo transfer procedure should only be performed by physicians trained in Assisted Reproductive Techniques.

OBTURATOR PRODUCT RECOMMENDATIONS

	Guardia™ Access and Guardia™ AccessET	Guardia™ Access Nano	Soft-Pass™	Pivet
K-JET-7002	X			X
J-SP-4420			X	

INSTRUCTIONS FOR USE

1. Place patient in a lithotomy or knee/chest position.
2. Insert vaginal speculum to expose cervix.
3. Use malleable/flexible obturator with guide catheter to traverse cervix if necessary. **NOTE:** Some obturators have markings to aid proper placement.
4. Remove malleable/flexible obturator prior to placement of trial/mock or transfer catheter. (May be performed under ultrasound guidance if desired.)
5. If so desired, move on to trial/mock procedure or embryo transfer using approved catheter and refer to instructions for use for that approved catheter.
6. Upon completion of procedure, remove the catheter/catheters, and discard.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2020
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd