- EN**
3 **Soft-Pass™ Embryo Mock Transfer Catheter Sets**
Instructions for Use
- CS**
4 **Soupravy katetrů Soft-Pass™ k provedení imitačního transferu embryí**
Návod k použití
- DA**
5 **Soft-Pass™ mock-/embryotransfereringskatetersæt**
Brugsanvisning
- DE**
6 **Soft-Pass™ Embryo-Probetransferkatheter-Sets**
Gebrauchsanweisung
- EL**
8 **Σετ καθετήρων εικονικής μεταφοράς εμβρύων Soft-Pass™**
Οδηγίες χρήσης
- ES**
9 **Equipos de catéter de transferencia de embriones de práctica Soft-Pass™**
Instrucciones de uso
- FR**
10 **Sets de cathéter de transfert pour simulation d'embryons Soft-Pass™**
Mode d'emploi
- HU**
12 **Soft-Pass™ embrióátültető gyakorló katéterkészletek**
Használati utasítás
- IT**
13 **Set con catetere per embrio transfer di prova Soft-Pass™**
Istruzioni per l'uso
- NL**
14 **Soft-Pass™ embryo proeftransferkathetersets**
Gebruiksaanwijzing
- NO**
15 **Soft-Pass™ forsøkskatetersett for embryooverføring**
Bruksanvisning
- PL**
17 **Zestawy cewnika do przenoszenia próbnego zarodków Soft-Pass™**
Instrukcja użycia
- PT**
18 **Conjuntos de cateteres para transferência de simulação de embriões Soft-Pass™**
Instruções de utilização
- SV**
19 **Soft-Pass™ Embryoöverföringskateterset för simulering**
Bruksanvisning

SOFT-PASS™ EMBRYO MOCK TRANSFER CATHETER SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Soft-Pass™ Embryo Mock Transfer Catheter Sets consist of a matched set guide and mock transfer catheter. The mock transfer catheters have a closed tip and are available with or without an EchoTip band. Both guide and mock transfer catheters have markings to assist with placement.

INTENDED USE / INDICATION FOR USE

The Soft-Pass™ Embryo Mock Transfer Catheter Sets are used for mock (simulated) embryo transfer procedures to aid in verifying catheter access and positioning in advance of an embryo transfer procedure.

CONTRAINDICATIONS

This device should not be used in the presence of:

- Active pelvic infection
- Sexually transmitted disease
- Pregnancy
- Intrauterine device
- Uterine perforation

WARNINGS

None known

PRECAUTIONS

- Infection may occur due to bacterial contamination of the device during vaginal manipulation. This may result in urinary tract infection (UTI), pelvic inflammatory disease (PID), or uterine infection. Recommendations to minimize occurrence include use of only embryo compatible materials, flushing catheters and other accessories with sterile, compatible culture media, and closely adhering to sterile techniques.
- Bleeding may occur as a result of trauma due to insertion of the catheter through the cervix and has been reported to be associated with a lower pregnancy rate. A simple and atraumatic transfer method has been noted to offer the best conditions for success.
- Embryo transfer procedure should only be performed by physicians trained in Assisted Reproductive Techniques.

INSTRUCTIONS FOR USE

1. Place patient in a lithotomy or knee/chest position.
2. Insert vaginal speculum to expose cervix.
3. If desired, place obturator (4.4 French or less in diameter) into guide catheter to aid in traversing the cervix.
4. Insert guide catheter through cervix. **NOTE:** The guide catheter has markings to aid proper placement.
5. If applicable, remove obturator prior to placement of the mock transfer catheter. **NOTE:** Ultrasound can be used to locate the tip of the transfer catheter. Some catheters are equipped with EchoTip bands for enhanced visualization under ultrasound.
6. Insert mock transfer catheter into the guide catheter. **NOTE:** The mock transfer catheter has markings to aid proper placement. Ink marks near hub of delivery/transfer catheter can be used to determine when delivery/transfer is flush with tip of guide catheter.
7. Once the mock embryo transfer procedure is complete, remove both the mock transfer and guide catheters and discard.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.