COOK
MEDICAL

CE
0123

EN
3

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Bruksanvisning

T _ J - S S G _ R E V O

SELECTIVE SALPINGOGRAPHY CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Selective Salpingography (SSG) Catheter has an angled tip with a radiopaque band.

INTENDED USE / INDICATIONS FOR USE

The Selective Salpingography Catheter is intended to be used through a uterine access device and be positioned at the tubal ostium for injection of contrast medium for fluoroscopic evaluation of tubal patency.

CONTRAINDICATIONS

Do not use this device in the presence of:

- Active uterine bleeding
- Active pelvic infection
- Sexually transmitted disease/sexually transmitted infection
- Known or suspected pregnancy
- Confirmed or suspected intrauterine device
- Gynecological malignancy
- Known allergy to dye or contrast

WARNINGS

- These catheter systems are only designed to access the tubal ostium of the fallopian tube from the uterine cavity. More distal advancement of the catheter may result in adverse effects.
- To reduce the risk of perforation, if any resistance is felt while inserting the catheter, do not force the catheter against the resistance.

PRECAUTIONS

- This device is intended for use by physicians or properly licensed practitioners trained and experienced in diagnostic and interventional techniques. Standard techniques for HSG and SSG procedures should be employed.
- The SSG Catheter is not intended for complete catheterization of the fallopian tubes or tubal recanalization. There are no data available showing any clinical benefit for use in this manner.
- Use of the SSG Catheter in postpartum patients, patients with recent cesarean section, or patients with a recent uterine perforation should be determined at the physician's discretion. Risk of uterine wall injury may be increased in these patients.

POTENTIAL ADVERSE EVENTS

- Discomfort, pain/cramping
- Bleeding
- Tubal injury/perforation
- Hypersensitivity to contrast media
- Extravasation
- Ectopic pregnancy
- Infection

INSTRUCTIONS FOR USE

1. Follow the instructions for use provided with the uterine access catheter.

NOTE: The uterine access catheter must have at least a 5.5 French working channel.

2. Under fluoroscopic guidance, inject 5-10 mL of contrast medium through the uterine access catheter to visualize the uterine cavity and tubal ostium. This volume of fluid serves to outline the basic anatomy and provide guidance in preparation for the SSG catheter placement.
3. Insert the SSG catheter through the working channel of the uterine access catheter.
4. While scanning with fluoroscopy, rotate the catheter to address the appropriate fallopian tube. The distal tip of the SSG catheter has a radiopaque band and will be visible under fluoroscopy. Advance the catheter towards the appropriate fallopian tube and place the tip of the catheter at the tubal ostium and no further. If access to the tubes is difficult, a wire guide (maximum diameter .035") may be introduced through the SSG catheter to facilitate positioning the catheter at the tubal ostium. If a wire guide was used to position the SSG catheter, remove it before moving to step 5.
5. Inject an appropriate contrast medium through the SSG catheter to assess tubal patency.

NOTE: In cases of bilateral observation, repeat steps 4-5 for other fallopian tube.

6. Once tubal evaluation has been performed, remove the SSG catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd