COOK

MEDICAL

CE
0123

EN
3

Dilator Sets

Instructions for Use

DA
4

Dilatatorsæt

Brugsanvisning

DE
6

Dilatator-Sets

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EL
8

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T _ J C D S _ R E V O

DILATOR SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

Dilator Sets include a set of serial dilators and may include an appropriately sized wire guide and access needle. The serial dilator components are designed with gradual tapers and may include hydrophilic coating, and are available in multiple inner and outer diameter and length combinations.

INTENDED USE

Intended to be used for dilating puncture sites or catheter tracts for percutaneous placement of devices for vascular and non-vascular applications such as in the venous, arterial, biliary and renal systems.

CONTRAINDICATIONS

None known

WARNINGS

- Do not use if any sign of product damage is visible.
- Do not alter the device in any way.
- This product is a delicate instrument. Avoid forceful angulation.
- Avoid manipulating or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape or shear material from the wire guide.
- Altering the tip's configuration or curve manually may damage the wire guide.

PRECAUTIONS

- This product is intended for use by clinicians trained and experienced in diagnostic and interventional techniques.
- The maximum diameter of the instrument or catheter to be introduced should be considered when selecting the dilator size.
- These wire guides come with two usable ends. Either end may be used, depending on the clinical situation as determined by an experienced operator.
- The wire guide (if present) in a Dilator Set should not be withdrawn through the needle (if present) once it has been placed beyond the needle tip.
- If resistance is encountered during dilator manipulation, stop and determine the cause before proceeding further.

POTENTIAL ADVERSE EVENTS

Potential complications include, but are not limited to:

- Air embolism
- Infection
- Hematoma
- Bleeding
- Perforation or laceration of the vessel wall

- Thrombus formation
- Pseudoaneurysm formation
- Wire guide embolization
- Vessel spasm
- Pneumothorax
- Risks normally associated with percutaneous diagnostic and/or interventional procedures

INSTRUCTION FOR USE

NOTE: If dilator has hydrophilic coating, immerse dilator in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

1. Prepare and drape the access site according to hospital protocol.
2. Introduce an access needle into target anatomy.
3. Introduce a wire guide through the needle and into the target anatomy.
4. Withdraw the needle, leaving the wire guide in place.
5. Smoothly advance the dilator over wire guide to the desired placement within the tract.
6. Remove dilator from tract while maintaining the positioning of the wire guide within the anatomy for advancement of next product to the selected place within the tract.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd