**COOK**

**MEDICAL**

**CE**  
0123

EN  
3

## **Dilator**

Instructions for Use

DA  
4

## **Dilatator**

Brugsanvisning

DE  
5

## **Dilatator**

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EL  
6

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T \_ J C D \_ R E V O

## DILATOR

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.**

### DEVICE DESCRIPTION

Dilators are designed with gradual tapers, may include hydrophilic coating, and are available in multiple inner and outer diameter and length combinations.

### INTENDED USE

Intended to be used for dilating puncture sites or catheter tracts for percutaneous placement of devices for vascular and non-vascular applications such as in the venous, arterial, biliary and renal systems.

### CONTRAINDICATIONS

None known

### WARNINGS

Do not alter the device in any way.

### PRECAUTIONS

- This product is intended for use by clinicians trained and experienced in diagnostic and interventional techniques.
- The maximum diameter of the instrument or catheter to be introduced should be considered when selecting the dilator size.
- If resistance is encountered during dilator manipulation, stop and determine the cause before proceeding further.

### POTENTIAL ADVERSE EVENTS

Potential complications include, but are not limited to:

- Air embolism
- Infection
- Hematoma
- Bleeding
- Perforation or laceration of the vessel wall
- Thrombus formation
- Pseudoaneurysm formation
- Wire guide embolization
- Vessel spasm
- Pneumothorax
- Risks normally associated with percutaneous diagnostic and/or interventional procedures

### INSTRUCTION FOR USE

**NOTE:** If dilator has hydrophilic coating, immerse dilator in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

1. Smoothly advance the dilator over a wire guide (not included) to the desired placement within the tract.
2. Remove dilator from tract while maintaining the positioning of the wire guide within the anatomy for advancement of next product to the selected place within the tract.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd