EN	Embryo Transfer Catheter
CS	Katetr pro embryotransfer
DA	Embryotransfereringskateter
DE	Embryotransferkatheter
EL	Καθετήρας μεταφοράς εμβρύων
ES	Catéter de transferencia de embriones
FR	Cathéter pour transfert d'embryons
HU	Embrióátültető katéter
IT	Catetere per embrio transfer
NL	Embryotransferkatheter
NO	Embryooverføringskateter
PL	Cewnik do przenoszenia zarodków
PT	Cateter para transferência de embriões
SV	Embryoöverföringskateter

ENGLISH

EMBRYO TRANSFER CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Guardia™ Access, Guardia™ AccessET, Guardia™ Access Nano, Pivet, Soft-Pass™, and Soft-Pass™ EchoTip. Embryo transfer catheter sets are comprised of a transfer catheter and a guide catheter.

INTENDED USE / INDICATION FOR USE

Used to place in vitro fertilized (IVF) embryos into the uterine cavity.

CONTRAINDICATIONS

This device should not be used in the presence of:

- Active pelvic infection
- Sexually transmitted disease
- Pregnancy
- Confirmed or suspected intrauterine device
- Uterine perforation

WARNINGS

None known

PRECAUTIONS

- Infection may occur due to bacterial contamination of the device during vaginal manipulation. This may result in urinary tract infection (UTI), pelvic inflammatory disease (PID), or uterine infection.
- Bleeding may occur as a result of trauma due to insertion of the catheter through the cervix.
- Embryo transfer procedure should only be performed by physicians trained in Assisted Reproductive Techniques.

INSTRUCTIONS FOR USE

One-step embryo transfer

- Place patient in a lithotomy or knee/chest position.
- Insert vaginal speculum to expose cervix.
- Insert transfer catheter into guide catheter.
- Attach syringe to luer-style fitting of transfer catheter and follow laboratory procedure for loading embryos into catheter.
- Retract transfer catheter until distal tip is level with, or inside, the distal tip of guide catheter.
- Insert catheters through the cervix.

NOTE: Some guide catheters have markings or positioners to aid proper placement. Ink marks near hub of delivery/transfer catheter can be used to determine when delivery/transfer catheter is flush with tip of guide catheter.
- Advance transfer catheter until the tip is in desired location within the uterine cavity.

NOTE: Ultrasound can be used to locate the tip of the transfer catheter. Some catheters are equipped with EchoTip® bands for enhanced visualization under ultrasound.
- Expel embryos into the uterine cavity.
- Remove catheters and check for embryo retention.

Two-step embryo transfer

- Place patient in a lithotomy or knee/chest position.
- Insert vaginal speculum to expose cervix.
- If applicable, place obturator into guide catheter to aid in traversing the cervix.
- Insert guide catheter through cervix.

NOTE: Some guide catheters have markings or positioners to aid proper placement.
- If applicable, remove obturator prior to placement of transfer catheter.
- Attach syringe to luer-style fitting of transfer catheter and follow laboratory procedure for loading embryos into catheter.
- Insert transfer catheter into guide catheter.
- Advance transfer catheter until the tip is in desired location within the uterine cavity.

NOTE: Ultrasound can be used to locate the tip of the transfer catheter. Some catheters are equipped with EchoTip bands for enhanced visualization under ultrasound.
- Expel embryos into the uterine cavity.
- Remove catheters and check for embryo retention.

NOTE: One cell MEA tested and passed with 80% or greater blastocyst development within 96 hrs. USP Endotoxin (LAL) tested and passed with 20 EU or less per device.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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