EN
2

Kumpe Access Catheter

Instructions for Use

CS
3

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DA
5

Kumpe adgangskateter

Brugsanvisning

DE
6

Kumpe Zugangskatheter

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EL
8

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T _ K A C _ R E V 1

KUMPE ACCESS CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

The Kumpe Access Catheter is composed of a polymer reinforced with a stainless steel braided wire. The proximal end has a Luer lock hub, and the distal end is formed into an angle shape. The catheter is 65 cm long and is 5.0 Fr in diameter.

NOTE: The Kumpe Access Catheter can accept up to .038 inch in diameter wire guides.

INTENDED USE

These devices are intended for access and catheterization of the urinary tract, including the following applications:

- Delivery of contrast media
- Navigation of a tortuous ureter
- Access, advancement, or exchange of wire guides

CONTRAINDICATIONS

This device is contraindicated for use in patients unsuitable for catheterization.

WARNINGS

Sterilized by Ethylene Oxide gas. Do not use if sterile barrier is damaged. The device is intended for one-time use.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in ureteral access techniques. Standard techniques should be employed.
- Avoid bending or kinking the catheter prior to placement. Doing so could damage the integrity of the catheter and result in patient injury.
- If you encounter resistance while advancing or withdrawing the catheter, stop. Determine the cause of the resistance before proceeding.
- Do not withdraw the catheter while it is deflected in the cystoscope/endoscope.

POTENTIAL ADVERSE EVENTS

Adverse events that may result from catheter placement include, but are not limited to:

- Extravasation
- Hemorrhage
- Sepsis
- Edema
- Loss of renal function
- Perforation or laceration of the kidney, renal pelvis, ureter, and/or bladder
- Peritonitis
- Urinary tract infection
- Abrasion

INSTRUCTIONS FOR USE

Placement with a Wire Guide

1. Place the cystoscope within the bladder to locate the ureteral orifice.
2. Flush the catheter of choice.
3. Advance the catheter through the cystoscope and into the bladder.

4. Advance the wire guide through the catheter. Slightly advance and retract the wire guide to assist in gaining access to the ureteral orifice and into the ureter.
5. Once the wire is advanced into the ureter, advance the catheter over the wire and to the desired location within the ureter.
6. Remove the wire guide.
7. Under fluoroscopy, inject contrast media through the catheter to visualize the ureter and to determine required intervention.

Placement without a Wire Guide

1. With the cystoscope in the bladder, locate the ureteral orifice to determine choice of catheter.
2. Flush the catheter.
3. Under direct vision, advance the catheter through the cystoscope and to the desired location within the ureter.
4. Under fluoroscopy, inject contrast media through the ureteral catheter to visualize the ureter and to determine required intervention.

Percutaneous Placement

1. Flush the catheter.
2. With a wire guide already in place within the renal pelvis, advance the catheter over the wire guide and into the renal pelvis.
3. Retract the wire guide to allow the tip of the catheter to form.
4. Advance the catheter to the desired location within the renal pelvis. Slightly advance and retract the wire guide to help navigate the catheter.
5. Once the catheter is in the desired location, deliver the working wire guide of choice through the catheter.
6. Maintain access with the wire guide while removing and discarding the catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd