EN
3

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Bruksanvisning

FLEXOR® CHECK-FLO® INTRODUCER SET MICROPUNCTURE® RADIAL ARTERY ACCESS WITH AQ® HYDROPHILIC COATING

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Flexor Check-Flo Introducer Set Micropuncture Radial Artery Access with AQ Hydrophilic Coating consists of:

- Micropuncture access needle
- Wire guide
- Flexor sheath with AQ Hydrophilic Coating and a coaxial dilator
- Syringes

INTENDED USE

The Flexor Check-Flo Introducer Set Micropuncture Radial Artery Access with AQ Hydrophilic Coating has been designed for the introduction of balloon, diagnostic and guiding catheters or other devices for diagnosis and intervention in radial artery access procedures.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in proper vascular access techniques. Standard techniques for the placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Possible allergic reactions or access site infections should always be considered.
- A sterile inflammatory response possibly associated with the use of the product in conjunction with latex and powdered non-latex based gloves has been reported with the use of this product.¹
- All catheters and instruments used with this product should move freely through the valve and sheath. Separation of the sheath may result when the fit is tight.
- Do not attempt to heat or reshape the curves as the products are made from heat sensitive material.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Prior to beginning procedure, perform the Allen Test to assess the presence and adequacy of a dual arterial circulation to the hand.
2. Apply appropriate local anesthetic agent at skin access site.
3. Using a surgical scalpel, make a shallow skin nick at the puncture site. (optional)
4. Insert the 21 gage Micropuncture needle into the vessel lumen.
5. Insert the .018 inch diameter nitinol wire guide through the cannula and into the vessel lumen.
6. Remove the needle over the wire guide.
7. Using heparinized saline, activate the hydrophilic coating of the coaxial sheath/dilator by wetting the outer surface of the device.
8. Insert the dilator/sheath combination over the wire guide and into the vessel lumen.
9. Slowly remove the dilator and wire guide together. The outer Flexor sheath will remain in the vessel.
10. The sheath may now be used for insertion and advancement of appropriately sized diagnostic or interventional devices.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

¹Kozak, et al., "Sterile Inflammation Associated with Transradial Catheterization and Hydrophilic Sheaths," *Catheterization and Cardiovascular Interventions*, June 2003; 59(number 2):207-213.

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Bloomington, IN 47404 U.S.A.