- EN**
3 **Koyle Diaper Stent**
Instructions for Use
- CS**
4 **Plenkový stent Koyle**
Návod k použití
- DA**
5 **Koyle blestent**
Brugsanvisning
- DE**
6 **Koyle Windel-Stent**
Gebrauchsanweisung
- EL**
7 **Ενδοπρόσθεση πάνας Koyle**
Οδηγίες χρήσης
- ES**
8 **Stent para pañal Koyle**
Instrucciones de uso
- FR**
9 **Endoprothèse de Koyle pour drainage entre deux couches**
Mode d'emploi
- HU**
10 **Koyle pelenkasztent**
Használati utasítás
- IT**
11 **Stent per pannolino Koyle**
Istruzioni per l'uso
- NL**
12 **Koyle luierstent**
Gebruiksaanwijzing
- NO**
13 **Koyle bleiestent**
Bruksanvisning
- PL**
14 **Stent pieluszkowy Koyle**
Instrukcja użycia
- PT**
15 **Fralda stent Koyle**
Instruções de utilização
- SV**
16 **Koyle blöjstent**
Bruksanvisning

KOYLE DIAPER STENT

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

This device is constructed of C-Flex tubing and a C-Flex suture funnel. The tubing is 50 cm in length and is available in diameters of 6 and 8 French. The tubing has four or eight distal sideports. The suture funnel is located in between the distal and proximal tip. One version of the 6 French stent has ink marks from the distal tip to the suture funnel.

INTENDED USE

This device is used for stenting the urethra during hypospadias or epispadias repair and to allow postoperative drainage of the bladder in infants (29 days to less than 2 years old) and children (2 years old to less than 12 years old).

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions that create unacceptable risk during catheterization.

WARNINGS

None known

PRECAUTIONS

This device is intended for use by physicians trained and experienced in urological techniques. Standard techniques should be employed.

POTENTIAL ADVERSE EVENTS

- Urinary tract infection
- Perforation or laceration of the urethra
- Perforation or laceration of the bladder
- Abrasions

INSTRUCTIONS FOR USE

1. Measure the urethra and cut the long segment of the stent to the correct length.
2. Advance the stent into the urethra.
3. Use absorbable suture to anchor the stent to the glans meatus at the funneled area.
4. Apply appropriate dressing and place the diaper on the patient, leaving the distal segment of the stent outside the diaper.
5. Place an additional, larger outside diaper over the inner diaper, and route the distal drainage segment of the stent into the outer diaper.

NOTE: The stent must not remain indwelling for more than 21 days. During the healing process, the sutures will be absorbed, and the stent should be expelled spontaneously during micturition.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2020
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd