EN **Intrauterine Access Balloon Catheter**

3 Instructions for Use

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24 Bruksanvisning

INTRAUTERINE ACCESS BALLOON CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Intrauterine Access Balloon Catheter is a dual lumen polymer catheter with a working channel, silicone balloon and packaged with shipping stylet and syringe.

INTENDED USE / INDICATIONS FOR USE

The Intrauterine Access Balloon Catheter is used for delivery of contrast medium into the uterine cavity during hysterosalpingography (HSG) for examination of the uterus and fallopian tubes, and evaluation of tubal patency. The catheter lumen accepts catheters up to 6.0 French in diameter to access the uterine cavity for assessment of tubal patency.

CONTRAINDICATIONS

This device should not be used in the presence of:

- Active uterine bleeding
- Active pelvic infection
- Sexually transmitted disease/sexually transmitted infection
- Known or suspected pregnancy
- Known allergy to dye or contrast medium
- Confirmed or suspected intrauterine device

WARNINGS

- Always inflate the balloon with sterile liquid; never inflate with air, carbon dioxide or any other gas.
- Do not exceed recommended balloon volume of 2 mL. Over-inflation of this device can cause the balloon to rupture.
- To reduce the risk of perforation, if any resistance is felt while inserting the catheter, do not force the catheter against the resistance.

PRECAUTIONS

- This procedure should only be performed by a physician or properly licensed practitioner trained in diagnostic and interventional techniques. Standard techniques for HSG procedures should be employed.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.
- Use of the Intrauterine Access Balloon Catheter in postpartum patients, patients with recent cesarean section, or patients with a recent uterine perforation should be determined at the physician's discretion. Risk of uterine wall injury may be increased in these patients.

POTENTIAL ADVERSE EVENTS

- Discomfort, pain/cramping
- Bleeding
- Infection
- Uterine perforation
- Hypersensitivity to contrast media

INSTRUCTIONS FOR USE

1. Remove the shipping stylet from catheter.

2. Pretest the catheter balloon by first priming the enclosed syringe and then filling it with a maximum of 2 mL of medium. Instill medium, check to ensure the balloon inflates, and then deflate the balloon.
3. Fill the instillation syringe with contrast medium and eliminate air from the syringe.
4. Prime the instillation lumen of the catheter with contrast medium.
5. Place patient in dorsal lithotomy position.
6. Position a vaginal speculum within the vagina so the external cervical os can be visualized during catheter insertion.
7. Cleanse the cervix with an appropriate cleansing solution.
8. Insert the catheter tip through cervix and into the uterus.
9. Inflate the balloon by instilling sterile saline, making sure not to exceed 2 mL.
10. Pull back gently on the catheter shaft to seat the inflated balloon against the internal cervical os.
11. Inject sufficient contrast medium through the catheter into the uterus.
12. Perform the fluoroscopy procedure.

NOTE: If spillage of contrast medium from the fimbriae portion of each tube into the peritoneal cavity is observed, this confirms tubal patency, and the procedure is complete.
13. If spillage of contrast medium from the fimbriae portion of the tube into the peritoneal cavity is not observed in 1 or both of the tubes the working channel of the catheter may be used to insert a catheter that is less than 6 French in diameter for further assessment of tubal patency, i.e., selective salpingography.
14. Deflate balloon and remove the catheter from the patient when the examination is complete and discard.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd