
- EN**  
5 **Modified Novy Cornual Cannulation Set**  
Instructions for Use
- CS**  
8 **Upravená souprava Novy pro kanylaci děložních rohů**  
Návod k použití
- DA**  
12 **Modifieret Novy cornualt kanyleringsset**  
Brugsanvisning
- DE**  
15 **Modifiziertes Novy Set zur Kanulierung der Cornua uteri**  
Gebrauchsanweisung
- EL**  
19 **Τροποποιημένο σετ καθετηριασμού Novy Cornual**  
Οδηγίες χρήσης
- ES**  
23 **Equipo de canulación cornual Novy modificado**  
Instrucciones de uso
- FR**  
27 **Set de Novy pour cathétérisme des cornes**  
Mode d'emploi
- HU**  
30 **Módosított Novy cornualis kanülálókészlet**  
Használati utasítás
- IT**  
34 **Set per incannulamento del corno uterino Novy modificato**  
Istruzioni per l'uso
- NL**  
38 **Aangepaste Novy uterushoorncanulatieset**  
Gebruiksaanwijzing
- NO**  
41 **Modifisert Novy kornualt kannuleringssett**  
Bruksanvisning
- PL**  
45 **Zmodyfikowany zestaw Novy do kaniulacji rogu macicy**  
Instrukcja użycia
- PT**  
49 **Conjunto de canulação Novy Cornual modificado**  
Instruções de utilização
- SV**  
52 **Modifierat Novy kanyleringsset för tubarhörn**  
Bruksanvisning



1


2

3

4

## MODIFIED NOVY CORNUAL CANNULATION SET

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### DEVICE DESCRIPTION

Set contains:

- Introducing catheter with side-arm fitting and Tuohy-Borst adapter
- Stiffening obturator
- Inner catheter with 1 cm markings and Tuohy-Borst adapter
- Wire guide with removable silicone safety cap
- Luer lock cap

### INTENDED USE/INDICATION FOR USE

The Modified Novy Cornual Cannulation Set is intended for use through the operating channel of a hysteroscope/laparoscope or other uterine access device, for hysteroscopic or fluoroscopic selective catheterization and cannulation of the proximal fallopian tube(s), followed by the introduction of chromotubation solution or contrast medium, in the evaluation of tubal patency.

**The wire guide in the Modified Novy Cornual Cannulation Set is intended only to facilitate placement of the inner catheter. It is not intended for tubal recanalization and should not be advanced beyond the tubal isthmus.**

### CONTRAINDICATIONS

- Pregnancy or suspected pregnancy
- Active or recent pelvic infection
- Unexplained or severe uterine bleeding
- Gynecologic malignancy (suspected or known)
- Known allergy to dye or contrast medium

### WARNINGS

- This device should not be used in the presence of hemorrhage, active pelvic infection, profuse bleeding, or pregnancy.
- This catheter system is only designed to access the fallopian tube from the uterine cavity. It should only be advanced as far as necessary to achieve a meaningful diagnostic study. Possible adverse effects of more distal advancement have not yet been determined. The catheter system is not intended for tubal recanalization, and there is no data available, to date, showing any clinical benefit when used this way.
- This catheter system is not designed nor intended for embryo transfer or other assisted reproduction procedures.

### PRECAUTIONS

- The product is intended for use by physicians trained and experienced in gynecological techniques. Standard techniques should be employed.
- The Modified Novy Cornual Cannulation Set is not intended for complete catheterization of the fallopian tube or tubal recanalization.

### POTENTIAL ADVERSE EVENTS

The following adverse events have been reported or are of concern with the use of this type of device. Patients must be advised of the following risks.

- Ectopic pregnancy
- Damage to normal tubes
- Tubal dissection
- Pain or discomfort
- Extravasation
- Tubal perforation
- Bleeding
- Infection

## **INSTRUCTIONS FOR USE**

### **With Hysteroscopy**

**NOTE:** The Modified Novy Cornual Cannulation Set should be loaded into the hysteroscope before beginning the examination, ensuring that the introducing catheter may be moved freely within the hysteroscope.

1. Advance the introducing catheter, with the obturator in place, into the operating channel of the hysteroscope.
2. Remove the stiffening obturator from the introducing catheter and occlude the straight arm of the fitting with the Luer lock cap provided. This port is used for irrigation or suction and injection of chromotubation solution or contrast medium into the uterine cavity.
3. Loosen the Tuohy-Borst adapter on the inner catheter and introduce the wire guide into it, positioning the distal tip of the wire guide evenly with the tip of the catheter. Tighten the Tuohy-Borst adapter to fix the position of the wire guide.
4. Loosen the Tuohy-Borst adapter on the side-arm of the introducing catheter fitting and advance the inner catheter, with the wire guide in place, through the introducing catheter until both catheter tips are even. For future reference, note the position of the proximal distance markings on the inner catheter with respect to the Tuohy-Borst adapter on the side-arm fitting. The first marking should be flush with the Tuohy-Borst adapter.
5. Remove the inner catheter and wire guide from the introducing catheter, and tighten the Tuohy-Borst adapter on the side-arm fitting to seal it and prevent leakage of fluid or gas.
6. Position the introducing catheter in the operating channel of the hysteroscope so it does not extend beyond the telescope lens.
7. Introduce the hysteroscope sheath, with its blunt obturator in place, through the dilated cervix and into the uterine cavity.
8. Remove the obturator from the hysteroscope sheath.
9. Introduce the hysteroscope, with the introducing catheter in place, into the hysteroscope sheath.
10. Under direct vision, guide the hysteroscope toward the tubal ostium.
11. After the tubal ostium is visualized, direct the introducing catheter toward the tubal orifice and wedge the tip of the catheter in the tubal ostium.
12. Remove the Luer lock cap from the straight arm of the side-arm fitting and attach a 3 mL syringe filled with chromotubation solution or dilute contrast medium (25% solution).
13. Inject the chromotubation solution or contrast medium into the tubal ostium while observing laparoscopically for spillage from the fimbrial portion of the tube. If such spillage is observed, tubal patency has been confirmed and the procedure is complete.
 

**NOTE:** If spillage of chromotubation solution or contrast medium from the fimbrial portion of the tube into the peritoneal cavity is not observed, it will be necessary to introduce the inner catheter into the proximal fallopian tube to confirm tubal blockage or rule out a false positive result due to temporary tubal spasm. Please proceed with step 14.
14. Detach the 3 mL syringe and occlude the straight arm of the side-arm fitting with the Luer lock cap.
15. Loosen the Tuohy-Borst adapter on the side-arm fitting and advance the inner catheter, with the wire guide in place, through the introducing catheter until the first proximal marking on the inner catheter is flush with the Tuohy-Borst adapter. The distal tip of the inner catheter should now be even with that of the introducing catheter.
16. Taking care to allow movement of the inner catheter, tighten the Tuohy-Borst adapter on the side-arm fitting to form a gas-tight seal around the inner catheter.

17. Begin cannulation of the intramural portion of the fallopian tube with the wire guide and inner catheter by first advancing the wire guide (**Fig. 1**). The flexibility of the wire guide tip will increase as the length exiting the catheter tip increases.
18. Advance the wire guide beyond the utero-tubal junction and into the tubal isthmus unless significant resistance is met (**Fig. 2**). Advancement of the wire guide through the utero-tubal junction is observed directly by laparoscopy and video display of the proximal oviduct or by characteristic arching of the tubal isthmus. Advancement may be facilitated by displacing or straightening the fallopian tube with an instrument under laparoscopic guidance.
19. Coaxially advance the inner catheter over the positioned wire guide for a short distance or until resistance is felt (**Figs. 3 and 4**). The distance the catheter passes into the tubal isthmus may be determined by observing the distance markings on either the proximal or distal portion of the catheter.

**WARNING: If significant resistance is felt, do not attempt to advance the catheter.**

**WARNING: Do not attempt to advance the catheter or wire guide beyond the tubal isthmus.**

20. Remove the wire guide and attach a 3 mL syringe filled with chromotubation solution or dilute contrast medium (25% solution).
21. Inject the chromotubation solution or contrast medium into the tubal ostium while observing laparoscopically for spillage from the fimbrial portion of the tube. If such spillage is observed, tubal patency is confirmed.

### With Fluoroscopy

1. Remove the stiffening obturator from the introducing catheter and occlude the straight arm of the fitting with the Luer lock cap provided. This port is used for irrigation or suction and injection of contrast medium into the uterine cavity.
  2. Loosen the Tuohy-Borst adapter on the inner catheter and introduce the wire guide into it, positioning the distal tip of the wire guide evenly with the tip of the catheter. Tighten the Tuohy-Borst adapter to fix the position of the wire guide.
  3. Loosen the Tuohy-Borst adapter on the side-arm of the introducing catheter fitting and advance the inner catheter, with the wire guide in place, through the introducing catheter until both catheter tips are even. For future reference, note the position of the proximal distance markings on the inner catheter with respect to the Tuohy-Borst adapter on the side-arm fitting. The first marking should be flush with the Tuohy-Borst adapter.
  4. Remove the inner catheter and wire guide from the introducing catheter, and tighten the Tuohy-Borst adapter on the side-arm fitting to seal it and prevent leakage.
  5. Introduce the appropriate uterine access device into the cervix and affix it in position either on the cervix or within the uterine cavity.
  6. Insert the introducing catheter through the uterine access device and advance the catheter until its distal tip is positioned in the lower portion of the uterine cavity.
  7. Remove the Luer lock cap from the straight arm of the side-arm fitting and attach a syringe filled with half-strength contrast medium. With fluoroscopic visualization, inject contrast medium through the introducing catheter to aid visualization of the uterine cavity and tubal ostium.
  8. Detach the syringe and introduce a .035 inch Safe-T-J<sup>®</sup> wire guide through the introducing catheter and into the uterine cavity. With the wire guide leading and under fluoroscopic guidance, advance the wire guide and catheter together toward the tubal ostium and wedge the tip of the catheter in the tubal ostium.
  9. Remove the wire guide from the introducing catheter and attach a syringe filled with full-strength contrast medium. Inject 2 to 10 mL of contrast medium into the tubal ostium while observing fluoroscopically for spillage from the fimbrial portion of the tube. If such spillage is observed, tubal patency has been confirmed and the procedure is complete.
- NOTE:** If spillage of contrast medium from the fimbrial portion of the tube into the peritoneal cavity is not observed, it will be necessary to introduce the inner catheter into the proximal fallopian tube to confirm tubal blockage or rule out a false positive result due to temporary tubal spasm. Please proceed with step 10.
10. Detach the syringe and occlude the straight arm of the side-arm fitting with the Luer lock cap.
  11. Loosen the Tuohy-Borst adapter on the side-arm fitting and advance the inner catheter, with the wire guide in place, through the introducing catheter until the first proximal marking on the inner catheter is flush with the Tuohy-Borst adapter. The distal tip of the inner catheter should now be even with that of the introducing catheter.

12. Taking care to allow movement of the inner catheter, tighten the Tuohy-Borst adapter on the side-arm fitting to form a fluid-tight seal around the inner catheter.
13. Begin cannulation of the intramural portion of the fallopian tube with the wire guide and inner catheter by first advancing the wire guide (**Fig. 1**). The flexibility of the wire guide tip will increase as the length exiting the catheter tip increases.
14. Advance the wire guide beyond the utero-tubal junction and into the tubal isthmus unless significant resistance is met (**Fig. 2**). Advancement of the wire guide through the utero-tubal junction may be observed fluoroscopically.
15. Coaxially advance the inner catheter over the positioned wire guide for a short distance or until resistance is felt (**Figs. 3 and 4**). The distance the catheter passes into the tubal isthmus may be determined by observing the distance markings on the proximal portion of the catheter.  
**WARNING: If significant resistance is felt, do not attempt to advance the catheter.**  
**WARNING: Do not attempt to advance the catheter or wire guide beyond the tubal isthmus.**
16. Remove the wire guide and attach a syringe filled with full-strength contrast medium.
17. Inject the contrast medium into the tubal ostium while observing fluoroscopically for spillage from the fimbrial portion of the tube. If such spillage is observed, tubal patency is confirmed.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

M. J. Novy, A. S. Thurmond, P. Patton, B. T. Uchida, J. Rösch: "Diagnosis of Cornual Obstruction by Transcervical Fallopian Tube Cannulation," *Fertility and Sterility*, 50 (1988), 434-440.5.

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Bloomington, IN 47404 U.S.A.